RESPIRONICS

Instructions for Use

AF531 Oro-Nasal Single-Use Leak 1

Intended Use

Medium and Large Size: The AF531 Oro-Nasal EE Leak 1 Face Mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single use in the hospital/institutional environment only. The mask is to be used on patients (>30 kg) for whom CPAP or bi-level therapy has been prescribed.

Small Size: The AF531 Oro-Nasal EE Leak 1 Face Mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single use in the hospital/institutional environment only. The mask is to be used on patients 7 years or older (>20 kg) for whom CPAP or bi-level therapy has been prescribed.

 **Note:** This mask is not made with natural rubber latex or DEHP.

 **Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.

Contraindications

This mask may not be suitable for patients with the following conditions: hiatal hernia, impaired cardiac sphincter function, excessive reflux, impaired cough reflex; or patients who are unable to remove the mask by themselves. This mask should also not be used on patients who are taking a prescription drug that may cause vomiting or on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.

Warnings:


- This mask is not suitable for providing life support ventilation.
- An exhalation port is not built into the AF531 Oro-Nasal EE Leak 1 Face Mask. A separate exhalation device must be used with this mask.
- Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. Do not block or otherwise try to seal the exhaust vent on the exhalation port.
 - * **Explanation of the Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation.
- A minimum of 3 cm H₂O (hPa) must be maintained when using this mask. Some rebreathing may occur.
 - * At low CPAP or EPAP pressures, the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing.
 - * If an additional exhalation device is added to the patient circuit, you may need to adjust the pressure level to compensate for the additional leak of the exhalation device.
- This mask, circuit and filter are for single use only and are not intended to be disinfected or cleaned after patient use. Dispose of the mask and circuit components after patient use. Respironics cannot guarantee the performance specifications if the mask, circuit, or filter are reused.

PHILIPS

- Check the patient regularly for adequate ventilation:
 - * Before instituting mechanical ventilation, confirm proper connection of the circuit and verify ventilator operation. Patient, ventilator, and circuit must be monitored on a regular basis per established standards of care.
 - * Do not block or otherwise try to seal the exhaust vent on the exhalation port. The exhaust vent on the exhalation port is designed to exhaust CO₂ from the patient circuit and continuous flow is required for safe operation. Inspect and verify the proper operation of the exhalation port regularly during use.
 - * To prevent increased expiratory resistance due to an occluded exhalation port filter, which may be caused by environmental conditions or use of active humidification, monitor the filter closely and replace it per your institution's protocol.
- Monitor or intervene if the patient experiences skin redness, irritation, or discomfort.
- The patient's physician should be contacted if the patient experiences the following symptoms while using the mask or after removing it:
 - * Unusual chest discomfort,
 - * Shortness of breath,
 - * Stomach distension,
 - * Belching, or severe headache;
 - * Drying of the eyes, eye pain, or eye infections; blurred vision. (Consult an ophthalmologist if symptoms persist.)
- Monitor or intervene if any of the following symptoms occur: tooth, gum, or jaw soreness. Use of a mask may aggravate a patient's existing dental conditions.
- Do not pull or stretch the tubing. This could result in circuit leaks.
- The mask contains small parts which could result in a choking hazard.


Mask Parts and Mask Options

Mask Frame and Elbow


- Ⓐ Mask Frame
- Ⓑ Forehead Adjuster
- Ⓒ NIVO Nebulizer Port  **Note: Consult Nebulizer Instructions for Use.**
- Ⓓ Entrainment Elbow with Anti-Asphyxia Valve
- Ⓔ Headgear Clips

CapStrap Headgear


- (F) CapStrap Headgear Crown Strap
- (G) CapStrap Forehead Pad
- (H) Headgear Top Straps
- (I) Headgear Bottom Straps
- (J) Headgear Clips

Four Point Headgear

- (K) Four Point Headgear Crown Strap
- (L) Headgear Top Straps
- (M) Headgear Bottom Straps
- (N) Headgear Clips

Optional Circuit (showing AF531 connection)

- ① Bacteria Filter
- ② Proximal Pressure Line
- ③ Smooth Inner Lumen Tubing
- ④ Disposable Filterable Exhalation Port (FEP)
- ⑤ Proximal Pressure Port

Before Use Read and Understand the Instructions Completely

⚠ Caution: Any deviation from these instructions may impact the performance of the product.

⚠ Caution: Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.

⚠ Caution: If the mask becomes soiled, please discard and replace.

- Wash the patient's face.
- ⚠ Caution:** Do not use moisturizer/lotion on your hands or the patient's face.
- Verify that the mask and headgear are the correct size.
- If a nasogastric tube or similar device is in place, use the optional nasogastric tube sealing pad. Position the pad so that its flat surface is against the patient's face and the C-shaped opening surrounds the tube.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Check the Entrainment Valve.
- Verify therapy device pressure(s).

Entrainment Valve

Locate the air entrainment valve flapper inside the mask elbow and verify it functions. Visually inspect the anti-asphyxia valve (flapper). The flapper inside the elbow must lay flat when the device is shut off so that room air can flow through the large opening in the valve (Figure 1). This allows the patient to breathe fresh air. Replace the mask if the flapper does not lay flat when the device is shut off.

⚠ Warning: Do not block the opening on the entrainment valve or the exhalation device. Check the valve before each use for malfunction or any blockage by patient secretions; ensure that the flapper is dry. To ensure proper air entrainment valve function, verify that the flapper is free-moving.


Figure 1

Leak Symbol and Port Settings

The leak symbol  printed on your elbow may be used with ventilators that are equipped with Interface Selection settings. On this type of ventilator, enter the leak symbol value (1) that corresponds with the symbol shown on the elbow. This value represents the intentional leak characteristics of the mask interface and enables the ventilator to calculate and monitor patient leak (intentional or unintentional). If using the mask and circuit, set the exhalation port selection to DEP.

Achieving the Right Fit

To get started, use the sizing gauge located on the AF531 package to select the proper mask size.

1. Hold the mask lightly against the patient's face and slide the headgear over their head.
 - If using the CapStrap Headgear, connect the bottom clips to the mask with a gentle push (Figure 2).
 - If using the Four Point Headgear strap, connect the bottom clips to the mask with a gentle push (Figure 3).
 - Disengage the headgear tabs from the bottom clips and tighten the straps to achieve a comfortable fit with minimal leak.
2. Adjust the Four Position Forehead Arm to ensure minimal mask contact pressure on the nasal bridge while achieving minimal mask cushion leak.
 - The lowest position provides the highest contact pressure at the bridge of the nose and provides the best mask cushion seal.
 - The highest position provides the lowest contact pressure at the bridge of the nose and provides lower mask cushion seal.
3. Connect the patient circuit (exhalation device and flexible tubing) to the entrainment valve.
4. Unfasten the headgear tabs on the top strap and gradually adjust to ensure the mask is tight enough to achieve minimal mask cushion leak and a comfortable fit. **DO NOT OVERTIGHTEN.** Adjust the crown strap as necessary to achieve proper mask support.


Figure 2

- 1 CapStrap Headgear
- 2 Four Position Forehead Arm
- 3 Entrainment Valve (EE)
- 4 Mask Faceplate
- 5 Top Headstrap


Figure 3

- 1 Four Point Headgear Strap
- 2 Four Position Forehead Arm
- 3 Entrainment Valve (EE)
- 4 Mask Faceplate
- 5 Connect to Forehead Bracket
- 6 Crown Strap - starts at top of head and around back of head

5. Turn on the noninvasive ventilator. Instruct the patient to breathe normally.
6. Make final adjustments to all adjustment points to lessen leaks and to provide a comfortable fit.

 **Note:** Re-adjust the strap tension if leaks occur as the patient's position changes. Do not overtighten the straps. Overtightening can cause or worsen leaks.

Connecting and Disconnecting the Mask and Headgear

The AF531 Oro-Nasal EE Leak 1 Face Mask can be used with either the CapStrap Headgear or the Four Point Headstrap. Use the following instructions to disconnect and reconnect the headgear to the mask.

Removing the CapStrap Headgear

First, disconnect one or both of the bottom clips. Next, gently twist the CapStrap hook so that it disconnects from the forehead support bracket (Figure 4).

Replacing the CapStrap Headgear

Press the CapStrap hook into the back of the slot on the forehead support bracket and press until they snap together. On the lower straps, press the bottom clips into place on the faceplate.

Removing the Four Point Headstrap

First, undo the bottom clips. Next, open the headgear tabs and pull them out of the slots on the forehead support arm.

Replacing the Four Point Headstrap

On the top strap, thread the headgear tabs through the slots in the forehead support arm (Figure 5). Fold the tabs back so that they catch onto the headgear straps.

Removing the Mask and Headgear from the Patient

Disconnect the bottom clips and swing the mask upward (Figure 6). This allows easy access to the patient. To remove completely, pull the entire assembly over the patient's head.

Replacing the Mask Elbow

To remove: Grasp and squeeze the elbow tabs and pull the elbow away from the mask faceplate.

To replace: Push the elbow onto the mask faceplate until the tabs snap into place (audible click). Do not use excessive force or any tools.

Make sure the elbow is firmly in place and will not separate from the hub during normal use.

Connecting the Nebulizer to the Mask Elbow

Disconnect the bottom clips and swing the mask upward (Figure 6). This allows easy access to the patient. To remove completely, pull the entire assembly over the patient's head. Remove the Nebulizing Port Cap on the top of the Nebulizing Elbow and attach the NIVO nebulizer. When the nebulizer is not attached, close the Nebulizing Port Cap.

 **Note:** Refer to the NIVO nebulizer *Instructions for Use*.


 Figure 4

1 CapStrap Hook


 Figure 5



 Figure 6

Specifications

⚠ Warning: The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

Pressure Flow Curve

Resistance with Anti-Asphyxia Valve Closed to Atmosphere
Drop in Pressure at

	50 SLPM	100 SLPM
All sizes	0.3 cm H ₂ O	0.8 cm H ₂ O

Closed to Atmosphere: A minimum of 3 cm H₂O must be maintained to ensure that the anti-asphyxia valve closes.

Open to Atmosphere: In the absence of positive airway pressure (0 cm H₂O), the anti-asphyxia valve will open.

Inspiratory and Expiratory Resistance with Anti-Asphyxia Valve Open to Atmosphere at 50 SLPM

Inspiratory Resistance	1.5 cm H ₂ O
Expiratory Resistance	1.0 cm H ₂ O

Deadspace

Small	304 mL
Medium	422 mL
Large	466 mL

Sound Levels at 10 cm H₂O

A-weighted Sound Power Level	17.3 dBA
A-weighted Sound Pressure Level at 1m	10.2 dBA

Disposal

Dispose of in accordance with local regulations.

The useful life of the AF531 Oro-Nasal face mask is single use for one week. Dispose of within 7 days, or before, as needed.

Storage Conditions

Temperature: -20° to +60° C
Relative Humidity: 15% to 95%

Symbols Glossary

	Does not contain phthalates		Does not contain DEHP		Not made with natural rubber latex
	Warning/Caution		Operator's manual; operating instructions Consult instructions for use.		Note
	Do not reuse		Packaging unit To indicate the number of pieces in the package.		Authorized representative in the European Community
	Prescription device.		Reorder number Indicates the manufacturer's catalogue number so the medical device can be identified.		Mask Leak Symbol and Value
	Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified.		Temperature Limits		Unique Device Identifier Indicates the Unique Device Identifier information.
	Humidity Limits		Importer - Indicates the entity importing the medical device into the locale.		
	Manufacturer - Indicates the medical device manufacturer		Date of manufacture - Indicates the date when the medical device was manufactured. Country of manufacture - Indicates the country of manufacture of the product. Note: When applied to the label, "CC" is replaced by the country code.		

Please refer to <http://symbols.philips.com> for a description of the symbols used on this device and its packaging.

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