REF 45FK12

Specimen collection and handling

1. Specimen collection

Throat swab specimens should be collected using standard clinical methods. To collect a throat swab specimen, hold down the tongue with a depressor and rub the swab on the tonsils or any areas of inflammation in the back of the throat with pus or redness. Be careful not to touch the tongue or sides of the mouth with the swab. Use sterile swabs on solid plastic shafts to collect throat specimen; do not use swabs with calcium alginate or cotton tip or with wooden shafts. The swab's extraction reagents will kill the bacteria, making it impossible to culture a specimen on the swab. If bacterial culture is required, gently streak the swab on a 5 % sheep blood agar plate before testing. Alternatively, a dual swab procedure or a subsequent second swab specimen may be collected for a culture.

2. Specimen transport and storage

It is recommended that swab specimens be processed as soon as possible after collection.

Swabs can be held in any clean, dry plastic tube or sleeve for up to 4 hours at 15 – 30 °C, or for up to 24 hours if refrigerated (2 – 8 °C) before processing. The use of charcoal or agar medium is not recommended.

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About the test

[Introduction]

Streptococcus is a group of gram-positive bacteria, some species of which cause disease. Streptococci are spherical and divide by fission but remain attached and thus grow in beadlike chains. The incidence and severity of streptococcal diseases decreased dramatically after the introduction of antibiotics (penicillin, erythromycin, and select cephalosporins are effective against the organisms), but in the late 1980s the medical community was shaken by the emergence of several severe forms of streptococcal infection and drug-resistant strains. Group A Streptococcus (GAS) causes numerous infections, chiefly of the respiratory tract, but can also affect other organ systems such as the skin and lungs. The throat infections are rarely serious and, as with most respiratory infections, are typically easily handled by the immune system. Sometimes, however, the immune system forms antibodies that can cross react with the blood vessels in various organs of the body, a condition known as vasculitis. The vasculitis characteristic of a streptococcal infection is scarlet fever or scarlatina.

[Test principle]

The Bioline™ Strep A Test strip has 2 pre-coated lines on the surface of the strip: "T" (test line) and "C" (control line). Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and should always appear if the test procedure is performed properly and the control reagents are working. Both the test line and the control line will be visible in the result window if enough GAS antigens are present in the specimen. If GAS antigens are not present or are present at very low levels in the specimen, no color will appear at the test line.

[Intended use]

The Bioline™ Strep A Test is a chromatographic immunoassay for the qualitative detection of GAS antigens in throat swab specimens or confirmation of presumptive GAS colonies recovered from a culture. The assay detects both viable and nonviable organisms in throat swabs and culture colonies within 5-10 minutes. The Bioline™ Strep A test strips are immobilized with goat anti-Streptococcus Group A and specially selected antibodies are used as detector materials. This test is intended for professional use.

ENGLISH

Materials provided and active ingredients of main components

1. The Bioline™ Strep A Test kit contains the following items to perform the assay:

• 25 Test strips with desiccant in individual foil pouches

• Extraction reagent A (1 x 5 ml/vial) 

• Extraction reagent B (1 x 5 ml/vial)

• Positive control (1 x 1 ml/vial)

• Negative control (1 x 1 ml/vial)

• 25 Sterile throat swabs

• 25 Disposable test tubes

• 1 Instructions for use

• 2 Active ingredients of main components

• 1 test strip includes:

- Gold Conjugates: Rabbit anti-Strep A-gold colloid (0.126±0.025 µg), Mouse IgG-gold colloid (0.042±0.008 µg)

- Test line: Goat anti-Strep A (1.0±0.2 µg)

- Control line: Goat anti-mouse IgG (0.5±0.1 µg)

• Extraction reagent A contains 2 M sodium nitrite (q.s.)

• Extraction reagent B contains 0.2 M acetic acid (q.s.)

• Positive control: heat-inactivated Group A Streptococcus (1.25x10⁷ CFU/ml)

• Negative control: heat-inactivated Group B Streptococcus (1.25x10⁷ CFU/ml)

• Classification according to Regulation (EC) No 1272/2008:

ENGLISH

Materials provided and active ingredients of main components

1. The Bioline™ Strep A Test kit contains the following items to perform the assay:

• 25 Test strips with desiccant in individual foil pouches

• Extraction reagent A (1 x 5 ml/vial) 

• Extraction reagent B (1 x 5 ml/vial)

• Positive control (1 x 1 ml/vial)

• Negative control (1 x 1 ml/vial)

• 25 Sterile throat swabs

• 25 Disposable test tubes

• 1 Instructions for use

• 2 Active ingredients of main components

• 1 test strip includes:

- Gold Conjugates: Rabbit anti-Strep A-gold colloid (0.126±0.025 µg), Mouse IgG-gold colloid (0.042±0.008 µg)

- Test line: Goat anti-Strep A (1.0±0.2 µg)

- Control line: Goat anti-mouse IgG (0.5±0.1 µg)

• Extraction reagent A contains 2 M sodium nitrite (q.s.)

• Extraction reagent B contains 0.2 M acetic acid (q.s.)

• Positive control: heat-inactivated Group A Streptococcus (1.25x10⁷ CFU/ml)

• Negative control: heat-inactivated Group B Streptococcus (1.25x10⁷ CFU/ml)

• Classification according to Regulation (EC) No 1272/2008:

ENGLISH

Materials provided and active ingredients of main components

1. The Bioline™ Strep A Test kit contains the following items to perform the assay:

• 25 Test strips with desiccant in individual foil pouches

• Extraction reagent A (1 x 5 ml/vial) 

• Extraction reagent B (1 x 5 ml/vial)

• Positive control (1 x 1 ml/vial)

• Negative control (1 x 1 ml/vial)

• 25 Sterile throat swabs

• 25 Disposable test tubes

• 1 Instructions for use

• 2 Active ingredients of main components

• 1 test strip includes:

- Gold Conjugates: Rabbit anti-Strep A-gold colloid (0.126±0.025 µg), Mouse IgG-gold colloid (0.042±0.008 µg)

- Test line: Goat anti-Strep A (1.0±0.2 µg)

- Control line: Goat anti-mouse IgG (0.5±0.1 µg)

• Extraction reagent A contains 2 M sodium nitrite (q.s.)

• Extraction reagent B contains 0.2 M acetic acid (q.s.)

• Positive control: heat-inactivated Group A Streptococcus (1.25x10⁷ CFU/ml)

• Negative control: heat-inactivated Group B Streptococcus (1.25x10⁷ CFU/ml)

• Classification according to Regulation (EC) No 1272/2008:

ENGLISH

Materials provided and active ingredients of main components

1. The Bioline™ Strep A Test kit contains the following items to perform the assay:

• 25 Test strips with desiccant in individual foil pouches

• Extraction reagent A (1 x 5 ml/vial) 

• Extraction reagent B (1 x 5 ml/vial)

• Positive control (1 x 1 ml/vial)

• Negative control (1 x 1 ml/vial)

• 25 Sterile throat swabs

• 25 Disposable test tubes

• 1 Instructions for use

• 2 Active ingredients of main components

• 1 test strip includes:

- Gold Conjugates: Rabbit anti-Strep A-gold colloid (0.126±0.025 µg), Mouse IgG-gold colloid (0.042±0.008 µg)

- Test line: Goat anti-Strep A (1.0±0.2 µg)

- Control line: Goat anti-mouse IgG (0.5±0.1 µg)

• Extraction reagent A contains 2 M sodium nitrite (q.s.)

• Extraction reagent B contains 0.2 M acetic acid (q.s.)

• Positive control: heat-inactivated Group A Streptococcus (1.25x10⁷ CFU/ml)

• Negative control: heat-inactivated Group B Streptococcus (1.25x10⁷ CFU/ml)

• Classification according to Regulation (EC) No 1272/2008:

ENGLISH

Materials provided and active ingredients of main components

1. The Bioline™ Strep A Test kit contains the following items to perform the assay:

• 25 Test strips with desiccant in individual foil pouches

• Extraction reagent A (1 x 5 ml/vial) 

• Extraction reagent B (1 x 5 ml/vial)

• Positive control (1 x 1 ml/vial)

• Negative control (1 x 1 ml/vial)

• 25 Sterile throat swabs

• 25 Disposable test tubes

• 1 Instructions for use

• 2 Active ingredients of main components

• 1 test strip includes:

- Gold Conjugates: Rabbit anti-Strep A-gold colloid (0.126±0.025 µg), Mouse IgG-gold colloid (0.042±0.008 µg)

- Test line: Goat anti-Strep A (1.0±0.2 µg)

- Control line: Goat anti-mouse IgG (0.5±0.1 µg)

• Extraction reagent A contains 2 M sodium nitrite (q.s.)

• Extraction reagent B contains 0.2 M acetic acid (q.s.)

• Positive control: heat-inactivated Group A Streptococcus (1.25x10⁷ CFU/ml)

• Negative control: heat-inactivated Group B Streptococcus (1.25x10⁷ CFU/ml)

• Classification according to Regulation (EC) No 1272/2008:

Product identifier	Trade name	Bioline™ Strep A Extraction reagent A
Hazardous substance	Sodium nitrite (CAS No. 7632-00-0)	
Classification	Acute toxicity (oral) Category 4	Precautory statements
Hazard Pictogram		Preven-tion P264: Wash hands thoroughly after handling P270: Do not eat, drink or smoke when using this product Re-sponse P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell P330: Rinse mouth Signal Word Warning Disposal P501: Dispose of contents/container in accordance with local/regional/national/international regulation
Signal Word	Warning	
Hazard Statement	H302: Harmful if swallowed	

Materials required but not provided

• Protective gloves, Timer, Biohazard container.

Kit storage and stability

1. The Bioline™ Strep A Test strip should be stored at a temperature between 1 °C and 30 °C . Do not freeze the kit or its components.

2. The test strip is sensitive to both heat and humidity. Perform the test immediately after removing the test strip from the foil pouch.

3. Do not use the test kit beyond its expiration date.

4. The shelf life of the kit is as indicated on the outer package.

5. Do not use the test kit if the pouch is damaged or the seal is broken.

Warnings

1. For in vitro diagnostic use only. Do not reuse the test strip.

2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.

3. Do not eat or smoke while handling specimens.

4. All specimens should be handled as being potentially infectious.

5. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.

6. Avoid splashing or aerosol formation of specimen and extraction reagent.

7. Clean up spills thoroughly using an appropriate disinfectant.

8. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.

9. Do not mix or interchange with different specimens.

10. The components in the test kit have been quality control tested as a standard batch unit. Do not mix components with different lot numbers.

11. To avoid cross contamination, do not allow the tip of the reagent bottles to come in contact with specimen swabs.

12. Do not drink extraction reagent.

13. The extraction reagent contains a proprietary anti-microbial agent which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with extraction reagent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

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