DETERMINE™ HIV CONTROLS

Key to symbols used / Schlüssel für die verwendeten Symbole /
Clave de los símbolos utilizados / Légende des symboles utilisés /
Legenda dei simboli utilizzati / Descrição dos símbolos utilizados

2°C

8°C

Store at 2-8°C /
Bei 2-8°C lagern /
Guardar a temperaturas entre 2 y 8°C /
Conservar entre 2-8°C /
Conservare a 2-8°C /
Armazenar entre 2-8°C

HIV

CONTROL

-

HIV Negative Control/
HIV-Negativkontrolle/
Control negativo del VIH/
Solution de contrôle négatif VIH/
Controllo Negativo HIV/
Controlo negativo do VIH

HIV-1 Ab

CONTROL

+

HIV-1 Antibody Positive Control/
HIV-1-Antikörper-Positivkontrolle/
Control positivo de anticuerpos del VIH-1/
Solution de contrôle positif à l'égard des anticorps anti-VIH-1/
Controllo Positivo per gli Anticorpi HIV-1/
Controlo positivo para anticorpos do VIH-1

HIV-2 Ab

CONTROL

+

HIV-2 Antibody Positive Control/
HIV-2-Antikörper-Positivkontrolle/
Control positivo de anticuerpos del VIH-2/
Solution de contrôle positif à l'égard des anticorps anti-VIH-2/
Controllo Positivo per gli Anticorpi HIV-2/
Controlo positivo para anticorpos do VIH-2

p24 Ag

CONTROL

+

HIV-1 p24 Antigen Positive Control/
HIV-1-p24-Antigen-Positivkontrolle/
Control positivo de antígenos p24 del VIH-1/
Solution de contrôle positif à l'égard de l'antigène p24 VIH-1/
Controllo Positivo per Antigeni p24 di HIV-1/
Controlo positivo para antígenos p24 do VIH-1

BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA

1. Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
2. EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1-5; RI-R3, AI-A24.
3. CDC. Guideline for Isolation Precaution; Preventing Transmission of Infectious Agents in Healthcare Settings 2007.

Determine™ HIV Controls is for use with **Determine™ HIV Ultra** or **Determine™ HIV Early Detect**.

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
Advice Line
For further information, please contact your distributor, or call to one of the following Abbott Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+ (44) 161 483 9032	EMEpductsupport@abbott.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@abbott.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCISproductsupport@abbott.com
Latin America	+ (57) 2 661 8797	LAProductsupport@abbott.com

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EN

This instruction for use must be read carefully prior to use. Instruction for use must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from this instruction for use.

NAME AND INTENDED USE
Determine™ HIV Controls are quality control reagents for use with **Determine™ HIV Ultra** or **Determine™ HIV Early Detect** only. Perform a test using **Determine HIV Controls** under the following circumstances:

- For each new operator prior to performing tests on patient specimens
- In compliance with local laboratory guidelines
- For training purposes

It is the responsibility of each laboratory using **Determine HIV Ultra** or **Determine HIV Early Detect** to establish an adequate quality assurance program to ensure the performance of the device under its specified conditions.

SUMMARY AND EXPLANATION
Determine HIV Controls are specially formulated to ensure the test is performing correctly and to verify the user's ability to perform the assay properly and interpret the results. The Negative Control will produce a negative result. The HIV-1 Antibody Positive Control and the HIV-2 Antibody Positive Control will produce a positive antibody result. The HIV-1 p24 Antigen Positive Control will produce a positive result of antigen.

- CONTENTS**
Determine HIV Controls contain an instruction for use and four vials, HIV Negative Control, HIV-1 Antibody Positive Control, HIV-2 Antibody Positive Control, and HIV-1 p24 Antigen Positive Control as described below.
- **HIV Negative Control** [HIV CONTROL -]
One white capped vial containing 1.0 mL of processed normal human plasma negative for HIV-1 antibody, HIV-2 antibody, HIV-1 RNA, HBsAg and anti-HCV antibody. Contains ProClin™ 300 as a preservative.
 - **HIV-1 Antibody Positive Control** [HIV-1 Ab CONTROL +]
One red capped vial containing 1.0 mL of heat inactivated human serum or plasma positive for HIV-1 antibody, diluted in the HIV Negative Control.
 - **HIV-2 Antibody Positive Control** [HIV-2 Ab CONTROL +]
One yellow capped vial containing 1.0 mL of heat inactivated human serum or plasma positive for HIV-2 antibody, diluted in the HIV Negative Control.
 - **HIV-1 p24 Antigen Positive Control** [p24 Ag CONTROL +]
One blue capped vial containing 1.0 mL of recombinant HIV-1 p24 antigen diluted in the HIV Negative Control.

MATERIALS (required but not provided)
• **Determine HIV Ultra (7D2846, 7D2847) or Determine HIV Early Detect (7D2842, 7D2843)**
• Micropipette capable of delivering a 50 µl sample
• Personal protective equipment
• Timing device

WARNINGS AND PRECAUTIONS
For In Vitro diagnostic use
For professional use only
Safety data sheet available for professional user on request.

CAUTION
Contains 0.1% ProClin™ 300 (mixture of 5-Chloro-2-methyl-4-isothiazolin-3-one&2-methyl-4-isothiazolin-3-one) as preservative

Warning

- May cause an allergic skin reaction.
- Wear protective gloves/protective clothing/eye protection/face protection.
Safety glasses, disposable gloves, laboratory coat and face mask are recommended.
- If skin irritation or rash occurs: Get medical advice/attention.
- Wash contaminated clothing before reuse.

Contains materials of human origin
This product should be handled as potentially infectious material, as no known test method/inactivation procedure can offer complete assurance that products derived from human blood will not transmit infectious agents. Refer to 2000/54/EC for information on handling biohazardous materials. Avoid contact with the eyes, skin and mucous membranes. Wash hand before breaks and after work. Clean work areas with hypochlorite or other disinfecting agent.

Appropriate biosafety practices^{1,2,3} should be used when handling reagents. These precautions include, but are not limited to the following:

- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of reagents using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.^{1,2,3}
- Decontaminate and dispose of all reagents and other potentially contaminated materials in accordance with local regulations.^{1,2,3}

STORAGE

- **Determine HIV Controls** must be stored at 2 to 8 °C.
- When handled and stored as directed, kit components are stable until the expiration date.
- Do not use beyond the indicated expiration date.
- Open the vial only when performing the test. Close the cap and store the vial in the original container at 2 to 8 °C immediately after use.
- Do not expose the controls to temperatures exceeding 30 °C during testing.
- Store vials upright position to prevent leakage.
- Do not use controls that are visibly turbid.

TEST PROCEDURE
All components of **Determine HIV Controls** are ready to use as supplied.

- a. Read **Determine HIV Ultra** or **Determine HIV Early Detect** instruction for use prior to using **Determine HIV Controls**.
- b. Use four controls, HIV Negative Control, HIV-1 Antibody Positive Control, HIV-2 Antibody Positive Control and HIV-1 p24 Antigen Positive Control.
- c. Remove the protective foil cover from each test strip.
- d. Mix each control vial and open one vial at a time.
- e. Apply 50 µL of control (micropipette) to the sample pad (marked by the arrow symbol). Use a new pipette tip with each control.
- f. Close the same color cap tightly and store the vial in the original container at 2 to 8 °C.
- g. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read results.

EXPECTED RESULTS

HIV Negative Control (One Bar) [HIV CONTROL -]
One red bar appears in the control window of the strip (labeled "C"), and no red bar appears in the Ab and Ag windows of the strip (labeled "AG" and "AB").

HIV-1 Antibody Positive Control (Two Bars - Control and Ab Bars) [HIV-1 Ab CONTROL +]
Red bars appear in both the control window (labeled "C") and in the Ab window (labeled "AB") of the strip. Any visible red bar in the Ab window should be interpreted as positive.

HIV-2 Antibody Positive Control (Two Bars - Control and Ab Bars) [HIV-2 Ab CONTROL +]
Red bars appear in both the control window (labeled "C") and in the Ab window (labeled "AB") of the strip. Any visible red bar in the Ab window should be interpreted as positive.

HIV-1 p24 Antigen Positive Control (Two Bars - Control and Ag Bars) [p24 Ag CONTROL +]
Red bars appear in both the control window (labeled "C") and in the Ag window (labeled "AG") of the strip. Any visible red bar in the Ag window should be interpreted as positive.

Invalid (No Control Bar)
If there is no red bar in the control window of the strip, and even if a red bar appears in the Ab or Ag windows of the strip, the result is invalid and the test should be repeated.

NOTES:

- The test result is positive even if the Ag and/or Ab bars appear lighter or darker than the control bar.
- If an invalid result occurs repeatedly, contact your local distributor or call Technical Support.
- If HIV Negative Control, HIV-1 Antibody Positive Control, HIV-2 Antibody Positive Control or HIV-1 p24 Antigen Positive Control does not produce the expected results, please confirm the testing procedure. If the result is repeated, contact your local distributor or call Technical Support.

LIMITATIONS

- **Determine HIV Controls** are quality control reagents for use ONLY with **Determine HIV Ultra** or **Determine HIV Early Detect**. No other quality control reagent has been verified for use with **Determine HIV Ultra** or **Determine HIV Early Detect**.
- **Determine HIV Controls** have not been verified with other HIV detection kits.
- NOT for use with **Determine™ HIV-1/2**.

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