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ENGLISH

About the test

[Introduction] Noroviruses are a major cause of gastroenteritis worldwide with an estimated 23 million cases a year in the USA. They are transmitted directly from person to person and indirectly via contaminated water and food. They are frequently involved in outbreaks in communal facilities, such as nursing homes, hospitals, day nurseries and prisons, and on cruise ships.

[Test principle] The Bioline® Norovirus test has 2 pre-encoded lines on the surface of the strip: "T" (test line) and "C" (control line). Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural verification. The test always appears if the test procedure is performed properly and the control line reagents are working. Purple test line will be visible in the result window if there is enough norovirus present in the specimen. When a specimen is added to the test, noroviruses in the specimen reacts with colloidal gold-labeled norovirus specific antibodies and forms a complex of antigen-antibody colloidal gold conjugates. As this complex migrates along the length of the result window by capillary action, the complex is captured by the norovirus specific monoclonal antibodies immobilized in the test line across the result window and generates a colored line. In the absence of norovirus in specimens, a complex is not formed and no colored line appears in the result window of test device.

[Intended Use] The Bioline® Norovirus test is a chromatographic immunoassay for the qualitative detection of norovirus genogroup I (GI) and genogroup II (GII) in human fecal samples.

Materials provided and active ingredients of main components

- 1. All test kit components contain following items to perform the assay:
 - 20 Test devices with desiccant in individual foil pouch
 - 20 Specimen collection tubes
 - 20 Dropping caps
 - Assay diluent (1 x 20 mL/vial)
 - 20 Specimen collection swabs
 - 20 Disposable droppers
 - 1 Instructions for use
- 2. Active ingredients of main components:
 - Testant: anti-norovirus IgG conjugated + Mouse monoclonal anti-norovirus GI gold colloid (0.067 ± 0.013 µg), Mouse monoclonal anti-norovirus GII gold colloid (0.067 ± 0.013 µg), Test line: Mouse monoclonal anti-norovirus GI (0.480 ± 0.096 µg), Mouse monoclonal anti-norovirus GII (0.480 ± 0.096 µg), Control line: Goat anti-mouse IgG (0.480 ± 0.096 µg), Tris buffer (200 mM), Triton X-100 (2 %), Sodium azide (0.01 %)

Materials required but not provided

- Protective gloves, Tiner, Biohazard container

Kit storage and stability

- The test device should be stored at a temperature between 1°C and 30 °C. Do not freeze the kit or its components.
- Assay diluent is stable until the expiration date, if kept between 1°C and 30 °C.
- The test device is sensitive to both heat and humidity. Perform the test immediately after opening the kit. The test results may become unstable and have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test kit.

Internal quality control

The Bioline® Norovirus test device has a line "T" and a control line "C" on the surface of the strip. Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but it does not guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Test limitations

- A negative result does not exclude the possibility of norovirus infection in the patient. Failure to detect norovirus may be a result of factors such as collection of a specimen at a time when too few viruses are present or improper specimen sampling or handling.
- Test results should be interpreted in conjunction with information available from epidemiological studies. Clinical assessment of the patient and other diagnostic procedures.
- Do not mix or interchange different specimens.
- Care should be taken to avoid contamination of the end of the bottle when dropping assay diluent into a specimen well.
- Do not drink assay diluent.
- Bioline® Norovirus assay diluent contains a proprietary antimicrobial agent that prevents no hazard to the user if normal laboratory safety precautions are followed.

- Specimen collection and handling (Refer extraction)
- Sampling should take place as soon as possible after symptoms (diarrhea and vomiting) appear since virus elimination is its maximum 1 - 3 days after the onset of the symptoms.
- Fecal specimen must be collected in clean containers without any additives and fecal specimen may be stored refrigerated (2 - 8 °C) for 72 hours.

Performance characteristics				
The performance of the Bioline® Norovirus test was confirmed by realtime RT-PCR method. The results are summarized in the following tables.				
*Norovirus Agt. sensitivity and specificity				
*The testing sites: Soonchunhyang University Cheonan Hospital				
Realtime RT-PCR	Positive	Negative	Total	
Positive	58	11	69	
Negative	11	269	280	
Sensitivity	84.1%(58/69)			
Specificity	96.1%(269/280)			
PPV	86.3% (69/80)			
NPV	96.2% (280/291)			

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- 1

EN Open the package and look for the following:

 1. Test device with desiccant in individual foil pouch
 2. Assay diluent
 3. Disposable dropper
 4. Specimen collection swab
 5. Specimen collection tube
 6. Dropping cap
 7. Instructions for use

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