SURESTEP™
URINE TEST
DRUG SCREEN CASSETTE

SURESTEP™ URINE TEST DRUG SCREEN CASSETTE

Package Insert

English

Package Insert for testing any combination of the following drugs in the table below:

Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Cotinine, Fentanyl, Ketamine, Marijuana, Methadone, Methadone metabolite, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Opiate, Oxycodone, Phencyclidine, Propoxyphene, Tramadol and Tricyclic Antidepressants.

A rapid, one step screening test for the simultaneous, qualitative detection of drugs and drug metabolites in human urine.

For medical and other professional in vitro diagnostic use only.

INTENDED USE & SUMMARY

Urine based tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method to screen urine for drugs of abuse.

The SureStep™ Urine Test Drug Screen Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of any combination of the following drugs and drug metabolites without the need of instruments:¹

Test	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP 300)	d-Amphetamine	300
Amphetamine (AMP 500)	d-Amphetamine	500
Amphetamine (AMP 1000)	d-Amphetamine	1,000
Barbiturates (BAR 300)	Secobarbital	300
Benzodiazepines (BZO 200)	Oxazepam	200
Benzodiazepines (BZO 300)	Oxazepam	300
Buprenorphine (BUP 10)	Buprenorphine	10
Cocaine (COC 150)	Benzoylgonine	150
Cocaine (COC 300)	Benzoylgonine	300
Cotinine (COT 100)	Cotinine	100
Fentanyl (FTY 20)	Norfentanyl	20
Ketamine (KET 1000)	Ketamine	1,000
Marijuana (THC 50)	11-nor- Δ^9 -THC-9 COOH	50
Methadone (MTD 300)	Methadone	300
Methadone metabolite (EDDP 100)	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	100
Methamphetamine (MET 500)	d-Methamphetamine	500
Methamphetamine (MET 1000)	d-Methamphetamine	1,000

Methylenedioxyamphetamine (MDMA 500)	d,l-Methylenedioxyamphetamine	500
Morphine (MOP 300)	Morphine	300
Opiate (OPI 2000)	Morphine	2,000
Oxycodone (OXY 100)	Oxycodone	100
Phencyclidine (PCP 25)	Phencyclidine	25
Propoxyphene (PPX 300)	Propoxyphene	300
Tramadol (TRA 100)	Tramadol	100
Tricyclic Antidepressants (TCA 1000)	Nortriptyline	1,000

This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

PRINCIPLE

The SureStep™ Urine Test Drug Screen Cassette is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test line region. The presence of drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test line region.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a colored line will always appear at the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

Each test contains specific drug antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in the control line.

PRECAUTIONS

- For medical and other professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test device should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test device should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed well before testing.

MATERIALS

Materials Provided

- Test Devices
- Droppers
- Package Insert


Materials Required But Not Provided

- Specimen Collection Container
- Timer

DIRECTIONS FOR USE

Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it as soon as possible.
2. Place the test device on a clean and level surface. Hold the dropper vertically and **transfer 3 full drops of urine** (approx. 100 µL) to the specimen well (S) of the test device, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
3. Wait for the colored line(s) to appear. **Read results at 5 minutes.** Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE: **Two lines appear.** One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T). This negative result indicates that the drug concentration is below the detectable level.

***NOTE:** The shade of color in the test line region (T) will vary, but it should always be considered as negative whenever there is even a faint colored line.

POSITIVE: **One colored line appears in the control line region (C).** No line appears in the test line region (T). This positive result indicates that the drug concentration exceeds the detectable level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The SureStep™ Urine Test Drug Screen Cassette provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method.^{2,3}
2. There is a possibility that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result does not indicate level or intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. The test does not distinguish between drugs of abuse and certain medications.
7. A positive result may be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the SureStep™ Urine Test Drug Screen Cassette and a commercially available drug rapid test. Testing was performed on approximately 300 specimens previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC-MS. Negative urine specimens were screened initially by Predicate test, 10% negative specimens were confirmed by GC-MS. The following results were tabulated:

% Agreement with Commercial Kit

Specimen	AMP 300	AMP 500	AMP 1000	BAR 300	BZO 200	BZO 300	BUP 10**	COC 150
Positive	>99%	*	96%	>99%	*	90%	84%	>99%
Negative	>99%	*	>99%	99%	*	97%	>99%	>99%
Total	>99%	*	98%	99%	*	94%	95%	>99%

Specimen	COC 300	COT 100	FTY 20	KET 1000	THC 50	MTD 300	EDDP 100	MET 500
Positive	95%	>99%	*	*	>99%	>99%	*	>99%
Negative	>99%	>99%	*	*	>99%	>99%	*	82%
Total	98%	>99%	*	*	>99%	>99%	*	89%

Specimen	MET 1000	MDMA 500	MOP 300	OPI 2000	OXY 100	PCP 25	PPX 300	TRA 100	TCA 1000
Positive	99%	>99%	>99%	99%	96%	97%	>99%	*	95%
Negative	>99%	99%	>99%	>99%	99%	>99%	>99%	*	>99%
Total	>99%	99%	>99%	>99%	98%	99%	>99%	*	99%

* NOTE: Commercial kit unavailable for comparison testing.

** NOTE: BUP was compared to the self-reported use of Buprenorphine.

% Agreement with GC-MS

Specimen	AMP 300	AMP 500	AMP 1000	BAR 300	BZO 200	BZO 300	BUP 10*	COC 150
Positive	>99%	97%	96%	92%	98%	96%	98%	99%
Negative	99%	99%	95%	98%	99%	96%	>99%	99%
Total	99%	98%	95%	95%	99%	96%	>99%	99%

Specimen	COC 300	COT 100*	FTY 20*	KET 1000	THC 50	MTD 300	EDDP 100	MET 500
Positive	96%	>99%	99%	>99%	97%	99%	>99%	99%
Negative	90%	>99%	89%	97%	96%	94%	>99%	98%
Total	93%	>99%	93%	97%	97%	96%	>99%	98%

Specimen	MET 1000	MDMA 500	MOP 300	OPI 2000	OXY 100	PCP 25	PPX 300	TRA 100*	TCA 1000**
Positive	99%	>99%	>99%	98%	98%	>99%	94%	96%	>99%
Negative	93%	98%	94%	97%	99%	97%	99%	97%	89%
Total	96%	99%	97%	98%	99%	98%	97%	97%	91%

* NOTE: BUP, COT, FTY and TRA were based on LC-MS data instead of GC-MS.

** NOTE: TCA was based on HPLC data instead of GC-MS.

Analytical Sensitivity

A drug-free urine pool was spiked with drugs to the concentrations at $\pm 50\%$ cut-off and $\pm 25\%$ cut-off. The results are summarized below.

Drug Conc. (Cut-off range)	AMP 300		AMP 500		AMP 1000		BAR 300		BZO 200		BZO 300		BUP 10	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	90	0	30	0	30	0	60	0	30	0	90	0
-50% Cut-off	30	0	90	0	30	0	30	0	60	0	30	0	90	0
-25% Cut-off	25	5	88	2	23	7	20	10	60	0	26	4	78	12
Cut-off	16	14	45	45	9	21	13	17	22	38	12	18	48	42
+25% Cut-off	4	26	1	89	1	29	8	22	2	58	3	27	24	66
+50% Cut-off	0	30	0	90	0	30	0	30	0	60	0	30	0	90

Drug Conc. (Cut-off range)	COC 150		COC 300		COT 100		FTY 20		KET 1000		THC 50	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	90	0	90	0	90	0	30	0
-50% Cut-off	30	0	30	0	90	0	90	0	90	0	30	0
-25% Cut-off	27	3	30	0	90	0	79	11	48	42	30	0
Cut-off	13	17	9	21	49	41	36	54	6	84	21	9
+25% Cut-off	7	23	7	23	4	86	7	83	0	90	17	13
+50% Cut-off	0	30	0	30	0	90	0	90	0	90	0	30

Drug Conc. (Cut-off range)	MTD 300		EDDP 100		MET 500		MET 1000		MDMA 500		MOP 300	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	90	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	0	90	0	30	0	30	0	30	0	30	0
-25% Cut-off	26	4	80	10	27	3	24	6	23	7	28	2
Cut-off	16	14	51	39	13	17	18	12	15	15	20	10
+25% Cut-off	4	26	3	87	7	23	1	29	6	24	3	27
+50% Cut-off	0	30	0	90	0	30	0	30	0	30	0	30

Drug Conc. (Cut-off range)	OPI 2000		OXY 100		PCP 25		PPX 300		TRA 100		TCA 1000	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	90	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	90	0	30	0
-25% Cut-off	24	6	30	0	26	4	26	4	90	0	26	4
Cut-off	10	20	21	9	11	19	19	11	58	32	14	16
+25% Cut-off	4	26	6	24	8	22	8	22	22	68	4	26
+50% Cut-off	0	30	0	30	0	30	0	30	2	88	0	30

Analytical Specificity

The following table lists the concentration of compounds (ng/mL) that are detected positive in urine by the SureStep™ Urine Test Drug Screen Cassette at 5 minutes.

Column 1	
AMPHETAMINE 300 (AMP 300)	
d-Amphetamine	300
d,l-Amphetamine	390
l-Amphetamine	50,000
p-Hydroxyamphetamine	1,560
p-Hydroxynorephedrine	100,000
3,4-Methylenedioxyamphetamine (MDA)	1,560
β-Phenylethylamine	100,000
Phenylpropanolamine (d,l-Norephedrine)	100,000
Tyramine	100,000
AMPHETAMINE 500 (AMP 500)	
d-Amphetamine	500
d,l-Amphetamine	1,500
Methamphetamine	780
p-Hydroxybuprenorphine	1,562
l-Amphetamine	25,000
AMPHETAMINE 1,000 (AMP 1000)	
d-Amphetamine	1,000
d,l-Amphetamine	3,000
l-Amphetamine	50,000
d,l-3,4-Methylenedioxyamphetamine (MDA)	2,000
Phentermine	3,000
BARBITURATES 300 (BAR 300)	
Secobarbital	300
Alphenal	150
Amobarbital	300
Aprobarbital	200
Butabarbital	75
Butalbital	2,500
Butethal	100
Cyclopentobarbital	600
Pentobarbital	300
Phenobarbital	100
BENZODIAZEPINES 200 (BZO 200)	
Oxazepam	200
Alprazolam	30

Column 2 (Continued from Column 1)	
7-Aminoclonazepam	4,000
7-Aminoflunitrazepam	390
7-Aminonitrazepam	625
Bromazepam	390
Chlordiazepoxide	300
Clobazam	48
Clorazepate	97
Desalkylflurazepam	1,560
Diazepam	97
Estazolam	125
Flunitrazepam	25,000
α-Hydroxylalprazolam	30
d-Lorazepam	3,125
Midazolam	195
Nitrazepam	3,125
Norchlordiazepoxide	780
Nordiazepam	780
Temazepam	33
Triazolam	150
BENZODIAZEPINES 300 (BZO 300)	
Oxazepam	300
Alprazolam	196
Bromazepam	1,562
Chlordiazepoxide	1,562
Clobazam	98
Clonazepam	781
Clorazepate	195
Delorazepam	1,562
Desalkylflurazepam	390
Diazepam	195
Estazolam	2,500
Flunitrazepam	390
α-Hydroxylalprazolam	1,262
d,l-Lorazepam	1,562
RS-Lorazepam glucuronide	156

Column 3 (Continued from Column 2)

Midazolam	12,500
Nitrazepam	98
Norchlordiazepoxide	195
Nordiazepam	390
Temazepam	98
Triazolam	2,500
BUPRENORPHINE 10 (BUP 10)	
Buprenorphine	10
Buprenorphine 3-D-glucuronide	15
Norbuprenorphine	20
Norbuprenorphine 3-D-glucuronide	200
COCAINE 150 (COC 150)	
Benzoyllecgonine	150
Cocaine	400
Cocaethylene	6,250
Ecgonine	12,500
Ecgonine methylester	50,000
COCAINE 300 (COC 300)	
Benzoyllecgonine	300
Cocaine	780
Cocaethylene	12,500
Ecgonine	32,000
COTININE 100 (COT 100)	
l-Cotinine	100
S-l-Nicotine	12,500
FENTANYL 20 (FTY 20)	
Norfentanyl	20
Alfentanyl	562,500
Buspirone	12,500
Fenfluramine	37,500
Fentanyl	100
Sufentanyl	57,500
Risperidone	1,000
9-hydroxyrisperidone	1,000
KETAMINE 1,000 (KET 1000)	
Ketamine	1,000
Pentobarbital	50,000
Secobarbital	100,000

Column 4 (Continued from Column 3)

Norketamine	50,000
METHADONE 300 (MTD 300)	
Methadone	300
Doxylamine	50,000
MARIJUANA 50 (THC 50)	
11-nor- Δ^9 -THC-9 COOH	50
Cannabinol	20,000
11-nor- Δ^8 -THC-9 COOH	30
Δ^8 -THC	15,000
Δ^9 -THC	15,000
METHADONE METABOLITE 100 (EDDP 100)	
2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	100
METHAMPHETAMINE 500 (MET 500)	
d-Methamphetamine	500
d,l-Amphetamine	75,000
d-Amphetamine	50,000
Chloroquine	12,500
(1R,2S)-l-Ephedrine	50,000
p-Hydroxymethamphetamine	15,000
Mephentermine	25,000
l-Methamphetamine	4,000
3,4-Methylenedioxyamphetamine (MDMA)	1,000
l-Phenylephrine	100,000
β -Phenylethylamine	75,000
METHAMPHETAMINE 1,000 (MET 1000)	
d-Methamphetamine	1,000
p-Hydroxymethamphetamine	30,000
Mephentermine	50,000
l-Methamphetamine	8,000
d,l-3,4-Methylenedioxyamphetamine (MDMA)	2,000
METHYLENEDIOXYMETHAMPHETAMINE 500 (MDMA 500)	
d,l-3,4-Methylenedioxyamphetamine (MDMA)	500
d,l-3,4-Methylenedioxyamphetamine (MDA)	3,000

Column 5 (Continued from Column 4)

3,4-Methylenedioxyethylamphetamine (MDEA)	300
MORPHINE 300 (MOP 300)	
Morphine	300
Codeine	300
Ethylmorphine	6,250
Hydrocodone	50,000
Hydromorphone	3,125
Levorphanol	1,500
6-Monoacetylmorphine (6-MAM)	400
Morphine 3-β-D-glucuronide	1,000
Norcodeine	6,250
Normorphine	100,000
Oxycodone	30,000
Oxymorphone	100,000
Procaine	150,000
Thebaine	6,250
OPIATE 2,000 (OPI 2000)	
Morphine	2,000
Codeine	2,000
Ethylmorphine	5,000
Hydrocodone	12,500
Hydromorphone	5,000
Levorphanol	75,000
6-Monoacetylmorphine (6-MAM)	5,000
Morphine 3-β-D-glucuronide	2,000
Norcodeine	12,500
Normorphine	50,000
Oxycodone	25,000
Oxymorphone	25,000
Procaine	150,000
Thebaine	100,000
OXYCODONE 100 (OXY 100)	
Oxycodone	100

Column 6 (Continued from Column 5)

Hydrocodone	6,250
Hydromorphone	50,000
Levorphanol	50,000
Naloxone	37,500
Naltrexone	37,500
Oxymorphone	200
PHENCYCLIDINE 25 (PCP 25)	
Phencyclidine	25
4-Hydroxyphencyclidine	12,500
PROPOXYPHENE 300 (PPX 300)	
d-Propoxyphene	300
d-Norpropoxyphene	300
TRAMADOL 100 (TRA 100)	
Cis-tramadol	100
n-Desmethyl-cis-tramadol	195
o-Desmethyl-cis-tramadol	6,250
Phencyclidine	100,000
Procyclidine	100,000
d,l-O-Desmethyl venlafaxine	25,000
TRICYCLIC ANTIDEPRESSANTS 1,000 (TCA 1000)	
Nortriptyline	1,000
Amitriptyline	1,500
Clomipramine	12,500
Desipramine	200
Doxepin	2,000
Imipramine	400
Maprotiline	2,000
Nordoxepin	1,000
Promazine	1,500
Promethazine	25,000
Trimipramine	3,000

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Cotinine, Fentanyl, Ketamine, Marijuana, Methadone, Methadone metabolite, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Opiate, Oxycodone, Phencyclidine, Propoxyphene, Tramadol and Tricyclic Antidepressants positive urine. The following compounds show no cross-reactivity when tested with the SureStep[™] Urine Test Drug Screen Cassette at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

4-Acetamidophenol	Diclofenac	Ketoprofen	Phenothiazine
Acetone	Dicyclomine	Labetalol	Prednisolone
Acetophenetidin	Diflunisal	Lidocaine	Prednisone
Acetylsalicylic acid	Digoxin	Lindane	d,l-Propanolol
Albumin	4-Dimethylaminoantipyrine	Lithium	Quinacrine
alpha-Naphthaleneacetic Acid	Diphenhydramine	Loperamide	Quinidine
Aminopyrine	5,5-Diphenylhydantoin	l-Thyroxine	Quinine
Amoxapine	EMDP	Meperidine	R(-) Deprenyl
Amoxicillin	Erythromycin	Meprobamate	Riboflavin
Ampicillin	β-Estradiol	Methaqualone	Salicylic acid
Apomorphine	Estrone-3-sulfate	Methoxyphenamine	Serotonin
Ascorbic acid	Ethyl alcohol	Methylphenidate	Seroquel
Aspartame	Ethyl-p-aminobenzoate	Metoprolol	Sertraline
Atropine	Etodolac	N-Acetylprocainamide	Sodium Chloride
Benzilic acid	Famprofazone	Nalidixic acid	Sulfamethazine
Benzoic acid	Fenoprofen	Naproxen	Sulindac
Benzydamine	Fluoxetine	Niacinamide	Tetracycline
Brompheniramine	Furosemide	Nifedipine	Tetrahydrozoline
Caffeine	Gentisic acid	Nimesulide	Theophylline
Cannabidiol	d-Glucose	Norethindrone	Thiamine
Chloral Hydrate	Guaiacol Glyceryl Ether	Noscapine	Thioridazine
Chloramphenicol	Hemoglobin	d,l-Octopamine	Tolbutamide
Chlorothiazide	Hydralazine	Orphenadrine	Trans-2-phenylcyclopropylamine
Chlorpromazine	Hydrochlorothiazide	Oxalic acid	Trazodone
Chlorprothixene	Hydrocortisone	Oxolinic acid	Triamterene
Cholesterol	o-Hydroxyhippuric acid	Oxymetazoline	Trifluoperazine
Cimetidine	3-Hydroxytyramine	Papaverine	Trimethoprim
Clonidine	Ibuprofen	Pemoline	d,l-Tryptophan
Cortisone	Iproniazid	Penicillin	d,l-Tyrosine
Creatinine	Isoproterenol	Pentazocine	Uric acid
Deoxycorticosterone	Isoxsuprine	Phenelzine	Verapamil
Dextromethorphan	Kanamycin	Pheniramine	Zomepirac

	Consult instructions for use
	<i>In vitro</i> diagnostic medical device
	Temperature limit
	Authorized representative in the European Community
	Unique Device Identifier

Index of Symbols

	Contains sufficient for <n> tests
	Use by Date
	Batch code
	Do Not Get Wet

	Manufacturer
	Do not re-use
	Catalogue Number
	Keep Out of Sunlight

ASSISTANCE

For product inquiries, please contact the following via phone or email based on your region:

Europe & Middle East

Phone: +44 161 483 9032

Email: EMEprouductsupport@abbott.com

Asia Pacific

Phone: +61 7 3363 7711

Email: APprouductsupport@abbott.com

SURESTEP™ URINE TEST DRUG SCREEN CASSETTE

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Europa & Médio Oriente

Telefone: +44 161 483 9032

E-mail: EMProductsupport@abbott.com

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Telefone: +61 7 3363 7711

E-mail: APProductsupport@abbott.com

BIBLIOGRAPHY/ LITERATUR/ BIBLIOGRAFÍA/ BIBLIOGRAPHIE/ BIBLIOGRAFIA/ BIBLIOGRAFIA

1. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man, 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
3. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

Innovacon, Inc.
9942 Mesa Rim Rd.
San Diego, CA 92121, USA
www.toxicology.abbott

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© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Revision date: 2022-09-29

Number: 1156261601

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MDSS GmbH