SURESTEP™
URINE TEST
DRUG SCREEN STRIP

SURESTEP™ URINE TEST DRUG SCREEN STRIP

Package Insert

English

Package Insert for testing any combination of the following drugs in the table below:

Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Fentanyl, Heroin (6-acetylmorphine), Ketamine, Marijuana, Methadone, Methadone metabolite, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Opiate, Oxycodone, Phencyclidine, Propoxyphene, Tramadol and Tricyclic Antidepressants.

A rapid, one step screen test for the qualitative detection of drugs and metabolites in human urine.

For medical and other professional in vitro diagnostic use only.

INTENDED USE & SUMMARY

Urine based tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method to screen urine for drugs of abuse.

The SureStep™ Urine Test Drug Screen Strip is a lateral flow chromatographic immunoassay for the qualitative detection of any combination of the following drugs and drug metabolites without the need of instruments:¹

Test	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP 300)	d-Amphetamine	300
Amphetamine (AMP 500)	d-Amphetamine	500
Amphetamine (AMP 1000)	d-Amphetamine	1,000
Barbiturates (BAR 300)	Secobarbital	300
Benzodiazepines (BZO 200)	Oxazepam	200
Benzodiazepines (BZO 300)	Oxazepam	300
Buprenorphine (BUP 10)	Buprenorphine	10
Cocaine (COC 150)	Benzoyllecgonine	150
Cocaine (COC 300)	Benzoyllecgonine	300
Fentanyl (FTY 20)	Norfentanyl	20
Heroin (HRN 10)	6-Acetylmorphine	10
Ketamine (KET 1000)	Ketamine	1,000
Marijuana (THC 20)	11-nor- Δ^9 -THC-9 COOH	20
Marijuana (THC 50)	11-nor- Δ^9 -THC-9 COOH	50
Marijuana (THC 150)	11-nor- Δ^9 -THC-9 COOH	150
Methadone (MTD 300)	Methadone	300
Methadone metabolite (EDDP 100)	2-Ethylidene-1,5-dimethyl-3,3-dipheylpyrrolidine (EDDP)	100

Methamphetamine (MET 500)	d-Methamphetamine	500
Methamphetamine (MET 1000)	d-Methamphetamine	1,000
Methylenedioxymethamphetamine (MDMA 500)	d,l-Methylenedioxymethamphetamine	500
Morphine (MOP 300)	Morphine	300
Opiate (OPI 2000)	Morphine	2,000
Oxycodone (OXY 100)	Oxycodone	100
Phencyclidine (PCP 25)	Phencyclidine	25
Propoxyphene (PPX 300)	Propoxyphene	300
Tramadol (TRA 100)	Tramadol	100
Tricyclic Antidepressants (TCA 1000)	Nortriptyline	1,000

This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

PRINCIPLE

The SureStep™ Urine Test Drug Screen Strip is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test line region. The presence of drug above the cut-off concentration will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test line region.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a colored line will always appear at the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

Each test contains specific drug antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in the control line.

PRECAUTIONS

- For medical and other professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed well before testing.

MATERIALS

Materials Provided

- Test Strips
- Package Insert


Materials Required But Not Provided

- Specimen Collection Container
- Timer

DIRECTIONS FOR USE

Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test strip from the sealed pouch and use it as soon as possible.
2. With arrows pointing toward the urine specimen, **immerse the test strip vertically in the urine specimen for at least 10-15 seconds.** Do not pass the maximum line (MAX) on the test strip when immersing the strip. See the illustration below.
3. Place the test strip on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear. **Read results at 5 minutes.** Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE: Two lines appear. One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T). This negative result indicates that the drug concentration is below the detectable level.

***NOTE:** The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: One colored line appears in the control line region (C). No line appears in the test line region (T). This positive result indicates that the drug concentration exceeds the detectable level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The SureStep™ Urine Test Drug Screen Strip provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method.^{2,3}
2. There is a possibility that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result does not indicate level or intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. The test does not distinguish between drugs of abuse and certain medications.
7. A positive result may be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the SureStep™ Urine Test Drug Screen Strip and a commercially available drug rapid test. Testing was performed on a minimum of 200 specimens previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC-MS. Negative urine specimens were screened initially by Predicate test, 10% negative specimens were confirmed by GC-MS. The following results were tabulated:

% Agreement with Commercial Kit

Specimen	AMP 300	AMP 500	AMP 1000	BAR 300	BZO 200	BZO 300	BUP 10**	COC 150	COC 300
Positive	>99%	*	97%	>99%	*	90%	84%	>99%	95%
Negative	>99%	*	>99%	99%	*	97%	>99%	>99%	>99%
Total	>99%	*	98%	99%	*	94%	95%	>99%	98%

Specimen	FTY 20	HRN 10	KET 1000	THC 20	THC 50	THC 150	MTD 300	EDDP 100	MET 500
Positive	*	*	*	*	98%	*	>99%	*	>99%
Negative	*	*	*	*	>99%	*	>99%	*	80%
Total	*	*	*	*	99%	*	>99%	*	87%

Specimen	MET 1000	MDMA 500	MOP 300	OPI 2000	OXY 100	PCP 25	PPX 300	TRA 100	TCA 1000
Positive	98%	>99%	>99%	99%	96%	98%	>99%	*	95%
Negative	>99%	99%	>99%	>99%	99%	>99%	>99%	*	>99%
Total	99%	99%	>99%	>99%	98%	>99%	>99%	*	99%

* NOTE: Commercial kit unavailable for comparison testing.

** NOTE: BUP was compared to the self-reported use of Buprenorphine.

% Agreement with GC-MS

Specimen	AMP 300	AMP 500	AMP 1000	BAR 300	BZO 200	BZO 300	BUP 10*	COC 150	COC 300
Positive	>99%	95%	97%	92%	98%	97%	98%	99%	96%
Negative	99%	>99%	95%	98%	99%	95%	>99%	>99%	90%
Total	99%	98%	96%	95%	99%	96%	>99%	99%	93%

Specimen	FTY 20*	HRN 10	KET 1000	THC 20	THC 50	THC 150	MTD 300	EDDP 100
Positive	99%	99%	>99%	87%	96%	91%	99%	98%
Negative	90%	>99%	95%	99%	97%	96%	94%	>99%
Total	93%	>99%	95%	95%	96%	96%	96%	99%

Specimen	MET 500	MET 1000	MDMA 500	MOP 300	OPI 2000	OXY 100	PCP 25	PPX 300	TRA 100*	TCA 1000**
Positive	>99%	99%	>99%	>99%	98%	98%	>99%	94%	99%	>99%
Negative	96%	94%	98%	94%	97%	99%	96%	99%	96%	89%
Total	98%	96%	99%	97%	98%	99%	97%	97%	97%	91%

* NOTE: BUP, FTY and TRA were based on LC-MS data instead of GC-MS.

** NOTE: TCA was based on HPLC data instead of GC-MS.

Analytical Sensitivity

A drug-free urine pool was spiked with drugs to the concentrations at $\pm 50\%$ cut-off and $\pm 25\%$ cut-off. The results are summarized below.

Drug Conc. (Cut-off range)	AMP 300		AMP 500		AMP 1000		BAR 300		BZO 200		BZO 300		BUP 10	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	60	0	30	0	90	0
-50% Cut-off	30	0	30	0	30	0	30	0	60	0	30	0	90	0
-25% Cut-off	27	3	25	5	22	8	27	3	60	0	27	3	75	15
Cut-off	13	17	11	19	12	18	22	8	22	38	11	19	60	30
+25% Cut-off	4	26	5	25	2	28	8	22	2	58	5	25	31	59
+50% Cut-off	0	30	0	30	0	30	2	28	0	60	0	30	0	90

Drug Conc. (Cut-off range)	COC 150		COC 300		FTY 20		HRN 10		KET 1000		THC 20		THC 50	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	90	0	90	0	90	0	30	0	30	0
-50% Cut-off	30	0	30	0	90	0	90	0	90	0	30	0	30	0
-25% Cut-off	24	6	30	0	82	8	N/A	N/A	61	29	27	3	12	18
Cut-off	14	16	4	26	48	42	N/A	N/A	20	70	24	6	1	29
+25% Cut-off	7	23	0	30	11	79	N/A	N/A	2	88	17	13	1	29
+50% Cut-off	0	30	0	30	0	90	0	90	0	90	5	25	0	30

Drug Conc. (Cut-off range)	THC 150		MTD 300		EDDP 100		MET 500		MET 1000		MDMA 500	
	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	90	0	30	0	90	0	30	0	30	0	30	0
-50% Cut-off	90	0	29	1	90	0	30	0	30	0	30	0
-25% Cut-off	90	0	24	6	90	0	23	7	30	0	26	4
Cut-off	46	44	21	9	37	53	13	17	18	12	17	13
+25% Cut-off	5	85	2	28	8	82	8	22	1	29	4	26
+50% Cut-off	0	90	0	30	0	90	0	30	0	30	0	30

Drug Conc. (Cut-off range)	MOP 300		OPI 2000		OXY 100		PCP 25		PPX 300		TRA 100		TCA 1000	
	-	+	-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	0	30	0	30	0	30	0	30	0	90	0	30	0
-50% Cut-off	30	0	30	0	30	0	30	0	30	0	90	0	30	0
-25% Cut-off	25	5	25	5	30	0	19	11	24	6	90	0	29	1
Cut-off	17	13	15	15	18	12	16	14	17	13	61	29	18	12
+25% Cut-off	1	29	6	24	6	24	6	24	7	23	21	69	5	25
+50% Cut-off	0	30	0	30	0	30	0	30	0	30	2	88	0	30

Analytical Specificity

The following table lists the concentration of compounds (ng/mL) that are detected positive in urine by the SureStep™ Urine Test Drug Screen Strip at 5 minutes.

Column 1	
AMPHETAMINE 300 (AMP 300)	
d-Amphetamine	300
d,l-Amphetamine	390
l-Amphetamine	50,000
p-Hydroxyamphetamine	1,560
p-Hydroxynorephedrine	100,000
3,4-Methylenedioxyamphetamine (MDA)	1,560
β-Phenylethylamine	100,000
Phenylpropanolamine (d,l-Norephedrine)	100,000
Tyramine	100,000
AMPHETAMINE 500 (AMP 500)	
d-Amphetamine	500
d,l-Amphetamine	1,500
3,4-Methylenedioxyamphetamine (MDA)	800
Phentermine	1,500
β-Phenylethylamine	50,000
Tryptamine	50,000
Tyramine	25,000
AMPHETAMINE 1,000 (AMP 1000)	
d-Amphetamine	1,000
d,l-Amphetamine	3,000
l-Amphetamine	50,000
d,l-3,4-Methylenedioxyamphetamine (MDA)	2,000
Phentermine	3,000
BARBITURATES 300 (BAR 300)	
Secobarbital	300
Alphenal	150
Amobarbital	300
Aprobarbital	200
Butabarbital	75
Butalbital	2,500
Butethal	100

Column 2 (Continued from Column 1)	
Cyclopentobarbital	600
Pentobarbital	300
Phenobarbital	100
BENZODIAZEPINES 200 (BZO 200)	
Oxazepam	200
Alprazolam	30
7-Aminoclonazepam	4,000
7-Aminoflunitrazepam	390
7-Aminonitrazepam	625
Bromazepam	390
Chlordiazepoxide	300
Clobazam	48
Clorazepate	97
Desalkylflurazepam	1,560
Diazepam	97
Estazolam	125
Flunitrazepam	25,000
α-Hydroxyalprazolam	30
d-Lorazepam	3,125
Midazolam	195
Nitrazepam	3,125
Norchlordiazepoxide	780
Nordiazepam	780
Temazepam	33
Triazolam	150
BENZODIAZEPINES 300 (BZO 300)	
Oxazepam	300
Alprazolam	196
Bromazepam	1,562
Chlordiazepoxide	1,562
Clobazam	98
Clonazepam	781
Clorazepate	195

Column 3 (Continued from Column 2)

Delorazepam	1,562
Desalkylflurazepam	390
Diazepam	195
Estazolam	2,500
Flunitrazepam	390
α -Hydroxyalprazolam	1,262
d,l-Lorazepam	1,562
RS-Lorazepam glucuronide	156
Midazolam	12,500
Nitrazepam	98
Norchlordiazepoxide	195
Nordiazepam	390
Temazepam	98
Triazolam	2,500

BUPRENORPHINE 10 (BUP 10)

Buprenorphine	10
Buprenorphine 3-D-glucuronide	15
Norbuprenorphine	20
Norbuprenorphine 3-D-glucuronide	200

COCAINE 150 (COC 150)

Benzoylcegonine	150
Cocaine	400
Cocaehtylene	6,250
Ecgonine	12,500
Ecgonine methylester	50,000

COCAINE 300 (COC 300)

Benzoylcegonine	300
Cocaine	780
Cocaehtylene	12,500
Ecgonine	32,000

METHADONE METABOLITE 100 (EDDP 100)

2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	100
--	-----

FENTANYL 20 (FTY 20)

Norfentanyl	20
Alfentanyl	562,500

Column 4 (Continued from Column 3)

Bupirone	12,500
Fenfluramine	37,500
Fentanyl	100
Sufentanyl	57,500
Risperidone	1,000
9-hydroxyrisperidone	1,000

HEROIN 10 (HRN 10)

6-Acetylmorphine	10
6-Acetylcodeine	1,562
Morphine	500,000

KETAMINE 1,000 (KET 1000)

Ketamine	1,000
Pentobarbital	50,000
Secobarbital	100,000
Norketamine	50,000

METHADONE 300 (MTD 300)

Methadone	300
Doxylamine	50,000

MARIJUANA 20 (THC 20)

11-nor- Δ^9 -THC-9 COOH	20
Cannabinol	12,500
11-nor- Δ^8 -THC-9 COOH	20
Δ^8 -THC	10,000
Δ^9 -THC	12,500

MARIJUANA 50 (THC 50)

11-nor- Δ^9 -THC-9 COOH	50
Cannabinol	20,000
11-nor- Δ^8 -THC-9 COOH	30
Δ^8 -THC	15,000
Δ^9 -THC	15,000

MARIJUANA 150 (THC 150)

11-nor- Δ^9 -THC-9 COOH	150
Cannabinol	25,000
11-nor- Δ^8 -THC-9 COOH	500
Δ^8 -THC	25,000
Δ^9 -THC	25,000

Column 5 (Continued from Column 4)

METHAMPHETAMINE 500 (MET 500)	
d-Methamphetamine	500
d,l-Amphetamine	75,000
d-Amphetamine	50,000
Chloroquine	12,500
(1R,2S)-l-Ephedrine	50,000
p-Hydroxymethamphetamine	15,000
Mephentermine	25,000
l-Methamphetamine	4,000
3,4-Methylenedioxymethamphetamine (MDMA)	1,000
l-Phenylephrine	100,000
β-Phenylethylamine	75,000
METHAMPHETAMINE 1,000 (MET 1,000)	
d-Methamphetamine	1,000
p-Hydroxymethamphetamine	30,000
Mephentermine	50,000
l-Methamphetamine	8,000
d,l-3,4-Methylenedioxymethamphetamine (MDMA)	2,000
METHYLENEDIOXYMETHAMPHETAMINE 500 (MDMA 500)	
d,l-3,4-Methylenedioxymethamphetamine (MDMA)	500
d,l-3,4-Methylenedioxyamphetamine (MDA)	3,000
3,4-Methylenedioxyethylamphetamine (MDEA)	300
MORPHINE 300 (MOP 300)	
Morphine	300
Codeine	300
Ethylmorphine	6,250
Hydrocodone	50,000
Hydromorphone	3,125
Levorphanol	1,500
6-Acetylmorphine	400

Column 6 (Continued from Column 5)

Morphine 3-β-D-glucuronide	1,000
Norcodeine	6,250
Normorphine	100,000
Oxycodone	30,000
Oxymorphone	100,000
Procaine	150,000
Thebaine	6,250
OPIATE 2,000 (OPI 2000)	
Morphine	2,000
Codeine	2,000
Ethylmorphine	5,000
Hydrocodone	12,500
Hydromorphone	5,000
Levorphanol	75,000
6-Acetylmorphine	5,000
Morphine 3-β-D-glucuronide	2,000
Norcodeine	12,500
Normorphine	50,000
Oxycodone	25,000
Oxymorphone	25,000
Procaine	150,000
Thebaine	100,000
OXYCODONE 100 (OXY 100)	
Oxycodone	100
Hydrocodone	6,250
Hydromorphone	50,000
Levorphanol	50,000
Naloxone	37,500
Naltrexone	37,500
Oxymorphone	200
PHENCYCLIDINE 25 (PCP 25)	
Phencyclidine	25
4-Hydroxyphencyclidine	12,500
PROPOXYPHENE 300 (PPX 300)	
d-Propoxyphene	300
d-Norpropoxyphene	300

Column 7 (Continued from Column 6)

TRAMADOL 100 (TRA 100)	
Cis-tramadol	100
n-Desmethyl-cis-tramadol	195
o-Desmethyl-cis-tramadol	6,250
Phencyclidine	100,000
Procyclidine	100,000
d,l-O-Desmethyl venlafaxine	25,000
TRICYCLIC ANTIDEPRESSANTS 1,000 (TCA 1000)	
Nortriptyline	1,000
Amitriptyline	1,500

Column 8 (Continued from Column 7)

Clomipramine	12,500
Desipramine	200
Doxepin	2,000
Imipramine	400
Maprotiline	2,000
Nordoxepin	1,000
Promazine	1,500
Promethazine	25,000
Trimipramine	3,000

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Cotinine, Fentanyl, Heroin (6-acetylmorphine), Ketamine, Marijuana, Methadone, Methadone metabolite, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Opiate, Oxycodone, Phencyclidine, Propoxyphene, Tramadol and Tricyclic Antidepressants positive urine. The following compounds show no cross-reactivity when tested with the SureStep™ Urine Test Drug Screen Strip at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

4-Acetamidophenol	Caffeine	Digoxin	Hemoglobin
Acetone	Cannabidiol	4-Dimethylaminoantipyrine	Hydralazine
Acetophenetidin	Chloral Hydrate	Diphenhydramine	Hydrochlorothiazide
Acetylsalicylic acid	Chloramphenicol	5,5-Diphenylhydantoin	Hydrocortisone
Albumin	Chlorothiazide	EMDP	o-Hydroxyhippuric acid
alpha-Naphthaleneacetic Acid	Chlorpromazine	Erythromycin	3-Hydroxytyramine
Aminopyrine	Chlorprothixene	β-Estradiol	Ibuprofen
Amoxapine	Cholesterol	Estrone-3-sulfate	Iproniazid
Amoxicillin	Cimetidine	Ethyl alcohol	Isoproterenol
Ampicillin	Clonidine	Ethyl-p-aminobenzoate	Isosuprine
Apomorphine	Cortisone	Etodolac	Kanamycin
Ascorbic acid	(-)-Cotinine	Famprofazone	Ketoprofen
Aspartame	Creatinine	Fenoprofen	Labetalol
Atropine	Deoxycorticosterone	Fluoxetine	Lidocaine
Benzoic acid	Dextromethorphan	Furosemide	Lindane
Benzoic acid	Diclofenac	Gentisic acid	Lithium
Benzydamine	Dicyclomine	d-Glucose	Loperamide
Brompheniramine	Diflunisal	Guaiaacol Glyceryl Ether	l-Thyroxine

Meperidine	Orphenadrine	Quinacrine	Theophylline
Meprobamate	Oxalic acid	Quinidine	Thiamine
Methaqualone	Oxolinic acid	Quinine	Thioridazine
Methoxyphenamine	Oxymetazoline	RC(-) Deprenyl	Tolbutamide
Metoprolol	Papaverine	Riboflavin	Trans-2-phenylcyclopropylamine
N-Acetylprocainamide	Pemoline	Salicylic acid	Trazodone
Nalidixic acid	Penicillin	Serotonin	Triamterene
Naproxen	Pentazocine	Seroquel	Trifluoperazine
Niacinamide	Phenelzine	Sertraline	Trimethoprim
Nifedipine	Pheniramine	Sodium Chloride	d,l-Tryptophan
Nimesulide	Phenothiazine	Sulfamethazine	d,l-Tyrosine
Norethindrone	Prednisolone	Uric acid	Verapamil
Noscapine	Prednisone	Tetracycline	Zomepirac
d,l-Octopamine	d,l-Propanolol	Tetrahydrozoline	

Index of Symbols

	Consult instructions for use		Contains sufficient for <n> tests		Manufacturer
	In vitro diagnostic medical device		Use by Date		Do not re-use
	Temperature limit		Batch code		Catalogue Number
	Authorized representative in the European Community		Do Not Get Wet		Keep Out of Sunlight
	Unique Device Identifier				

ASSISTANCE

For product inquiries, please contact the following via phone or email based on your region:

Europe & Middle East

Phone: +44 161 483 9032

Email: EMEproductsupport@abbott.com

Asia Pacific

Phone: +61 7 3363 7711

Email: APproductsupport@abbott.com

SURESTEP™ URINE TEST DRUG SCREEN STRIP

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E-mail: EMEdproductsupport@abbott.com

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Telefone: +61 7 3363 7711

E-mail: APproductsupport@abbott.com

BIBLIOGRAPHY/ LITERATUR/ BIBLIOGRAFÍA/ BIBLIOGRAPHIE/ BIBLIOGRAFIA/ BIBLIOGRAFIA

1. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
3. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

Innovacon, Inc.
9942 Mesa Rim Rd.
San Diego, CA 92121, USA
www.toxicology.abbott

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© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Revision date: 2022-09-22

Number: 1156261901

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MDSS GmbH