QuantisalTM

Oral Fluid Collection Device

INSTRUCTIONS FOR USE

TABLE OF CONTENTS

<u>ENGLISH</u>	<u>4</u>
<u>FRANÇAIS</u>	<u>11</u>
<u>ITALIANO</u>	<u>19</u>
<u>DEUTSCH</u>	<u>27</u>
<u>ESPAÑOL</u>	<u>35</u>
<u>PORTUGUÊS</u>	<u>43</u>
<u>NORSK</u>	<u>51</u>
<u>SVENSKA</u>	<u>58</u>
<u>DANSK</u>	<u>66</u>
<u>NEDERLANDS</u>	<u>74</u>

REF

QSI-0025NINT

25 devices per kit

GTIN: 00840937112311

QSI-0500NINT

500 devices per kit

GTIN: 00840937112328

MD

INTENDED PURPOSE

The Quantisal™ Oral Fluid Collection Device is intended for the collection, preservation, and transport of oral fluid specimens for analytical testing of drugs or drug metabolites. This device is for use only under observed collections.

SUMMARY AND EXPLANATION OF FUNCTION

Collection: An oral fluid specimen is collected by placing a cellulose pad affixed to a polypropylene stem (Collector) under the tongue of an individual until a defined volume of oral fluid has saturated the cellulose pad. The defined volume taken up by the cellulose pads is indicated by coloration (blue) in a window on the stem (volume adequacy).

Preservation: The collector is then transferred into a provided polypropylene tube containing a specific volume of preservative buffer. The tube is stoppered with the provided cap. The specimens are ready for storage or transport.

PRINCIPLE OF OPERATION

The Quantisal™ device is designed to collect 1 mL ± 10% of oral fluid on the collection pad. The cellulose pad placed under the tongue for collection is extended into a cellulose tail enclosed in the interior of the polypropylene stem. Blue dye is placed on the cellulose tail below the window in the stem. Upon saturating the cellulose pad with specimen, the oral fluid aqueous medium migrates along the cellulose tail by capillary action and dissolves the dye. Upon further migration of liquid, the dye becomes visible in the window of the stem. Specimen collection is complete when the indicator turns completely blue.

MATERIALS PROVIDED

- 25 or 500 pouched Quantisal™ Oral Fluid Collection Devices consisting of:
 - 1 pouched Collector comprised of a cellulose pad affixed to polypropylene stem with volume adequacy indicator; and
 - 1 Transport Tube with Cap containing 3 mL of Transportation Buffer with preservative.
- 1 Instructions for Use

WARNINGS AND PRECAUTIONS

Warning

- Preservative Buffer contains a hazardous preservative: Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone.
 - H317 – May cause an allergic skin reaction.
 - H412 – Harmful to aquatic life with long lasting effects.
 - P261 – Avoid breathing vapor, mist, spray.
 - P272 – Contaminated work clothing should not be allowed out of the workplace.
 - P280 – Wear eye protection, face protection, face shield, protective clothing, protective gloves.

- P302+P352 – IF ON SKIN: Wash with plenty of soap and water.
 - P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.
 - P362+P364 – Take off contaminated clothing and wash it before reuse.
 - P501 – Dispose of contents and container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.
- Wear protective gloves classified under EN 374 or ATSM F1296 when handling the Transportation Tube when prolonged or repeated contact with the Preservative Buffer could occur.
 - Refer to Safety Data Sheet for additional information on the safe handling of the product.
 - The Preservative Buffer is not intended to contact the mouth. Do not allow the specimen donor to ingest the liquid.
 - The sample collection should be observed and instructed by the operator to ensure consistent and complete collection of oral fluid specimens.
 - Quantisal™ is a single use only device. Do not re-use
 - The device is designed to collect oral fluid specimens. Use of this product may pose a choking hazard.
 - Do not chew or suck on pad or stem, remove the stem cover, separate the collector into split pad/stem or talk before collection is completed. Misuse of this product may cause separation of the plastic stem and lead to ingestion or lodging of pad or plastic in the windpipe.
 - Samples collected with the Quantisal™ device have been diluted 1 to 4. Test results should be corrected for the dilution factor. However, Immunoanalysis enzyme immunoassays intended to be used with the Quantisal™ device are designed to factor in this dilution, and no additional calculation or correction to the test result is needed.
 - Do not use device beyond the expiration date printed on the package.

STORAGE

Store unused collection devices between 15-30°C. The devices are stable until the expiry date on the label.

SPECIMEN COLLECTION PROCEDURE

The following instructions are written for the operator to administer the collection procedure. The operator must instruct and observe the donor using the Quantisal™ device. Refer to the illustrations on the pouched Quantisal™ device when instructing the donor on proper use of the product.

A video on the collection procedure is available at <https://immunoanalysis.com/resources/training-videos/>

Avoid exposure to human oral fluid specimens. When prolonged or repeated contact with the Preservative Buffer could occur, it is advised to use gloves to prevent skin irritation.

1. Ensure donor has refrained from eating, drinking, or smoking for 10 minutes prior to specimen collection.
2. Check the expiration date of the device before use.
3. Instruct the donor to check packaging for integrity. Do not use the device if the package or device is broken or the Transport Tube leaks.
4. Record donor information and collection date and time on tube labels and/or the custody and control form or requisition form associated with the specimen.
5. Instruct the donor to remove the Collector from the pouch.

6. Instruct the donor to position the Collector under the tongue and to close mouth. The Donor should keep the head down to allow gravity to help with saliva collection.
WARNING: THE DONOR MUST NOT CHEW OR SUCK ON PADS OR STEM, REMOVE THE STEM COVER OFF THE COLLECTOR, TALK, PULL-ON OR REMOVE COLLECTOR FROM MOUTH UNTIL INDICATOR COMPLETELY TURNS BLUE.
7. The donor should keep the pad under the tongue until the indicator turns completely blue. The blue color indicates that the Collector has collected 1 mL $\pm$ 10% of oral fluid. The collection time may take from 2-10 minutes. If the indicator has not turned blue within 10 minutes, the pad should be removed from the mouth and discarded. Re-collection with a new device may begin immediately but only after oral fluid has accumulated in the mouth.
8. When the indicator completely turns blue, instruct the donor to remove the Collector from mouth.
9. Hold the Transport Tube in an upright position and remove the cap by pushing up with thumb(s).
10. Insert the saturated collector into the tube. Do not spill or empty liquid from tube.
WARNING: DO NOT INGEST LIQUID. DO NOT PLACE COLLECTOR INTO MOUTH AFTER IT HAS BEEN IN THE PRESERVATIVE BUFFER.
11. Carefully place cap over the top of the collector stem in the tube. **FORCEFULLY PUSH CAP DOWNWARD UNTIL CAP "SNAPS" FLUSH WITH TOP OF TUBE.**
12. Ship the oral fluid specimen at ambient temperature via overnight courier service to a laboratory for testing as soon as possible. Follow all lab-provided instructions and applicable regional regulations for labeling and packaging the Transport Tube for shipment. Studies have been performed to demonstrate that samples can withstand temperature excursions down to -20°C and up to 40°C, not to exceed 24 hours. Specimen must be received at the laboratory, stored, and tested within the timeframe detailed in the "Instructions for the Laboratory" section.

INSTRUCTIONS FOR THE LABORATORY

Check the custody and control form or requisition form associated with the specimen.

Check the collection date and time on the transport tube label. Samples should not be tested until a minimum of 4 hours has passed after specimen collection to allow sufficient time for any drug(s) present to be extracted from the pad into the preservation buffer.

Samples must be tested within the timeframe listed in **Table 2**.

Prior to testing, inspect the sample tube for signs of leakage. If the tube is broken or shows signs of leakage, the sample should not be tested. **Insufficient preservative buffer volume due to leaky or broken transport tube may cause erroneous results.**

1. Open the red cap of the Transport Tube. Hold Transport Tube in an upright position and uncap by pushing up with thumb(s). Do not use if the tube leaks.
2. Aliquot the sample from the transport tube for drug testing.
3. When ready to discard, safely dispose of the sample per local guidelines.

Samples collected with the Quantisal™ Oral Fluid Collection Device have been diluted 1 to 4. Immunalysis enzyme immunoassays intended to be used with the Quantisal™ Oral Fluid Collection Device are designed to factor this dilution. However, results from alternative chromatographic or mass spectral testing methods should calculate this dilution factor into test results.

PERFORMANCE CHARACTERISTICS

Device performance was established for the drugs listed in **Table 1** using mass spectrometry.

Table 1. Drugs and Concentrations

Drugs	Concentration (ng/mL)
THC	4
Benzoylecgonine	15
Cocaine	15
Morphine	30
Codeine	30
Oxycodone	30
Hydrocodone	30
6-acetylmorphine	4
Phencyclidine	10
Amphetamine	50
Methamphetamine	50
Buprenorphine	3
Methadone	20
Benzodiazepines	5
Tramadol	50

Sample Volume and Collection Time: Seventy-five oral fluid samples from known drug users were collected using the Quantisal™ device. Sample volume and collection time were recorded. All (75/75) drug users provided oral fluid volume of 1 mL ± 10%. All (75/75) samples were collected within 10 minutes.

Drug Stability: The stability of the drugs in oral fluid specimens collected with the Quantisal™ Oral Fluid Collection Device was evaluated with low positive (above and near the concentration listed in **Table 1**) samples at room temperature (8-25°C) and refrigerated (2-8°C). The stability of drugs in oral fluid specimens was measured by comparing the concentration of the specimen over time to the initial concentration. Results were within ±10% of the initial concentration and are listed in **Table 2**.

Table 2. Oral Fluid Specimens Storage in the Quantisal™ Transport Tube

Drugs	Initial Concentration (ng/mL)	Stability at 8-25°C	Stability at 2-8°C
THC	6.0	10 days	2 months
Benzoylecgonine	22	10 days	12 months
Cocaine	22	5 days	1 month
Morphine	45	10 days	12 months
Codeine	46	10 days	12 months
Oxycodone	47	10 days	12 months
Hydrocodone	46	10 days	12 months
6-acetylmorphine	5.8	10 days	12 months

Table 2. Oral Fluid Specimens Storage in the Quantisal™ Transport Tube continued

Drugs	Initial Concentration (ng/mL)	Stability at 8-25°C	Stability at 2-8°C
Phencyclidine	14	10 days	12 months
Amphetamine	75	10 days	12 months
Methamphetamine	74	10 days	12 months
Buprenorphine	4.5	10 days	12 months
Methadone	29	10 days	12 months
Benzodiazepines	7.5	10 days	12 months
Tramadol	75	10 days	12 months

Drug Recovery: Three Quantisal™ devices were introduced independently into test solutions prepared by spiking the drug into drug free negative oral fluid at $\pm 25\%$ and $+50\%$ of the concentrations listed in **Table 1**. The collectors were then placed into the Transport Tubes, sealed, and stored overnight at room temperature. The next day the samples were tested by mass spectrometry to demonstrate that the percentage recovery was greater than 80% of the initial concentration.

Clinical Specimen Study: At least sixty expectorated samples for each drug collected from self-reported drug users were directly processed through the Quantisal™ device. The samples were analyzed by mass spectrometry and the results were compared with those obtained from expectorated samples. The results showed 899/900 samples met the criteria that drug concentration following the introduction of Quantisal™ devices into the expectorated oral fluid was within $\pm 20\%$ of expectoration concentration. This study was designed for clinical study purposes only. Follow the intended use and sample collection instructions for the correct use of the device.

An additional forty deidentified, unaltered drug free clinical oral fluid samples collected by expectoration and the Quantisal™ device were obtained from clinical research facilities, and analyzed using mass spectrometry. The results were negative for all forty samples collected by expectoration and Quantisal™ device.

DISPOSAL OF DEVICE

Dispose of all used devices in accordance with national and local disposal procedures.

Dispose of contents and container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

TECHNICAL ASSISTANCE

If technical assistance is required when using this product, please contact Immunalysis Corporation by post, telephone or email at the contact details below.

Region	Group Emails
Europe & Middle East	Phone: +44 161 483 9032 Option 3 Email: EME.TechSupport@abbott.com
Asia Pacific	Phone: +61 7 3363 7100 Email: AP.TechSupport@abbott.com
Africa	Phone: +27 10 500 9700 Option 3 Email: arcis.techsupport@abbott.com

Region	Group Emails
Russia, and the Commonwealth of Independent States (RCIS)	Phone: +44 161 483 9032 Option 3 Email: arcis.techsupport@abbott.com
Latin America	Phone: +57 1 482 4033 Email: LA.TechSupport@abbott.com
Canada	Phone: 1-800-818-8335 Email: Canproductsupport@abbott.com
US	Phone: +1 877 441 7440, Option 2 Email: USproductsupport@abbott.com

SERIOUS INCIDENT NOTICE

For users in the European Union and in countries with identical regulatory regimes (Regulation (EU) 2017/745 on Medical Devices): if a serious incident occurs during the use of this device or as a result of its use, please report it to the manufacturer and/or its authorized representative and to your national regulatory authority.

GLOSSARY OF SYMBOLS


Symbol	Description	ISO 15223-1 ¹ Reference #
	Authorized representative in the European Community	5.1.2
	Batch code	5.1.5
	Catalogue number	5.1.6
	CE Marking of Conformity	Not Applicable
	Consult instructions for use	5.4.3
	Do not re-use	5.4.2
	EC Importer: Registered Importer for products distributed within the European Community and the European Free Trade Association only	Not Applicable
	Keep away from sunlight	5.3.2
	Manufacturer	5.1.1

¹ ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

GLOSSARY OF SYMBOLS CONTINUED

Symbol	Description	ISO 15223-1 ¹ Reference #
	Medical device	Not Applicable
	QR Code for Unique Device Identifier (UDI)	Not Applicable
	Temperature limit	5.3.7
	Use by date	5.1.4

¹ ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

Abbott Rapid Dx International Limited
Parkmore East Business Park,
Ballybrit, Galway, H91 VK7E, Ireland
+353 (0)91 429 900

IMMUNALYSIS Corporation
829 Towne Center Drive, Pomona, CA 91767 USA
+1 (909)-482-0840
www.immunalysis.com

FOR EXPORT ONLY (from USA)

Indication of Change:

Revision AD	Reformatted size of the Instructions for Use.
Previous Revision AC	Added GTIN, Updated SKUs, Intended purpose, Warnings and Precautions, Specimen Stability, Specimen Collection, Instructions for Laboratory, Disposal of Device, technical assistance, Serious Incident Notice, Glossary of Symbols. Removed User requirements, Limitations. Added QSI-0025NINT, QSI-0500NINT. Updated intended purpose. Added H412 and P272 statements. Added shipping temperature limitation. Revised disposal, technical assistance and serious incident section. Removed ISO 7000 references. Added distributor, EC Importer symbol. Updated authorized representative to Abbott Rapid Dx International Limited. Added indications for change section.
Previous Revision AB	Added language headers. Updated contents section. Added French, Italian, German, Spanish, Portuguese, Swedish, Norwegian, Dutch languages.
Previous Revision AA	Created Instructions for Use

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© 2022 Immunalysis Corporation. Tous droits réservés. Toutes les marques mentionnées sont