SURESTEP™

URINE TEST

DRUG SCREEN PANEL

PGB Rapid Test Panel (Urine)

Package Insert

REF AB-DPG-114 **English**

*A rapid test for the qualitative detection of Pregabalin in human urine.
For medical and other professional in vitro diagnostic use only.*

【INTENDED USE】

The PGB Rapid Test Panel (Urine) is a rapid chromatographic immunoassay for the detection of Pregabalin in urine at a cut-off concentration of 500ng/ml. This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/ mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result.

【SUMMARY】

Pregabalin, also known as β-isobutyl-γ-amino butyric acid (beta-isobutyl-GABA), is a medication used to treat epilepsy, neuropathic pain, fibromyalgia, and generalized anxiety disorder.^[3] Common side effects include: sleepiness, confusion, trouble with memory, poor coordination, dry mouth, problem with vision, and weight gain. Potentially serious side effects include angioedema, drug misuse, and an increased suicide risk.^[2]

Pregabalin is eliminated from the systemic circulation primarily by renal excretion as unchanged drug. The Pregabalin is predominantly excreted unchanged in urine (≥ 98%)^[3]. Pregabalin mean elimination half-life is 6.3 hours.^[4] 50% would be expected to have negative urine specimens within 3 days and a total of 5 days would be needed to achieve negative urine specimens in the subject with the maximum urinary concentration measured.^[5]

The PGB Rapid Test Panel (Urine) is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Pregabalin in urine. The PGB Rapid Test Panel (Urine) yields a positive result when Pregabalin in urine exceed 500ng/mL.

【PRINCIPLE】

The PGB Rapid Test Panel (Urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Pregabalin, if present in the urine specimen below the cut-off level, will not saturate the binding sites of the antibody in the test. The antibody coated particles will then be captured by immobilized Pregabalin-protein conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Pregabalin level exceeds the cut-off level, because it will saturate all the binding sites of anti- Pregabalin antibody.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

【PRECAUTIONS】

- For medical and other professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

【STORAGE AND STABILITY】

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

【SPECIMEN COLLECTION AND PREPARATION】

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

【MATERIALS】

Materials Provided

- Test Panels

Materials Required But Not Provided

- Specimen collection containers
- Timer

【DIRECTIONS FOR USE】

Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

- Remove the test Panel from the sealed pouch and use it within one hour.
- Remove the cap.
- With the arrow pointing toward the urine specimen, immerse the test Panel vertically in the urine specimen for **at least 10 to 15 seconds**. Immerse the panel to at least the level of the wavy lines, but do not touch the plastic device.
- Replace the cap and place the test Panel on a non-absorbent flat surface.
- Start the timer and wait for the colored line(s) to appear.
- The result should be read at 5 minutes.** Do not interpret the result after 10 minutes.

【INTERPRETATION OF RESULTS】

(Please refer to the illustration)


NEGATIVE: Two colored lines appear. One colored line should be in the control line region (C), and another colored line should be in the test line region (T). This negative result indicates that the Pregabalin concentration is below the detectable cut-off level.

***NOTE:** The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: One colored line appears in the control line region (C). No line appears in the test line

region (T). This positive result indicates that the Pregabalin concentration exceeds the detectable cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

【LIMITATIONS】

- The PGB Rapid Test Panel (Urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/ mass spectrometry (GC/MS) is the preferred confirmatory method.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.
- The PGB Rapid Test Panel can only be used with urine specimen.

【EXPECTED VALUES】

This negative result indicates that the Pregabalin concentration is below the detectable level of 500ng/ml. Positive result means the concentration of Pregabalin is above the level of 500ng/ml. The PGB Rapid Test Panel has a sensitivity of 500ng/ml.

【PERFORMANCE CHARACTERISTICS】

A side-by-side comparison was conducted using the PGB Rapid Test Panel and GC/MS at the cut-off of 500ng/ml. Testing was performed on 75 clinical specimens previously collected from subjects present for Drug Screen Testing. The following results were tabulated:

Method	GC/MS		Total Results
	Results		
	Positive	Negative	
PGB Rapid Test Panel	Positive	20	22
	Negative	1	53
	Total Results	21	75
% Agreement	95.2%	96.3%	96.0%

Analytical Sensitivity

A drug-free urine pool was spiked with Pregabalin at the following concentrations: 0 ng/ml, 250ng/ml, 375ng/ml, 500ng/ml, 625ng/ml, 750ng/ml and 1500ng/ml. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Pregabalin Concentration (ng/ml)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
250	-50%	30	30	0
375	-25%	30	25	5
500	Cut-off	30	15	15
625	+25%	30	6	24
750	+50%	30	0	30
1500	3X	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the PGB Rapid Test Panel (Urine) at 5 minutes.

Compound	Concentration (ng/mL)
Pregabalin	500

Precision

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Pregabalin, 25% Pregabalin above and below the cut-off, and 50% Pregabalin above and below the 500ng/ml cut-off was provided to each site. The results are given below:

Pregabalin Concentration (ng/ml)	n per Site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
250	10	10	0	10	0	10	0

375	10	9	1	8	2	8	2
625	10	2	8	2	8	2	8
750	10	0	10	0	10	0	10
1500	10	0	10	0	10	0	10

Effect of Urinary Specific Gravity

Fifteen urine specimens of normal, high, and low specific gravity ranges were spiked with 250ng/ml and 750ng/ml of Pregabalin. The PGB Rapid Test Panel (Urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Pregabalin to 250ng/ml and 750ng/ml. The spiked, pH-adjusted urine was tested with The PGB Rapid Test Panel (Urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Pregabalin positive urine. The following compounds show no cross-reactivity when tested with The PGB Rapid Test Panel (Urine) at a concentration of 100µg/ml.

Non Cross-Reacting Compounds	
Acetaminophen	d,l-Chloropheniramine
N-Acetylprocainamide	Chloroquine
Ampicillin	Trans-2-phenylcyclopropylamine
	hydrochloride
Apomorphine	l-Cotinine
Atropine	Deoxycorticosterone
Benzoic acid	Diclofenac
d,l-Brompheniramine	Digoxin
Chloral-hydrate	l -Ψ-Ephedrine
Chlorothiazide	Estrone-3-sulfate
Chlorpromazine	l(-)-Epinephrine
Cholesterol	Fenoprofen
Cortisone	Genistic acid
Creatinine	Hydralazine
Dextromethorphan	Hydrocortisone
Diflunisal	p-Hydroxytyramine
Diphenhydramine	Iproniazid
β-Estradiol	Isosuxiprine
Ethyl-p-aminobenzoate	Isotamine
Erythromycin	Labetalol
Furosemide	Meperidine
Hemoglobin	Methylphenidate
Hydrochlorothiazide	Naproxen
o-Hydroxyhippuric acid	Nifedipine
lbuopren	d-Norpropoxyphene
d,l-Isoproterenol	d,l-Octopamine
Acetophenetidin	Oxolinic acid
Acetylsalicylic acid	Papaverine
Amoxicillin	Pentazocine hydrochloride
	Phenelzine
l-Ascorbic acid	Phenylpropanolamine
Aspartame	Prednisone
Benzoic acid	d-Propoxyphene
Benzphetamine	Quinacrine
Caffeine	Quindine
Chloramphenicol	Zomepirac
Salicylic acid	

【BIBLIOGRAPHY】

- Frampton, JE (September 2014). "Pregabalin: a review of its use in adults with generalized anxiety disorder." CNSDrugs 28 (9): 835-54.
- "Pregabalin". The American Society of Health-System Pharmacists. Retrieved Oct 23,2015.
- D.R. Guay. Pregabalin in neuropathic pain: a more "pharmaceutically elegant" gabapentin? Am. J. Geriatr. Pharmacother. 3: 274-287 (2005).
- "Summary of product characteristics" (PDF). European Medicines Agency. 6 March 2013. Retrieved 6 May 2013.
- "Detection times of Pregabalin in urine after illicit use: when should a positive specimen be considered a new intake?" TherDrugMonit. 2013 Feb;35(1):137-40.

Index of Symbols			
	Consult Instructions For Use		Tests per kit
	For in vitro diagnostic use only		Use by
	Store between 2-30°C		Lot Number
	Do not use if package is damaged		Keep away from sunlight
	Authorized Representative		Keep dry

Hangzhou AllTest Biotech Co.,Ltd.
#550,Yinhai Street
Hangzhou Economic & Technological Development Area
Hangzhou - 310018,P.R.China
www.alltests.com.cn

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LITERATUR

- Frampton, JE (September 2014). "Pregabalin: a review of its use in adults with generalized anxiety disorder." CNS Drugs 28 (9): 835-54.
- "Pregabalin". The American Society of Health-System Pharmacists. Retrieved Oct 23,2015.
- D.R. Guay. Pregabalin in neuropathic pain: a more "pharmaceutically elegant" gabapentin? Am. J. Geriatr. Pharmacother. 3: 274-287 (2005).
- "Summary of product characteristics" (PDF). European Medicines Agency. 6 March 2013. Retrieved 6 May 2013.
- "Detection times of Pregabalin in urine after illicit use: when should a positive specimen be considered a new intake?" TherDrugMonit. 2013 Feb;35(1):137-40.

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[BIBLIOGRAFIA]
1. Frampton, JE (September 2014). "Pregabalin: a review of its use in adults with generalized anxiety disorder."
2. "Pregabalin". The American Society of Health-System Pharmacists. Retrieved Oct 23,2015.

Table with 2 columns: Compound Name, Pregabalin Concentration. Rows listing various drugs and their corresponding Pregabalin concentrations.

- 1. Frampton, JE (September 2014). "Pregabalin: a review of its use in adults with generalized anxiety disorder."
2. "Pregabalin". The American Society of Health-System Pharmacists. Retrieved Oct 23,2015.
3. D.R. Guay. Pregabalin in neuropathic pain: a more "pharmaceutically elegant" gabapentin?
4. "Summary of product characteristics" (PDF). European Medicines Agency. 6 March 2013. Retrieved 6 May 2013.
5. "Detection times of Pregabalin in urine after illicit use: when should a positive specimen be considered a new intake?"

Table with 2 columns: Symbol, Description. Rows for instructions like 'Consultare le istruzioni per l'uso', 'Solo per uso diagnostico in vitro', 'Conservare tra 2-30°C', etc.

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- Frampton, JE (September 2014). "Pregabalin: a review of its use in adults with generalized anxiety disorder." CNSDrugs 28 (9): 835–54.
- "Pregabalin". The American Society of Health-System Pharmacists. Retrieved Oct 23,2015.
- D.R. Guay. Pregabalin in neuropathic pain: a more "pharmaceutically elegant" gabapentin? Am. J. Geriatr. Pharmacother. 3: 274–287 (2005).
- "Summary of product characteristics" (PDF). European Medicines Agency. 6 March 2013. Retrieved 6 May 2013.
- "Detection times of Pregabalin in urine after illicit use: when should a positive specimen be considered a new intake?" TherDrugMonit. 2013 Feb;35(1):137-40.

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