
Air Hoses for Noninvasive Blood Pressure Cuffs

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




989803209741, M1596C, 989803209751, M1597C

989803209761, M1598B, 989803209771, M1599B

EN - Air Hoses for NBP Cuffs Instructions for Use

Associated Documents - NBP Cuffs Instructions for Use; NBP Cuffs and Air Hoses Care, Cleaning and Disinfection, and Multi Care NBP Cuffs Instructions for Use. **Intended Use** - Philips air hoses are to be used with compatible devices intended for use by, or under the supervision of, a licensed physician or other healthcare provider for noninvasive measurement of neonatal, pediatric and adult human blood pressure. **Intended Purpose** - In accordance with European Medical Device Regulation 2017/745, this accessory device enables a vital signs patient monitoring instrument to record and generate alarms for multiple physiological parameters when monitoring adults, pediatrics, and neonates. **Indications** - These devices are limited by the indications for use of the connected monitor in healthcare facilities. These devices are intended to interact with patient intact skin only. **Contraindications** - There are no known contraindications.

Before use, read all warnings and cautions listed in this IFU. Also refer to the IFU for the monitor or defibrillator for additional warnings and cautions.

Warnings - Use neonatal and infant #5 cuffs only with hoses designated for Neonatal use (see cover and *Selecting the Appropriate Hose*). • Use pediatric and adult cuffs only with hoses designated for Adult use (see cover and *Selecting the Appropriate Hose*). • Hoses are not designed to be sterilized. Do not autoclave. • Refer to the *NBP Cuffs and Air Hoses Care, Cleaning and Disinfection* instructions for validated cleaners and disinfectants and appropriate procedures to prevent cross-contamination between patient use. • Ensure that liquids do not enter the air hose since any liquid in the air hose can damage monitoring equipment. • Philips NBP air hoses use connectors which have the potential to misconnect with small-bore connectors used in non-NBP healthcare applications. When connecting Philips NBP air hose, ensure that connectors are properly connected to Philips NBP devices only. • When cleaning the hoses, use care to prevent liquid from entering the connectors. Do not immerse the hose in any liquid since this action may lead to cuff or monitor failure. • Modification of this equipment is not allowed.

Cautions - Federal (US) law restricts this device to sale by or on the order of a medical practitioner. • Do not substitute any component of the NBP product since a substitution might result in a measurement error.

Visual Inspection Prior to Use - Before using the hose on a patient, perform a visual inspection to determine if the hose is at its end of life. Visually inspect the hose for cracks, damaged connectors, or similar wear or damage. If damage or severe soiling occurs, see *Disposal*. **Selecting the Appropriate Hose** - Neonatal, Infant #5 cuff use 989803209741, M1596C, 989803209751, or M1597C; Infant, Pediatric or Small Adult, Adult, Large Adult, Thigh cuff use 989803209761, M1598B, 989803209771, or M1599B. **Compatibility** - Only use hoses on systems that list them as accessories. **Reordering Information** - The following part numbers are CE marked and available in the European Union: 989803209741 and M1596C — Neonatal NBP Air Hose 1.5 m; 989803209751 and M1597C — Neonatal NBP Air Hose 3.0 m; 989803209761 and M1598B — Adult NBP Air Hose 1.5 m; 989803209771 and M1599B — Adult NBP Air Hose 3.0 m. **Cleaning and Disinfection** - Refer to the *NBP Cuffs and Air Hoses Care, Cleaning and Disinfection* instructions for a list of validated cleaners and disinfectants and the appropriate procedures. When cleaning cuffs, use care to prevent liquid from entering the tubing. Do not immerse the tubing. **Disposal** - To determine if a product has reached its end of life, visually check for cracks, damaged connectors, or similar wear or damage. Whenever the inspection reveals the product is no longer suitable for use, then follow approved medical waste disposal methods as specified by your patient care facility or local regulations. **EU Incident Reporting** - Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. **Environmental Specifications - Storage**: Temperature: -20° to 60° C (-4° - 140° F); Humidity: 0 to 90% relative humidity at 60° C (140° F); **Operating**: Temperature: 0° to 50° C (32° - 122° F) Humidity: 0 to 90% relative humidity.

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CS - Air Hoses for NBP Cuffs Instructions for Use

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Copyright © 2022
Böblingen GmbH