





LoFlo Accessories - Instructions for Use

Инструкции за употреба на аксесоари LoFlo	LoFlo piederumu lietošanas instrukcija
Príslušenství LoFlo – návod k použití	LoFlo додатки - Упатство за употреба
LoFlo-tilbehør – Brugsanvisning	Gebruiksaanwijzing voor LoFlo-accessoires
Gebrauchsanweisung LoFlo-Zubehör	Bruksanvisning for LoFlo-tilbehør
Οδηγίες χρήσης εξαρτημάτων LoFlo	Instrukcja stosowania akcesoriów LoFlo
Instrucciones de uso de los accesorios LoFlo	Instruções de uso de acessórios LoFlo
Tootesarja LoFlo lisatarvikute kasutusjuhend	Instrucțiuni de utilizare pentru accesorii LoFlo
LoFlo-lisävarusteet – käyttöohjeet	Инструкции по применению аксессуаров LoFlo
Mode d'emploi des accessoires LoFlo	Návod na použitie príslušenstva LoFlo
Pribor LoFlo – upute za uporabu	Navodila za uporabo za dodatke LoFlo
A LoFlo tartozékainak használati útmutatója	Упутство за употребу за LoFlo додатке
Petunjuk Penggunaan Aksesoris LoFlo	Bruksanvisning för LoFlo-tilbehör
Istruzioni per l'uso degli accessori LoFlo	LoFlo Aksesuarları Kullanım Talimatları
LoFloアクセサリ使用方法	Інструкції із застосування аксесуарів LoFlo
LoFlo керек-жарақтары - Пайдалану нұсқаулары	Hướng dẫn sử dụng phụ kiện LoFlo
LoFlo 부속품 사용 지침	LoFlo 附件使用说明
„LoFlo“ priedų naudojimo instrukcijos	LoFlo 配件使用說明

LoFlo Accessories Instructions for Use

EN

English

Refer to the Instructions for Use accompanying your host monitor. Also refer to page 2 and 3 for a setup diagram.

NOTE: Any serious incident that has occurred in relation to this device should be reported to Respirationics Novamatrix and the competent authority of the country in which the user and/or patient is established.

Symbols

These symbols may be found on LoFlo accessories, labels, packaging, and other documentation.

	Do not reuse		Consult Instructions for Use
	Quantity		This device complies with Council Directive 93/42/EEC
	Storage humidity		This device complies with Regulation (EU) 2017/745
	Storage temperature		Adult only
	U.S. Federal law restricts this device to sale by or on the order of a physician.		Adult/Pediatric only
	ETT > 4.0mm Deadspace < 6ml Adult Pediatric		Pediatric only
	ETT ≤ 4.0mm Deadspace < 1ml Infant Neonatal		Neonatal/Infant only
	Importer		Medical device
	Electronic instructions for use		MR Unsafe

Cautions

- Single-use products are not designed to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy and/or system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization and/or reuse.
- Verify that the sampling cannula, airway adapter, or sampling line is intact and undamaged.
- For cannulas with an O₂ connector, verify oxygen flow prior to use.
- Cleaning and sterilization will compromise system performance.
- Remove sample cell from monitoring system when not in use.
- Do not use on patients that cannot tolerate the removal of 50ml/min from their total minute ventilation.
- Do not position cables or tubing in any manner that may cause entanglement or strangulation.

Accessory - LoFlo	Intended Use/Indications for Use:
CO ₂ Sampling Cannula, Nasal (3468ADU-00, 3468PED-00, 3468INF-00)	For use with the LoFlo Sidestream CO ₂ Sensor when monitoring patients that are not intubated.
CO ₂ Sampling Cannula w/ Dehumidification Tubing, Nasal (3468ADH-00, 3468PEH-00, 3468INH-00)	
CO ₂ w/O ₂ Cannula, Nasal (3469ADU-00, 3469PED-00, 3469INF-00)	
CO ₂ w/O ₂ Cannula w/ Dehumidification Tubing, Nasal (3469ADH-00, 3469PEH-00, 3469INH-00)	
CO ₂ Sampling Cannula, Oral/Nasal (3470ADU-00, 3470PED-00)	
CO ₂ Sampling Cannula w/ Dehumidification Tubing, Oral/Nasal (3470ADH-00, 3470PEH-00)	
CO ₂ w/O ₂ Cannula, Oral/Nasal (3471ADU-00, 3471PED-00)	
CO ₂ w/O ₂ Cannula w/ Dehumidification Tubing, Oral/Nasal (3471ADH-00, 3471PEH-00)	For use with the LoFlo Sidestream CO ₂ Sensor when monitoring patients with an endotracheal tube size larger than 4.0mm.
Airway Adapter Kit, ETT >4.0mm, Deadspace <6ml (3472ADU-00)	
Airway Adapter Kit w/ Dehumidification Tubing, ETT >4.0mm, Deadspace <6ml (3473ADU-00)	For use with the LoFlo Sidestream CO ₂ Sensor when monitoring patients with an endotracheal tube size of 4.0mm or smaller.
Airway Adapter Kit, ETT ≤4mm, Deadspace <1ml (3472INF-00)	
Airway Adapter Kit w/ Dehumidification Tubing, ETT ≤4mm, Deadspace <1ml (3473INF-00)	For use with the LoFlo CO ₂ Sensor with ventilator and anaesthesia circuits that have an integrated airway adapter with a female Luer-Lok connector.
Sampling Line Kit, Male Luer (3474-00)	
Sampling Line Kit w/Dehumidification Tubing, Male Luer (3475-00)	For adapting CapnoTrak accessories for use with a LoFlo Module.
Sample Line, Female Luer Connector (3476-00)	

Disposal

Dispose of in accordance with all applicable federal, state, and local medical / hazardous waste practices.

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Respironics Novamatrix, LLC
Respironics Deutschland GmbH & Co. KG
Philips Electronics Australia Ltd