
Microstream™ Advance

Adult-Pediatric Intubated CO₂ Filter Line

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


EN Instructions for Use

Intended Use

The Microstream™ Advance Intubated CO₂ Filter Line is intended to conduct a sample of the patient's breathing from a ventilator or anesthesia machine to a gas measurement device for measuring the percentage of CO₂ in the patient's exhaled breath. The set is intended for single patient use only. Intended population: Intubated Adult-Pediatric patients.

Intended Users

The Microstream™ Advance Intubated CO₂ Filter Line is intended to be used by licensed healthcare providers only.

Indications for Use

The Microstream™ Advance Intubated CO₂ Filter Line is indicated for measurement of exhaled carbon dioxide (CO₂) for monitoring purposes. These devices are limited by the indications for use of the connected Philips monitoring equipment in healthcare facilities. These devices are intended to interact with patient intact skin and nasal mucosa.

Contraindications

There are no known contraindications.

WARNINGS

- Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components, turning the sampling line connector clockwise into the monitor CO₂ port until it can no longer be turned. Check connections for leaks according to standard clinical procedures, and ensure CO₂ values appear.
- Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line. Re-use of single-use accessories could pose a cross-contamination risk to the patient or damage the functioning of the monitor.
- When used with a suction catheter, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter, and that the suction process does not damage the airway adapter. Damage of the airway adapter may harm the patient.
- Ensure no kinks are present in CO₂ tubing, as kinked tubing may cause inaccurate CO₂ sampling.
- Carefully route patient cabling (CO₂ sampling line) to reduce the possibility of patient entanglement or strangulation.

Caution

Federal (US) law restricts this device to sale by or on the order of a medical practitioner.

Required Compatible Components

These devices are to be used with Philips monitors using Microstream™ technology.

Airway adapters for adult to pediatric patients:

- Added dead space <6.6 cc.
- To be used with an endotracheal tube > 4.5mm bore.
- Suitable for use with 22mm and 15mm patient tubing.

Applying Adult/Pediatric Intubated CO₂ Filter Line

When performing the following procedure, reference the step-by-step patient application graphics on Page i:

- ① This graphic displays the in-line setup for an intubated CO₂ measurement.
- ② Firmly connect the small-end (male) of the Microstream™ airway adapter to the female-end of the Wye piece (ventilation source).
- ③ Firmly connect the patient's endotracheal tube connector into the large-end (female) of the Microstream™ airway adapter.
- ④ This graphic displays the setup of intubated CO₂ measurement with suction catheter. Be sure to place the airway adapter as shown, so that a suction catheter does not interfere with the functioning of the airway adapter
- ⑤ Endotracheal tube
- ⑥ Filter Line Set with airway adapter.
- ⑦ Wye piece (ventilation source).
- ⑧ Ventilator connector on closed suction system, with lavage.
- ⑨ Suction tubing
- ⑩ Attach CO₂ connector to monitor as shown in this graphic.

Operational Notes

- During setup, ensure that airway adapter can be easily attached and detached from the tubing before proceeding.
- During nebulization, lavage, or suction, follow monitor instructions for use (IFU) to suspend pump activity, in order to avoid moisture buildup and sampling line occlusion. In addition, especially if unable to suspend the CO₂ pump, disconnect the sampling line connector from the CO₂ port on the monitor.
- The response time of the 2m sampling line is ≤3.5 seconds, 4m sampling lines ≤5.9 seconds..
- Sampling lines with Extended Duration or High Humidity in their names include a moisture reduction component (Nafion®* or its equivalent) for use in higher humidity environments where long duration use of CO₂ sampling is required.
- Verify that CO₂ readings appear on monitor screen as expected after insertion.

TM Trademark of its respective owner.*

Reordering Information

Part Number	Description	Quantity
989803204521	Adult-Pediatric Intubated CO ₂ Filter Line, Extended Duration Note: Length: 2m	25
989803204301	Adult-Pediatric Intubated CO ₂ Filter Line, High Humidity Note: Length: 2m	25
989803204331	Adult-Pediatric Intubated CO ₂ Filter Line, Extended Duration Note: Length: 4m	25
989803204511	Adult-Pediatric Intubated CO ₂ Filter Line, Short-term use: Procedural/Emergency Note: Length: 2m	25
989803204321	Adult-Pediatric Intubated CO ₂ Filter Line, Short-term use: Procedural/Emergency Note: Length: 4m	25
989803204641	Adult-Pediatric Intubated CO ₂ Filter Line, Short-term use: Procedural/Emergency Note: Length: 2m	25
989803204651	Adult-Pediatric Intubated CO ₂ Filter Line, Extended Duration Note: Length: 2m	25
989803204671	Adult-Pediatric Intubated CO ₂ Filter Line, High Humidity Note: Length: 2m	25

Disposal

Discard any product that becomes damaged or deteriorated. Follow approved medical waste disposal methods specified by your facility or local regulations.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey, in which the user and/or patient is established.

Environmental Specifications

The following table provides the temperature and humidity ranges that must be maintained for these products during storage and transport:

Specification	Storage	Transport
Temperature	-20°C to 70°C (-4°F to 158°F)	-20°C to 70°C (-4°F to 158°F)
Humidity	10% RH to 90% RH	10% RH to 90% RH

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Copyright © 2020
Microstream™ is a trademark of Oridion
Nederland B.V.
Medical 1987 Ltd. Oridion Medical 1987 Ltd.