

## Expression Model **MR400**

### MRI Patient Monitoring System



The Expression Model MR400 MRI Patient Monitoring System is designed to assist clinicians in monitoring patient vital signs in the dynamic magnetic resonance environment. The Expression Model MR400 MRI Patient Monitoring System combines wireless communication, radio frequency shielding and digital signal processing to address the challenges associated with patient monitoring in the MR environment.

The Expression Model MR400 MRI Patient Monitoring System consists of these primary components:

- Traditional Roll-Around Cart
- Wireless ECG (wECG) module
- Wireless SpO<sub>2</sub> (wSpO<sub>2</sub>) module

## Features and Benefits

- Integral color, 15" LED widescreen display for high resolution patient information
- Intuitive touchscreen graphical user interface
- Colored waves and large numeric displays
- Bedside-quality "SINC" parameters
- Exclusive, advanced ECG solution for MRI
- 8-hour battery life and user-replaceable batteries for extended run time
- Simultaneous display of up to 13 parameters, and 6 waveforms and associated values
- Multi-priority visual and audible alarm signals, unique "alarm flag" messages, and pulse tones
- Gating, both digital pulse and analog waveform

## Optional Components

- Expression Information Portal (IP5)
- Wireless remote printing to IP5

## Ordering Information

### Standard Features, 866185

**A01:** Standard Accessories

### Options

**F01:** Basic (NBP, ECG, SpO<sub>2</sub>, CO<sub>2</sub>, RR)

**F02:** Basic + IBP (x2)

**F03:** Basic + Temp

**F04:** Basic + AA, O<sub>2</sub>

**F05:** Basic + AA, O<sub>2</sub>, IBP (x2)

**F06:** Basic + AA, O<sub>2</sub>, Temp

**F07:** Basic + AA, O<sub>2</sub>, IBP (x2), Temp

## System Parameters

The Expression Model MR400 MRI Patient Monitoring System can include the following vital sign parameters:

- Electrocardiogram (ECG), dual channel
- Blood oxygen saturation/pulse oximetry (SpO<sub>2</sub>)
- Invasive blood pressure (IBP)
- Non-invasive blood pressure (NIBP)
- End-tidal and inspired CO<sub>2</sub>
- Respiration from CO<sub>2</sub> or bellows
- Anesthetic Agents, including end-tidal and inspired N<sub>2</sub>O, inspired O<sub>2</sub>, and Total MAC
- Temperature

The system can include the ability to display these parameters:

- Alarms: High and low selectable limits for each patient parameter
- ECG: Waveform scale, dual channels displayed

- Heart rate: Factory-default derived from ECG; also from pulse oximetry or IBP
- Pulse oximeter: Pulse rate, pulse waveform, and percent saturation
- CO<sub>2</sub>: End-tidal and inspired
- IBP (two channels, P1 and P2): Systolic, mean, and diastolic pressures
- NIBP: Systolic, mean, and diastolic pressures
- Anesthetic Agents, including end-tidal and inspired N<sub>2</sub>O, and inspired O<sub>2</sub>
- Bellows respiration: Rate derived from pneumatic chest bellows
- Temperature
- Trends: Heart rate, respiration rate, IBP (systolic, diastolic, mean), NIBP (systolic, diastolic, mean), CO<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>O, SpO<sub>2</sub>, agents, and temperature
- Respiration: Rate derived from CO<sub>2</sub>
- Time: Battery-backed quartz clock

## Main Component

### Display Panel

**Type:** Liquid Crystal Display (LCD), touch screen, color

**Screen size:** 39.5 cm (15.6 inches) diagonal

**Drive type:** a-Si TFT active matrix

**Pixels:** 1366 (H) by 768 (V) pixels, color

**Area:** 344.2 (H) by 193.5 (V) mm

**Dot pitch:** 0.084 (H) by 0.252 (V) mm

**Pixel pitch:** 0.252 (H) by 0.252 (V) mm

**Contrast ratio:** 500:1 (typical)

**Backlight:** LED

**Polarizer surface:** Anti-glare

**Tilt:** Adjustable, 5° to 35°

**Sweep speeds for ECG, SpO<sub>2</sub>, and IBP:** 25 mm/second gives 9.1 seconds of display time, while 50 mm/second gives 4.6 seconds. Sweep speeds for respiration: 3.125, 6.25, 12.5 or 25 mm/second are provided.

**Waveform display mode:** Fixed trace, moving erase bar

**Waveform display length:** ≥ 228 mm

**Waveform display height:** ECG (single trace): ≥ 40 mm; ECG (dual trace): ≥ 20 mm; All other waveforms: ≥ 25 mm

**Speaker:** Audio

**Control of monitored parameters provided by:** Power switch; Touchscreen

### User Interface

**Four groups of data are displayed:** Informational, Vital sign traces, Vital sign numeric values, System status

## Application Features

### Trends

- Data at 1-minute intervals, up to 12 hours of data
- Automatically can store the parameter trend information for heart rate, IBP, NIBP, SpO<sub>2</sub>, CO<sub>2</sub>, N<sub>2</sub>O, O<sub>2</sub>, agents, respiration, and temperature

- Trend arrows graphically indicate an increasing, decreasing or stable parameter
- Graphical trends (with the IP5 option)

## Alarms

**Severity:** High, medium, and low

**Visual indicators:** Alarm light, flashing numeric values, alarm flags, icons

**Audible alarms:** User-configurable for volume, tone, and silence

**Limits:** Configurable; 1-Touch Alarms allow limits to be quickly adjusted

## Device Connections

**Input/output ports permit the connection of external**

**equipment:** USB port (system update use only); ECG and peripheral gating output port

## Specifications

### Safety Standards

- Conforms to ANSI/AAMI ES 60601-1: 2012. Certified to CAN/CSA C22.2 No. 60601-1-08; IEC 60601-1-2
- Conforms to 93/42/EEC as amended by 2007/47/EEC, *Medical Device Directive*
- Defibrillator protection up to 5 KV

### Physical Specifications

Cart

**Size:** 127.3 x 47.5 x 55.9 cm (50.1 x 18.7 x 22 in.)

**Weight:** 46.9 kg (103.3 Lbs)

Wireless ECG module

**Size:** 18.2 x 6.7 x 3.1 cm (7.17 x 2.65 x 1.24 in.)

**Weight:** 340 g (12 oz.)

Wireless SpO<sub>2</sub> module

**Size:** 13.0 x 6.5 x 3.1 cm (5.13 x 2.55 x 1.24 in.)

**Weight:** 204 g (7.2 oz.)

### Electrical Specifications

Power Requirements

**Operating voltage range:** 100 to 240 VAC

**Frequency range:** 50 to 60 Hz

**Current:** 1.4 A @ 100 VAC / 0.7 A @ 240 VAC

**Maximum power consumption:** ≤ 65 Watts

Battery Type

**Cart:** Lithium-Ion

**Module:** Lithium polymer

Battery Operation Time - Cart

**All displays, alarms, and monitoring functions:** continuously for 8 hours

**ECG & SpO<sub>2</sub>:** 8 hours

**CO<sub>2</sub>** (with or without AGENT): continuously for 6 hours

**Invasive Pressures:** continuously for 6 hours

**AGENT analysis:** continuously for 6 hours

**Temperature:** continuously for 6 hours

**NIBP:** readings every 5 minutes for 6 hours

Battery Operation Time - Module

**Duration:** Approximately 8 hours

Battery Capacity

**Cart:** 75 Wh/battery (300 Wh capacity with 4 batteries installed)

**Module:** 3.1 Wh

### Environmental Specifications

**Operating temperature range:** 10 to 35°C (50 to 95°F)

Relative humidity range

**Cart, bellows, SpO<sub>2</sub> accessories, LoFlo accessories (optional), gating cables (optional), and IP5 (optional):** 15 to 80 %, non-condensing

**Wireless modules, ECG cables, and FlexTEMP II sensor:** 5 to 80 %, non-condensing

**IBP transducer kits and cable (optional):** 15 to 85 %, non-condensing

**Module batteries:** 15 to 90 %, non-condensing

Storage and transport temperature range

**Cart:** -20 to 50°C (-4 to 122°F)

**Cart batteries:** 0 to 40°C (32 to 104°F)

**Wireless modules, and all other accessories not specified below:** -20 to 60°C (-4 to 140°F)

**Quadtrodes:** 10 to 32°C (50 to 90°F)

**Optional Transducer kits/cable:** -15 to 60°C (5 to 140°F)

**O<sub>2</sub> sensor (AGENT option), storage:** 5 to 25 °C (41 to 77 °F);

**transport:** -40 to 50 °C (-40 to 122 °F)

**ECG skin prep gel:** Follow instructions on tube.

**NOTE:** When storing or transporting in temperatures beyond the ranges specified above, remove the designated component and store or move it appropriately.

**Operating pressure range:** Up to 3,000 m (9,842 feet) above sea level (708 mbar or 531 mmHg)

**Storage and transport pressure range:** 708 to 1020 mbar

### MRI Rating

**MR Conditional:** 1.5 and 3.0 Tesla, 5000 gauss, at RF power levels not exceeding 4W/kg SAR and 7.2 μT B<sub>1rms</sub> in all orientations

### Measurement Specifications

Electrocardiogram Channel (ECG)

**ECG Amplifier:** Protected against defibrillator and electro-surgery potentials

**Standard lead configurations:** I, II, III, AVR, AVL, AVF

**Lead Fail:** Passive, sensing signal imbalance

**ECG input impedance:** > 2.5MΩ (according to IEC 60601-2-27, 50.102.3)

**Electrode contact impedance:**  $\leq 20\text{K ohms @ } 10\text{ Hz}$

#### Heart Rate

**Range:** 30 to 250 bpm (adult); 30 to 300 bpm (pediatric and neonate)

**Resolution:** 1 bpm

**Accuracy:**  $\pm 1\%$  or  $\pm 1\text{ bpm}$ , whichever is greater

#### Cardiotach

**Sensitivity (Monitor filter):**  $> 200\text{ }\mu\text{V}$  (Adult ECG mode);  $> 100\text{ }\mu\text{V}$  (Neonate/Pediatric ECG mode)

**Bandwidth:** Monitor (0.5 to 40 Hz)

**QRS Duration:** 70 to 120 ms (adult); 40 to 120 ms (Neonate/Pediatric)

**Baseline Offset:** Automatically removed

**Tall T-wave rejection capability for heart rate indication:** 2 mV with a 1 mV QRS amplitude

**Leads-off sensing:** Detection by DC current waveform of  $< 100\text{ nA}$

#### Alarm Limits

| Lower*        | Upper*        |
|---------------|---------------|
| 30 to 250 bpm | 60 to 250 bpm |

\* Off is an option for each alarm limit

#### Test/Calibrations

**Square wave test signal:** 60 bpm  $\pm 1\text{ bpm}$ , 1 mV  $\pm 10\%$

#### Pulse Oximeter

*(Pitch of pulse tone is modulated by saturation value.)*

**Saturation range:** 1 to 100 %

**Saturation value resolution:** 1 %

**Saturation accuracy:**  $\pm 3\%$  at 70 to 100 %

**Pulse accuracy:**  $\pm 2\%$  or  $\pm 1\text{ bpm}$ , whichever is greater

**Pulse rate range:** 30 to 250 bpm

**Pulse rate resolution:** 1 bpm

**Perfusion Index:** Decimal Number

**Data update period:** 5, 10, or 15 seconds (according to the SpO<sub>2</sub> Averaging Time setting)

**Data update period during alarm:** 9, 14, or 19 seconds, maximum (4 seconds plus the SpO<sub>2</sub> Averaging Time setting of 5, 10, or 15 seconds)

**Wavelength range:** 500 to 1000 nm (Information about wavelength range can be especially useful to clinicians)

**Emitted light energy:**  $< 15\text{ mW}$

**Pulse oximeter calibration range:** 70 to 100 %

#### Alarm Limits

|  | Lower*        | Upper*        |
|--|---------------|---------------|
| SpO <sub>2</sub>                           | 50 to 100 %   | 70 to 100 %   |
| When "HR" is derived from SpO <sub>2</sub> | 30 to 250 bpm | 60 to 250 bpm |

\* Off is an option for each alarm limit

#### CO<sub>2</sub> (Optional LoFlo)

*(Side stream non-dispersive infrared absorption technique, including multiple water trap filtration system and microprocessor control of sample handling and calibration.)*

**Method for determining end tidal CO<sub>2</sub> measurement:**

Measures peak of the expired CO<sub>2</sub> waveform every 20 seconds.

**Output:** CO<sub>2</sub> waveform, EtCO<sub>2</sub> and FiCO<sub>2</sub> numeric values, and respiration rate

**Initialization time:** Waveform displayed in less than 20 seconds, at an ambient temperature of 25°C (77°F); full specifications attained within 2 minutes

**Zero calibration interval:** Automatic or user requested  
**CO<sub>2</sub> unit of measure:** Millimeters of mercury (mmHg) or kilopascals\* (kPa)

\*For kilopascals (kPa), allow  $\pm 1$  least significant digit to accommodate round-off error for calculated values.

**CO<sub>2</sub> resolution:** 1 mmHg (0.1 kPa)

**Flow rate:** 50 mL per minute  $\pm 10\text{ mL per minute}$

**Data sample rate:** 100 Hz

**EtCO<sub>2</sub> measurement range** (in which the accuracy specification is met: 0 to 76 mmHg (0 to 10.1 kPa) for respiration rates ranging from 4 to 60 breaths per minute, inclusive

**Inspired CO<sub>2</sub> (FiCO<sub>2</sub>) measurement range:** 0 to 50 mmHg (0 to 6.7 kPa) (method: lowest reading of the CO<sub>2</sub> waveform in the previous 20 seconds)

**CO<sub>2</sub> accuracy:**  $\pm 4\text{ mmHg}$  ( $\pm 0.5\text{ kPa}$ ) or  $\pm 12\%$ , whichever is greater

**Short term drift:** Not to exceed 0.8 mmHg (0.1 kPa) over a 4-hour period

**Long term drift:** Accuracy specification maintained over a 120-hour period

**Respiration accuracy:**  $\pm 1$  breath or  $\pm 3\%$ , whichever is greater

**Respiration resolution:** 1 rpm

**Respiration rate range** (in which the respiration accuracy specification is met: 4 to 100 rpm, inclusive

**Accessory usage:** Functional without changing accessories for a minimum of 6 hours

**Response and rise times:**

*(As measured from the patient gas input of the complete pneumatic circuit, including tubing, from 10 to 90 percent of the measured CO<sub>2</sub> levels)*

|                 | System Response (sec) | Rise Time (sec) |
|-----------------|-----------------------|-----------------|
| Airway Adaptor  | 10.89                 | 0.94            |
| Cannula         | 12.44                 | 1.12            |
| Divided Cannula | 16.17                 | 2.01            |

#### Compensations

*(Automatic CO<sub>2</sub> ambient pressure compensation 523 to 760 mmHg [69.7 to 101.3 kPa]; for expired O<sub>2</sub> balance gas (N<sub>2</sub>, N<sub>2</sub>O, O, He) and anesthetic agents; Uses gas compensation information to correct the raw carbon dioxide value)*

**Anesthetic agent effects (MAC levels)**

**Uncompensated Sensitivity:** Accuracy maintained for halogenated anesthetic agents present at accepted Minimum Alveolar Concentration clinical levels

**Compensated Sensitivity:** Testing at agent levels defined by accepted regulatory standards (80601-2-55:2011)

**Cross-sensitivity compensation error:**

*(Additional worst case error when compensation for O<sub>2</sub>, N<sub>2</sub>O, anesthetic agents, or helium is correctly selected for the actual fractional gas constituents present)*

**0 to 40 mmHg:**  $\pm 1\text{ mmHg}$  (0 to 5.3 kPa:  $\pm 0.1\text{ kPa}$ )

**41 to 70 mmHg:**  $\pm 2.5\text{ mmHg}$  (5.5 to 9.3 kPa:  $\pm 0.3\text{ kPa}$ )

**71 to 100 mmHg:**  $\pm 4\text{ mmHg}$  (9.5 to 13.3 kPa:  $\pm 0.5\text{ kPa}$ )

**101 to 150 mmHg:**  $\pm 5\text{ mmHg}$  (13.5 to 20 kPa:  $\pm 0.6\text{ kPa}$ )

**Quantitative effects of gas sample humidity or condensate\*\*:**

**0 to 40 mmHg:**  $\pm 2\text{ mmHg}$  (0 to 5.3 kPa:  $\pm 0.2\text{ kPa}$ )

**41 to 70 mmHg:**  $\pm 5\%$  (5.5 to 9.3 kPa:  $\pm 5\%$ )

**71 to 100 mmHg:**  $\pm 8\%$  (9.5 to 13.3 kPa:  $\pm 8\%$ )

**101 to 150 mmHg:**  $\pm 10\%$  (13.5 to 20 kPa:  $\pm 10\%$ )

*\*\*With appropriate compensations applied*

#### Alarm Limits

|                   | Lower*                          | Upper*                           |
|-------------------|---------------------------------|----------------------------------|
| EtCO <sub>2</sub> | 5 to 60 mmHg;<br>0.6 to 8.0 kPa | 5 to 76 mmHg;<br>0.7 to 10.1 kPa |
| FiCO <sub>2</sub> | No low alarm<br>limit           | 0 to 20 mmHg;<br>0 to 2.7 kPa    |
| Respiration       | 4 to 40 rpm                     | 20 to 100 rpm                    |

\* Off is an option for each alarm limit

#### Invasive Pressure (Optional)

##### Pressure Amplifier

**Isolation voltage:** 5 KVDC

**Signal range:** -30 to 250 mmHg

**Sensitivity:** 5  $\mu\text{V/V/mmHg}$

**Gain accuracy:**  $\pm 0.5\%$  percent

**Bandwidth:** 0 to 10 Hz (-3 dB)

**Offset range:**  $\pm 300$  mmHg

##### Transducer (REF 989803179721)

**Operating pressure:** -50 to 300 mmHg

**Overpressure limits:** -400 to 5000 mmHg

**Sensitivity:** 5  $\mu\text{V/V/mmHg}$   $\pm 1$  @ 6 VDC and 22°C (71.6°F)

**Zero offset:**  $< \pm 25$  mmHg

**Zero drift:**  $< \pm 2$  mmHg in 8 hours

**Input impedance:** 300 to 350 ohms

**Output impedance:** 300 ohms  $\pm 30$  ohms

##### Auto Zero

**Range:** +300 mmHg

**Zero accuracy:**  $\pm 1.0$  mmHg

**Response time:** 1 second, notification upon completion

##### Pressure Wave Display

**Number of channels:** 0, 1 or 2

**ABP, PAP and LAP:** Numeric display of systolic, diastolic and mean pressures

**CVP and ICP:** Numeric display of mean pressure only

##### Pressure Scale

**Ranges:** (user selectable) 0 to 250 mmHg; 0 to 200 mmHg; 0 to 150 mmHg; 0 to 100 mmHg; 0 to 75 mmHg; 0 to 45 mmHg

##### Pulse Rate (When derived from P1 or P2)

**Range:** 30 to 250 bpm

**Accuracy:**  $\pm 2\%$  of full scale

**Resolution:** 1 bpm

##### Alarm Delay

**Transducer disconnect:** 6 seconds

**Pressure disconnect:** 6 seconds

**High and low pressure:** 10 seconds

#### Alarm Limits

|                                    | Lower*                               | Upper*                                |
|------------------------------------|--------------------------------------|---------------------------------------|
| When "HR" is derived from P1 or P2 | 30 to 250 bpm                        | 60 to 250 bpm                         |
| Systolic, Mean, Diastolic          | - 30 to 250 mmHg;<br>4.0 to 33.3 kPa | - 30 to 250 mmHg;<br>-4.0 to 33.3 kPa |

\* Off is an option for each alarm limit

#### Transducer Connector Pin Compatibility

**Pin A:** - Signal

**Pin B:** + Excitation

**Pin C:** + Signal

**Pin D:** - Excitation

**Pin E:** Shield

#### Anesthetic Agents (Optional)

*Side stream, non-dispersive infrared (NDIR) absorption technique, including water trap filtration system and microprocessor control of sample handling and calibration*

**Simultaneously measured gases** (any two of the following, inspired or expired, while also measuring CO<sub>2</sub>, N<sub>2</sub>O, and O<sub>2</sub>):  
Halothane, Isoflurane, Desflurane, Enflurane, Sevoflurane

#### Measurement Range

*(After maximum warm-up period)*

**Halothane:** 0 to 5.0 VOL%

**Isoflurane:** 0 to 5.0 VOL%

**Desflurane:** 0 to 18.0 VOL%

**Enflurane:** 0 to 5.0 VOL%

**Sevoflurane:** 0 to 8.0 VOL%

**Carbon dioxide:** 0 to 10.0 VOL%

**Nitrous oxide:** 0 to 100 VOL%

#### Accuracy

*(Includes stability and drift)*

**Halothane:** 0 to 1.00 VOL% ( $\pm 0.15$  VOL%); 1.00 to 5.00 VOL% ( $\pm 0.20$  VOL%); > 5.00 (unspecified)

**Isoflurane:** 0 to 1.00 VOL% ( $\pm 0.15$  VOL%); 1.00 to 5.00 VOL% ( $\pm 0.20$  VOL%); > 5.00 (unspecified)

**Desflurane:** 0 to 1.00 VOL% ( $\pm 0.15$  VOL%); 1.00 to 5.00 VOL% ( $\pm 0.20$  VOL%); 5.00 to 10.00 VOL% ( $\pm 0.40$  VOL%); 10.00 to 15.00 VOL% ( $\pm 0.60$  VOL%); 15.00 to 18.00 VOL% ( $\pm 1.00$  VOL%); > 18.00 (unspecified)

**Enflurane:** 0 to 1.00 VOL% ( $\pm 0.15$  VOL%); 1.00 to 5.00 VOL% ( $\pm 0.20$  VOL%); > 5.00 (Unspecified)

**Sevoflurane:** 0 to 1.00 VOL% ( $\pm 0.15$  VOL%); 1.00 to 5.00 VOL% ( $\pm 0.20$  VOL%); 5.00 to 8.00 VOL% ( $\pm 0.40$  VOL%); > 8.00 (Unspecified)

**Carbon dioxide:** 0 to 1.00 VOL% ( $\pm 0.10$  VOL%); 1.00 to 5.00 VOL% ( $\pm 0.20$  VOL%); 5.00 to 7.00 VOL% ( $\pm 0.30$  VOL%); 7.00 to 10.00 VOL% ( $\pm 0.50$  VOL%); > 10.00 (Unspecified)

**Nitrous oxide:** 0 to 20 VOL% ( $\pm 2.00$  VOL%); 20.0 to 100 VOL% ( $\pm 3.00$  VOL%)

#### Interference Gas

**CO<sub>2</sub>:** N<sub>2</sub>O, O<sub>2</sub>, any agent = 0.1%<sub>ABS</sub> inaccuracy allowance for each

**N<sub>2</sub>O:** CO<sub>2</sub>, O<sub>2</sub>, any agent = 0.1%<sub>ABS</sub> inaccuracy allowance for each

**Agents:** CO<sub>2</sub> = 0%<sub>ABS</sub> inaccuracy allowance; N<sub>2</sub>O, O<sub>2</sub>, second agent = 0.1%<sub>ABS</sub> inaccuracy allowance for each

#### Flow Rate

**Adult and pediatric:** 200  $\pm 20$  ml per min

**Neonate:** 150  $\pm 15$  ml per min

**Maximum specified interval for intervention of water** (hours at specified minimum sample flow rate): **AGENT mode** - Adult and pediatric is 17 hours @ 200 ml/min, 37°C, 100% RH;

neonate is 17 hours @ 120 ml/min, 37°C, 100% RH; **CO<sub>2</sub> mode**  
- 8 hours @ 50 ml/min ±10 ml/min

### System Response and Rise Times

*As measured from patient gas input of the complete pneumatic circuit, including tubing, from 10 to 90 percent of measured levels*

#### Adult Cannula:

|                 | System Response (sec) | Rise Time (sec) |
|-----------------|-----------------------|-----------------|
| Halothane       | 11.56                 | 5.77            |
| Isoflurane      | 6.71                  | 0.88            |
| Desflurane      | 6.33                  | 0.57            |
| Enflurane       | 7.55                  | 1.75            |
| Sevoflurane     | 6.45                  | 0.62            |
| CO <sub>2</sub> | 6.62                  | 0.61            |
| Oxygen          | 6.99                  | 1.02            |
| Nitrous Oxide   | 6.28                  | 0.25            |

#### Infant Cannula:

|                 | System Response (sec) | Rise Time (sec) |
|-----------------|-----------------------|-----------------|
| Halothane       | 15.95                 | 8.63            |
| Isoflurane      | 9.26                  | 1.70            |
| Desflurane      | 6.47                  | 0.61            |
| Enflurane       | 11.98                 | 4.75            |
| Sevoflurane     | 6.48                  | 0.62            |
| CO <sub>2</sub> | 6.51                  | 0.48            |
| Oxygen          | 8.61                  | 1.13            |
| Nitrous Oxide   | 7.95                  | 0.72            |

#### Adult Divided Cannula:

|                 | System Response (sec) | Rise Time (sec) |
|-----------------|-----------------------|-----------------|
| Halothane       | 20.81                 | 14.18           |
| Isoflurane      | 10.99                 | 3.91            |
| Desflurane      | 7.38                  | 0.64            |
| Enflurane       | 13.83                 | 7.11            |
| Sevoflurane     | 7.48                  | 0.78            |
| CO <sub>2</sub> | 7.57                  | 0.64            |
| Oxygen          | 8.02                  | 1.07            |
| Nitrous Oxide   | 7.16                  | 0.51            |

#### Infant Divided Cannula:

|                 | System Response (sec) | Rise Time (sec) |
|-----------------|-----------------------|-----------------|
| Halothane       | 9.98                  | 3.95            |
| Isoflurane      | 6.75                  | 0.89            |
| Desflurane      | 6.25                  | 0.60            |
| Enflurane       | 7.32                  | 1.37            |
| Sevoflurane     | 5.45                  | 0.67            |
| CO <sub>2</sub> | 5.49                  | 0.49            |
| Oxygen          | 7.25                  | 0.84            |
| Nitrous Oxide   | 6.51                  | 0.39            |

**Data sample rate:** 25 Hz

**Full accuracy respiration rate** (range permitting specified gas accuracy): 2 to 60 rpm

**Total respiration range:** 2 to 100 rpm; accuracy is unspecified from 60 to 100 rpm

**Relevant interference:** 0.5 mmHg equivalent with 37.5°C saturated with H<sub>2</sub>O (0.1 % relative max)

**Display resolution:** 0.1 % volume

**Maximum warm-up time:** 10 minutes; ISO accuracy achieved in less than 45 seconds of activation

**Auto ID threshold** (full accuracy mode): 0.15 % (Primary agent ID); 0.3 % (Secondary agent ID)

**Multiple agents alarm threshold:** 0.3 % (0.5 % during ISO accuracy mode) or 5%<sub>REL</sub> (10 % for isoflurane) of primary

agent if primary agent > 10 % (For halothane add 0.1 %<sub>ABS</sub> to threshold values)

**CO<sub>2</sub> ambient pressure compensation range:** 500 to 900 mmHg

**Pressure compensation:** Unaffected by cyclical pressures of up to 10 kPa as, apart from the described automatic pressure compensation, the pump automatically regulates flow so that not only gas readings but also gas sample flow is unaffected

**Calibration interval:** Calibration verification (as described in service instructions) must be performed at one year intervals

### Alarm Limits

|                    | Lower*                          | Upper*                           |
|--------------------|---------------------------------|----------------------------------|
| EtCO <sub>2</sub>  | 5 to 60 mmHg;<br>0.6 to 8.0 kPa | 5 to 76 mmHg;<br>0.7 to 10.1 kPa |
| FiCO <sub>2</sub>  | No limit                        | 0 to 20 mmHg; 0 to 2.7 kPa       |
| FiN <sub>2</sub> O | No limit                        | 0 to 80 %                        |
| Et Halothane       | 0.1 to 5.0 VOL%                 | 0.1 to 5.0 VOL%                  |
| Fi Halothane       |                                 |                                  |
| Et Isoflurane      |                                 |                                  |
| Fi Isoflurane      |                                 |                                  |
| Et Desflurane      | 0.1 to 18.0 VOL%                | 0.1 to 18.0 VOL%                 |
| Fi Desflurane      | 0.1 to 5.0 VOL%                 | 0.1 to 5.0 VOL%                  |
| Et Enflurane       |                                 |                                  |
| Fi Enflurane       | 0.1 to 8.0 VOL%                 | 0.1 to 8.0 VOL%                  |
| Et Sevoflurane     |                                 |                                  |
| Fi Sevoflurane     | 18 to 100 %                     | 20 to 100 %                      |
| FiO <sub>2</sub>   |                                 |                                  |

\* Off is an option for each alarm limit

### CO<sub>2</sub>

**Range:** 0 to 76 mmHg (0 to 10.1 kPa)

**Resolution:** 1 mmHg (0.1 kPa)

### O<sub>2</sub>

**Range:** 0 to 100 %

**Resolution:** 1 %

**Signal Output** (at constant temperature and pressure): 10 mV ±1.5 mV @ 20° C / 20.95 % O<sub>2</sub>

**Maximum response time** (21 to 100 % step change through patient sampling line as seen in WPU gas monitor window): < 7.3 seconds (adult and pediatric); < 8.2 seconds (neonate)

**Accuracy** (includes stability and drift, full scale gas measurement performance met after the maximum warm-up period): ±1 % (0 to 40 %); ±2 % (40 to 60 %); ±3 % (60 to 80 %); ±4 % (80 to 100 %)

**Offset:** ±1 %

### O<sub>2</sub> interfering gas effects

**N<sub>2</sub>O:** < 0.3 VOL% @ 80 VOL% N<sub>2</sub>O

**CO<sub>2</sub>:** < 0.3 VOL% @ 5 VOL% CO<sub>2</sub>

**Halothane:** < 0.3 VOL% @ 5 VOL% halothane

**Enflurane:** < 0.3 VOL% @ 5 VOL% enflurane

**Isoflurane:** < 0.3 VOL% @ 5 VOL% isoflurane

**Desflurane:** < 0.3 VOL% @ 18 VOL% desflurane

**Sevoflurane:** < 0.3 VOL% @ 8 VOL% sevoflurane

**Acetone:** < 0.3 VOL% @ 1 VOL% acetone

**Ethanol:** < 0.3 VOL% @ 0.1 VOL% ethanol

**Helium:** < 0.3 VOL% @ 80 VOL% helium

**Methane:** < 0.3 VOL% @ 0.1 VOL% methane

**Nitric oxide:** < 0.3 VOL% @ 50 ppm nitric oxide



## Oxygen Sensor

**Operating temperature:** 15 to 35 °C (59 to 95°F)

**Expected operating life:** Use by the expiration date printed on sensor and its packaging

## Respiration (Pneumatic)

*Displayed numerically by detecting the patient's abdominal or chest wall motion through a pneumatic bellows placed at the patient's chest. No user adjustable options, including alarms, as this parameter is not intended for vital sign monitoring.*

**Respiration rate measurement range:** 0 to 60 rpm

**Respiration rate resolution:** 1 rpm

**Respiration rate accuracy:**  $\pm 1$  rpm

## Temperature (Optional)

*For use with the FlexTEMP II Sensor*

**Channel:** One

**Units:** Celsius and Fahrenheit

**Range:** 20.0 to 44.0°C (68.0 to 111.2°F)

**Resolution:** 0.1°C (0.1°F)

**Accuracy:**  $\pm 0.5^\circ\text{C}$  ( $\pm 0.9^\circ\text{F}$ )

**Response time:** The measuring time to obtain a steady-state reading within the manufacturer's accuracy specifications is within 15 seconds, compliant to ISO 80601-2-56.

**Numeric display update time:** 2 seconds

**Sensor type:** Fiber-optic, multiple-use (when used with single-use sterilized jackets)

**Application site:** Axillary, esophageal, rectal

**Measurement mode:** Direct

## Alarm Limits

| Lower*                             | Upper*                             |
|------------------------------------|------------------------------------|
| 20.0 to 44.0°C;<br>68.0 to 111.2°F | 20.0 to 44.0°C;<br>68.0 to 111.2°F |

## Non-invasive Blood Pressure

*Oscillometric method (with inflatable cuff) determines systolic, diastolic and mean arterial pressures, and pulse rate.*

**Patient Types:** Adult, pediatric, and neonate

**Unit of measure:** Millimeters of mercury (mmHg) or kilopascals\*\* (kPa)

**Cuff inflation pressure:** Initially 165 mmHg (22 kPa) for Adult, 130 mmHg (17.3 kPa) for Pediatric, and 100 mmHg (13.3 kPa) for Neonate; all pressures are  $\pm 15$  mmHg (2 kPa); Subsequent inflation pressures determined by the previous NIBP measurement

**Overpressure protection:** release of cuff pressure if inflation pressure exceeds 300 mmHg (40 kPa) for Adult and Pediatric modes, and 150 mmHg (20 kPa) for Neonate mode

## Measurement Range

|           | Adult                           | Pediatric                       | Neonate                          |
|-----------|---------------------------------|---------------------------------|----------------------------------|
| Systolic  | 30 to 270 mmHg;<br>4 to 36 kPa  | 30 to 180 mmHg;<br>4 to 24 kPa  | 30 to 130 mmHg;<br>4 to 17.3 kPa |
| Diastolic | 10 to 245 mmHg; 1.3 to 32.7 kPa | 10 to 150 mmHg; 1.3 to 20 kPa   | 10 to 100 mmHg; 1.3 to 13.3 kPa  |
| Mean      | 20 to 255 mmHg; 2.7 to 34 kPa   | 20 to 160 mmHg; 2.7 to 21.3 kPa | 20 to 120 mmHg; 2.7 to 16 kPa    |

## Accuracy

**Pressure measurement accuracy:** Maximum mean error  $\pm 5$  mmHg ( $\pm 0.6$  kPa) with a standard deviation of less than 8 mmHg (1 kPa)

**Pressure measurement resolution:** 1 mmHg (0.1 kPa)

**Pressure transducer range:** 0 to 300 mmHg (0 to 40 kPa)

## Modes

**Manual:** Immediate upon operator command

**Automatic:** Determinations automatically made with selectable intervals of 1, 2, 3, 5, 10, 15, 20, and 30 minutes

## Alarm Limits

|           |           | Lower*                             | Upper*                             |
|-----------|-----------|------------------------------------|------------------------------------|
| Systolic  | Adult     | 30 to 270 mmHg;<br>4.0 to 36.0 kPa | 30 to 270 mmHg;<br>4.0 to 36.0 kPa |
|           | Pediatric | 30 to 180 mmHg;<br>4.0 to 24.0 kPa | 30 to 180 mmHg;<br>4.0 to 24.0 kPa |
|           | Neonate   | 30 to 130 mmHg;<br>4.0 to 17.3 kPa | 30 to 130 mmHg;<br>4.0 to 17.3 kPa |
| Diastolic | Adult     | 10 to 245 mmHg;<br>1.3 to 32.7 kPa | 10 to 245 mmHg;<br>1.3 to 32.7 kPa |
|           | Pediatric | 10 to 150 mmHg;<br>1.3 to 20.0 kPa | 10 to 150 mmHg;<br>1.3 to 20.0 kPa |
|           | Neonate   | 10 to 100 mmHg;<br>1.3 to 13.3 kPa | 10 to 100 mmHg;<br>1.3 to 13.3 kPa |
| Mean      | Adult     | 20 to 255 mmHg;<br>2.7 to 34.0 kPa | 20 to 255 mmHg;<br>2.7 to 34.0 kPa |
|           | Pediatric | 20 to 160 mmHg;<br>2.7 to 21.3 kPa | 20 to 160 mmHg;<br>2.7 to 21.3 kPa |
|           | Neonate   | 20 to 120 mmHg;<br>2.7 to 16.0 kPa | 20 to 120 mmHg;<br>2.7 to 16.0 kPa |

\* Off is an option for each alarm limit

\*\*For kilopascals (kPa), allow  $\pm 1$  least significant digit to accommodate round-off error for calculated values.

## Gating

*Parameter result outputs to the MRI system as data and discrete signals.*

### Digital pulses (parameter event-associated signals):

- ECG (3.3 to 5.0 V p-p signal, pulse duration 10 ms  $\pm 3$  ms)
- SpO<sub>2</sub> (3.3 to 5.0 V p-p signal, pulse duration 10 ms  $\pm 3$  ms)
- Negative pulses (-3.3 to -5.0 V p-p signals), other characteristics same as above

**Analog waveforms** (monitored parameter representative signals):

- ECG (1 mV/mV scaling, 5 mA maximum current, 20 mV maximum output voltage)
- ECG (1 V/mV scaling,  $\pm 5$  V maximum output voltage, 5 mA maximum current)
- IBP (200 mV maximum output voltage)
- Respiration ( $\pm 5$  V maximum output voltage, 5 mA maximum current, 1 V p-p signal voltage)
- SpO<sub>2</sub> IR/red (1 V/mV scaling, 40 mV maximum output voltage)
- SpO<sub>2</sub> IR/red (2 V maximum output voltage)

## Accessories

(Original part number in parenthesis where applicable.)

### AGENT

989803152561 (9012) - Cannula, Disp, Adult  
989803152601 (9016) - Cannula, Disp, Adult  
989803152621 (9016B) - Cannula, Disp, INT INF, (Divided)  
989803152631 (9016C) - Cannula, Disp, PED, (Divided)  
989803152611 (9016A) - Cannula, Disp, Infant, (Divided)  
989803152591 (9015) - Cannula, Disp, INT Infant  
989803152571 (9013) - Cannula, Disp, PED  
989803152581 (9014) - Cannula, Disp, Infant  
989803162051: Anesthetic Oxygen (O<sub>2</sub>) Sensor  
989803152671 (94012) – Kit, Disposable, Water Trap, 3160  
989803152661 (94018): Kit, Sample, Agents, 3160

### CO2 Cannulas

989803183241 - LoFlo Sample Line, Adult Cannula, Box 20  
989803183251 - LoFlo Sample Line, Ped. Cannula, Box 20  
989803183261 - LoFlo Sample Line, Neo. Cannula, Box 20  
989803183271 - LoFlo Line, Adu Dvd Cannula, Box 20  
989803183281 - LoFlo Line, Ped Dvd Cannula, Box 20  
989803183291 - LoFlo Line, Adu Airway Adpt, Box 20  
989803185331 - LoFlo Sample Line, Adult Cannula, Box 100  
989803185341 - LoFlo Sample Line, Ped Cannula, Box 100  
989803185351 - LoFlo Sample Line, Neo Cannula, Box 100  
989803185361 - LoFlo Line, Adu Dvd Cannula, Box 100  
989803185371 - LoFlo Line, Ped Dvd Cannula, Box 100  
989803185381 - LoFlo Line Adu Airway Adpt, Box 100

### ECG

989803152291 (9009) – Gel, ECG/EEG Skin Prep, Tube, 3-pack  
989803193721 - Expression MR ECG Leads, AAMI, CV  
989803193731 - Expression MR ECG Leads, AAMI, Standard  
989803193741 - Expression MR ECG Leads, AAMI, Neonatal  
989803193751 - Expression MR ECG Leads, IEC, CV  
989803193761 - Expression MR ECG Leads, IEC, Standard  
989803193771 - Expression MR ECG Leads, IEC, Neonatal  
989803179031 - QUADTRODE MRI ECG Pad, 25/box  
989803179041 - ELCTRD, MRI ECG, QUTRD.CV, 25/box  
989803179051 - ELCTRD, MRI, NEO.QUDTRD, 25/box  
989803192761 - Wireless ECG Patient Module (GEN 3) 1-5  
989803194341 - Wireless ECG Patient Module (GEN 3) 6-10

### Gating

989803152821 (9292) - CAB, Digital Gating, GE, 3160  
989803152831 (9291) - CAB, Gating, Siemens, 3160  
989803152851 (9293) - CAB, Dig. Gating, HIT/TOSH, 3160  
989803195521 – Universal Gating Interface

### Invasive Blood Pressure

989803194601 - Expression MR IBP Transducer Cable, 5 Ft  
989803194631 - Expression MR IBP DPT Kit, A/P, Box 20  
989803194641 - Expression MR IBP DPT Kit, I/N, Box 20

### Non-invasive Blood Pressure (NIBP) Cuffs

989803182611 - NIBP Cuff, Single Lumen, Infant  
989803182621 - NIBP Cuff, Single Lumen, Pediatric  
989803182631 - NIBP Cuff, Single Lumen, Small Adult  
989803182641 - NIBP Cuff, Single Lumen, Adult  
989803182651 - NIBP Cuff, Single Lumen, Adult-L  
989803182661 - NIBP Cuff, Single Lumen, LRG Adult

989803182671 - NIBP Cuff, Single Lumen, LRG Adult-L  
989803182681 - NIBP Cuff, Single Lumen, Thigh

### Disposable Non-invasive Blood Pressure (NIBP) Cuffs

989803182511 - NIBP Cuff, Single Lumen, Infant, Disp  
989803182521 - NIBP Cuff, Single Lumen, Pediatric, Disp  
989803182531 - NIBP Cuff, Single Lumen, Small Adult, Disp  
989803182541 - NIBP Cuff, Single Lumen, Adult, Disp  
989803182551 - NIBP Cuff, Single Lumen, Adult-L, Disp  
989803182561 - NIBP Cuff, Single Lumen, LRG Adult, Disp  
989803182571 - NIBP Cuff, Single Lumen, LRG Adult-L, Disp  
989803182581 - NIBP Cuff, Single Lumen, Thigh, Disp  
989803183171 - NIBP Cuff, Single Lumen, Neo #1, Disp  
989803183181 - NIBP Cuff, Single Lumen, Neo #2, Disp  
989803183191 - NIBP Cuff, Single Lumen, Neo #3, Disp  
989803183201 - NIBP Cuff, Single Lumen, Neo #4, Disp  
989803183211 - NIBP Cuff, Single Lumen, Infant #5, Disp

### Non-invasive Blood Pressure (NIBP) Hoses

989803183221 - Adult Pressure Interconnect Hose  
989803183231 - Neonatal Pressure Interconnect Hose

### Respiration (Pneumatic)

989803152791 (94023) - Pneumograph, Chest, NM, 3160

### SpO<sub>2</sub>

989803161991 - Quick Connect SpO<sub>2</sub> Probe, MRI  
989803166531 - Quick Connect SpO<sub>2</sub> Clip, Adult  
989803166541 - Quick Connect SpO<sub>2</sub> Clip, Pediatric  
989803166551 - Quick Connect SpO<sub>2</sub> Grip, Adult, 20/box  
989803166561 - Quick Connect SpO<sub>2</sub> Grip, Ped, 20/box  
989803166571 - Quick Connect SpO<sub>2</sub> Grip, Infant, 20/box  
989803166581 - Quick Connect SpO<sub>2</sub> Grip, Neo, 20/box  
989803192771 - Wireless SpO<sub>2</sub> Patient Module (Gen 3) 1-5  
989803194331 - Wireless SpO<sub>2</sub> Patient Module (Gen 3) 6-10

### System

989803191341 - Battery, Module (Gen 3)  
989803169491 - Battery, MRI, 14.8V, 5.08 AH, UL  
865471 - Expression Information Portal - IP5  
989803176521 - Advanced Communications Option  
453564177501 - European Line Cord  
989803168211 - North American Line Cord  
989803168221 - Cord, Jumper, 25 Feet  
989803173901 - Brazilian Power Cord, 3 meter  
989803174171 - UK Line Cord, 3 meter  
989803181291 - Power Cord, AUS/NZL, 3 meter  
989803181321 - Power Cord, S Africa, 3 meter  
989803181331 - Power Cord, Danish, 3 meter  
989803181341 - Power Cord, Israeli, 3 meter  
989803181351 - Power Cord, Argentina, 3 meter  
989803181361 - Power Cord, Swiss, 3 meter

### Temperature

989803194511 - FlexTEMP II Sensor  
(Esophageal/Rectal/Axillary, Direct Mode)  
989803168891 - Surgical Lubricant, 12 pack  
989803178181 - FlexTEMP System, Jacket (Box 10)



## **User Guides**

989803193191 - Danish\*  
989803193201 - Dutch\*  
989803193211 - English\*  
989803193221 - Finnish  
989803193231 - French\*  
989803193241 - German\*  
989803193251 - Indonesian  
989803193261 - Italian\*  
989803193271 - Japanese  
989803193281 - Korean  
989803193291 - Norwegian\*  
989803193301 - Polish  
989803193311 - Brazilian Portuguese\*  
989803193321 - Russian  
989803193331 - Spanish\*  
989803193341 - Swedish\*  
989803193351 - Traditional Chinese  
989803193361 - Turkish  
989803196981 - Czech  
989803198661 - Bulgarian  
989803198671 - Hungarian  
989803198681 - Latvian  
989803198691 - Lithuanian  
989803198701 - Romanian  
989803198711 - Slovakian  
989803198821 - Estonian

\* - User interface software also translated in these languages

## **Miscellaneous**

989803196881 - Quick Reference Guide  
989803195211 - Service Manual



For more information about the Expression Model MR400 MRI Patient Monitoring System or any of our complete solution products, please contact us. We are glad to hear from you.



**Invivo, a division of  
Philips Medical Systems**  
12151 Research Parkway  
Orlando, FL 32826, USA

**Phone**  
**USA:** 800-722-9377

**E-mail**  
Worldwide: [Info@invivocorp.com](mailto:Info@invivocorp.com)

**Websites**  
[www.ExpressionMR.com](http://www.ExpressionMR.com)  
[www.invivocorp.com](http://www.invivocorp.com)  
[www.philips.com](http://www.philips.com)

© 2016 Koninklijke Philips N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

REF A-866185-90058 Rev. C

