LBM[®] CAT Broth

Instructions for Use

CAT BROTH

Information Leaflet & User Guide

See the glossary of symbols at the end of the information leaflet.

INTENDED USE

CAT BROTH is a device for preservation and cultivation of *Trichomonas vaginalis*, *Candida albicans* and other yeasts present in urogenital samples.

SUMMARY AND PRINCIPLES

Among the pathogens responsible for infectious genital disease, those determining the onset of vaginitis are essentially bacteria, fungi and protozoa. Mycotic vaginitis is one of the most frequent genital infections caused by *Candida albicans* in about 70-80% of cases and by non albicans species for the remaining 20-30%. An increase in the yeast load - normally present in the vagina even if healthy - due to microbiological and/or physiological factors increases the probability of infection symptoms setting in, especially when pseudohyphae, the most virulent form, manifest themselves. *Trichomonas vaginalis*, a sexually transmitted flagellate protozoan, is the cause of infection of the genito-urinary tract, which may be asymptomatic but can also manifest itself with characteristic symptomatology and can extend to the uterus and the tubes, at times causing severe pelvic infection and in pregnant women an increased risk of premature birth. The symptoms are less evident in men but can cause urethritis and prostatitis. In any event, it would appear that it increases susceptibility to HIV. The most common and fastest method used to diagnose trichomoniasis is direct microscopic examination of vaginal secretion, while the most sensitive technique is examination of vaginal and urethral sample cultures. Using a culture medium that allows preserving and growing *Trichomonas* and *Candida* is a useful tool for microbiological diagnostic examination to study urethrovaginal infections in order to ascertain the microbial aetiology of the disease.

The protein base constituting the medium provides the source of carbon and amine nitrogen added to which is horse serum as growth supplement and yeast extract as source of vitamin B complex. The medium contains a wide-spectrum antibiotic supplement (chloramphenicol) as inhibiting agent for concurrent flora with an inhibiting action towards gram-negative and gram-positive bacteria while the presence of reducing agents favours the anaerobic growth conditions.

REAGENTS

CAT BROTH medium composition

Pancreatic Digest of Casein
Yeast extract
Glucose
Sodium chloride
L-Cysteine
Sodium thioglycollate
Horse serum
Chloramphenicol

WARNINGS

1. CAT BROTH is a single-use device intended for in vitro diagnostic use.
2. Observe the approved precautions for potentially hazardous biological materials and aseptic techniques.
3. This device may only be used by adequately trained and qualified personnel. The manufacturer declines all responsibility for use by unqualified or unauthorised persons.
4. The clinical samples may contain pathogenic microorganisms, including hepatitis or human immunodeficiency viruses. Handle all materials and items contaminated with blood and other biological fluids in compliance with the guidelines and "standard precautions" dictated by the institutional bodies.
5. All samples and materials used for product analysis are to be considered potentially infected and must therefore be handled in such a way as to prevent infection of laboratory personnel. After use, sterilize all potentially hazardous waste including samples, containers and transport and culture media.
6. Strictly follow the instructions.
7. The device is not suitable for transport of microorganisms not included in the intended use.
8. Do not ingest the culture medium.
9. Do not use the culture medium to pre-moisten the swab before collecting the sample or rinse or irrigate the sampling sites.
10. The product is not reusable; reuse may result in incorrect diagnosis.
11. Do not use the device if there are evident signs of damage.
12. Avoid skin contact with the culture medium.
13. If the sample was taken and transported using common sampling and transport systems and is to be processed with automatic plating systems, remove the swab from the test tube containing the enrichment medium before processing.
14. Consult the operating instructions supplied with the device or available in electronic format, and which can be identified by the e-IFU indicator on the packaging label.

STORAGE

CAT BROTH is ready to use and does not require additional preparation. It must be stored in its original package at 5-25°C until the time of use. Do not overheat. Do not incubate or freeze prior to use; inappropriate storage determines a loss of effectiveness. Do not use after the expiry date printed on the external package and on the label of the culture medium test tube.

PRODUCT DETERIORATION

Do not use CAT BROTH if it shows visible signs of damage or contamination, leakage or other signs of deterioration.

MATERIAL INCLUDED

Code	Description	Availability	Use in automation
094CE.A	12x80mm tapered-bottom test tube with screw cap containing 2 ml of CAT BROTH medium	1 box of 300 pcs (6x50 units)	Yes
094CE	12x80mm tapered-bottom test tube with screw cap containing 2 ml of CAT BROTH medium	1 box of 300 pcs (6x50 units)	No
092CN	16x100mm tapered-bottom test tube with screw cap containing 4 ml of CAT BROTH medium and fiber swab with breakpoint in a pouch.	1 box of 300 pcs (6x50 units)	Yes

REQUIRED MATERIAL BUT NOT INCLUDED

Ancillary culture plates, reagents, organisms for quality control, various laboratory instruments, Regular Flocked Swab with pre-moulded breakpoint at 80 mm (502CS01), Minitip Flocked Swab with pre-moulded breakpoint at 80 mm (501CS01).

NB: If the swab is processed with automatic plating systems for culture samples, the Regular Flocked Swab with pre-moulded breakpoint at 80 mm (502CS01) is validated for use. Any other type of swab – including the Minitip Flocked Swab (code 501CS01) – must be removed from the tube before loading the sample on the automatic system.

INSTRUCTIONS FOR USE

Sample collection

The samples to be analysed to search for *Trichomonas vaginalis*, *Candida albicans* and other yeasts may be in various forms: vaginal exudate, prostatic secretion, urethral secretion.

- Loosen the cap of the test tube containing the culture medium.
- Collect the sample.
- Open the test tube and immediately insert the swab.
- Break the swab shaft inside the test tube at the breakpoint indicated.
- Securely screw the cap back on.
- Identify the test tube with the patient's data and send the sample to the laboratory for analysis.

Use in the laboratory

Culture examination for *T. vaginalis*: incubate the test tube at 35±2°C for 2-7 days and starting from the second day, observe the typical morphology under the microscope.

Culture examination for yeasts: the sample can directly be seeded and/or after test tube incubation at 35±2°C for 24 hours on a suitable culture medium plate (i.e. SDA plate) by means of:

- Automatic plating
- Direct seeding with swab

QUALITY CONTROL AND PERFORMANCE CHARACTERISTICS

The quality control procedures include checking system performance using pure microbial suspensions of ATCC® strains of the two target microbiological agents, *Trichomonas vaginalis* and *Candida albicans*, and checking the inhibiting capacity using *Escherichia coli*.

Procedures applied:

Growth test with *C. albicans*

1. From a pure culture of an ATCC® strain, prepare a 0.5 McF suspension corresponding to 1.5x10⁸ CFU/ml.
2. Prepare a dilution from the 0.5 McF bacterial suspension able to give a plate count of 30-300 CFU in 100 µl.
3. Retrieve 200 µl of medium from the thus prepared dilution and inoculate the CAT BROTH medium.
4. Close the system and vortex for 10 seconds at 2000/2500 RPM.
5. Plate 100 µl at time zero (immediately) on an appropriate culture medium plate (i.e. SDA plates).
6. Incubate the zero-time test tube and plates at 35 ±2°C for 24 hours.
7. After 24 hours incubation, retrieve the CAT BROTH device from the incubator.
8. Seed the plate (i.e. SDA plate) with 100 µl using a micropipette and a spreader.
9. Incubate the plates for a further 24 hours at 35 ±2°C.

Inhibition test with *E. coli*

1. From a pure culture of an ATCC® strain, prepare a 0.5 McF suspension corresponding to 1.5x10⁸ CFU/ml.
2. Prepare a dilution from the 0.5 McF bacterial suspension able to give a semiconfluent growth after seeding the plate with 100 µl (0.5 suspension diluted 1:100).
3. Retrieve 200 µl of medium from the thus prepared dilution and inoculate the CAT BROTH medium.
4. Close the system and vortex for 10 seconds at 2000/2500 RPM.
5. Plate 100 µl at time zero (immediately) on an appropriate culture medium plate (i.e. SDA plates).
6. Incubate the zero-time test tube and plates at 35 ±2°C for 24 hours.
7. After 24 hours incubation, retrieve the CAT BROTH device from the incubator.
8. Seed the plate (i.e. SDA plate) with 100 µl using a micropipette and a spreader.
9. Incubate the plates for a further 24 hours at 35 ±2°C.

Growth test with *T. vaginalis*

1. Immerse the bulb until the frozen material is completely covered in a bain-marie at 35±2°C for 2-3 minutes until it has fully thawed.
2. Immediately after thawing, aseptically transfer the strain to a test tube containing 13ml of reconstituted broth 2154 brought up to pH 6.0 by adding 150 µl of 1N HCl. Incubate the test tube at 35±2°C in a 15° slanted position.
3. Empirically determine the culture density peak by microscopic inspection evaluating the number of vital organisms (normally 3-4 days).
4. When the density peak has been reached, place the test tube in ice for 10 minutes. Gently turn the test tube upside down 10 times and aseptically transfer 200 µl of the suspension to a CAT BROTH test tube and immediately do a zero-time count by fresh microscopic examination.
5. Incubate the test tube at 35°C for 72 hours.
6. After 72 hours, retrieve 100 µl and check for the presence/number of protozoa by fresh microscopic examination.

RESULTS

ATCC® strain	Average of CFUs recovered at time 0 hr	Average of CFUs recovered at time 24 hrs
<i>C. albicans</i> ATCC® 10231	186	Confluent growth
<i>E. coli</i> ATCC® 25922	Semiconfluent growth	Inhibition (from partial to total)

ATCC® strain	Average count at 0 hr in 100 µl	Average count at 72 hrs in 100µl
<i>T. vaginalis</i> ATCC® 30001	155	> 155

The results of the tests were obtained with ATCC® strains in the laboratory. The results are expressed as the means obtained from three test-tubes at 24 hours for *C. albicans* and *E. coli*, and at 72 hours for *T. vaginalis*.

The results obtained largely depend on appropriate and adequate collection of the sample, as well as on prompt transport and laboratory analysis.

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BIBLIOGRAPHY

1. Krieger JN and Alderete JF. *Trichomonas vaginalis* and trichomoniasis. In: K. Holmes, P. Markh, P. Sparling et al (eds). *Sexually Transmitted Diseases*, 3rd Edition. New York: McGraw-Hill, 1999, 587-604.
2. Centers for Disease Control and Prevention. Sexually transmitted diseases treatment guidelines 2002. *MMWR. Morbidity and Mortality Weekly Report* 2002;51 (No.RR-6)
3. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. *Manual of Clinical Microbiology*. 7th ed. ASM, Washington, D.C..
4. National Committee for Clinical Laboratory Standards (NCCLS). 2004. *Quality Control for commercially Prepared Microbiological Culture Media*. Approved Standard – Third Edition M22-A3
5. Miller, J. M. 1999. *A Guide to Specimen Management in Clinical Microbiology*, 2nd ed. ASM, Washington, DC.
6. Isenberg, H. D., 2004. *Clinical Microbiology Procedures Handbook*, 2nd ed. ASM, Washington, DC.
7. Atlas, R.N.1993 *Handbook of Microbiological Media* p. 147-153. CRC Press, Boca Raton, FL

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Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com

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Copan Italia S.p.A.