

## Platinum-E Retroviral Packaging Cell Line, Ecotropic

**CATALOG NUMBER:** RV-101

**STORAGE:** Liquid nitrogen

*Note: For best results begin culture of cells immediately upon receipt. If this is not possible, store at -80°C until first culture. Store subsequent cultured cells long term in liquid nitrogen.*


**QUANTITY & CONCENTRATION:** 1.0 mL,  $>3 \times 10^6$  cells/mL in DMEM, 20% FBS, and 10% DMSO

### Background

Retroviruses are efficient tools for delivering heritable genes into the genome of dividing cells. However, conventional NIH-3T3 based retroviral packaging cell lines have limited stability and produce low viral yields, mainly due to poor expression level of the retroviral structure proteins (gag, pol, env) in the packaging cells.

The Platinum-E (Plat-E) Cell Line, a potent retrovirus packaging cell line based on the 293T cell line, was generated using novel packaging constructs with an EF1 $\alpha$  promoter to ensure longer stability and high-yield retroviral structure protein expression (gag, pol, ecotropic env). Plat-E cells can be kept in good condition for at least 4 months in the presence of drug selection, and can produce retroviruses with an average titer of  $1 \times 10^7$  infectious units/mL by transient transfection. In addition, replication competent retroviruses (RCR) are virtually nonexistent because only coding sequences of viral structural genes are used, avoiding any unnecessary retroviral sequences.

The Plat-E cell line is designed for rapid, transient production of high-titer ecotropic retrovirus.



**Figure 1.** NIH3T3 cells (left) and mouse ProB Ba/F3 cells (right) were infected with GFP retrovirus supernatant produced in Plat-E cells after transfection with pMX-GFP.

## **Medium**

1. Culture Medium: DMEM, 10% fetal bovine serum (FBS), 1 µg/mL puromycin, 10 µg/mL blasticidin, penicillin and streptomycin
2. Freeze Medium: 70% DMEM, 20% FBS, 10% DMSO

## **Methods**

### **I. Establishing Plat-E Cultures from Frozen Cells**

1. After quickly thawing the cells in a 37°C water bath, immediately transfer the thawed cell suspension into a 15 mL tube containing 10 mL of culture medium.
2. Centrifuge the tube for 5 min at 1300 to 1500 rpm.
3. Discard the supernatant and break the cell pellet by finger tapping.
4. Add a few drops of culture medium with gentle shaking and finger tap the tube a few times.
5. Add 2 mL of culture medium to the tube and gently pipet the cell suspension up and down twice.
6. Transfer the cell suspension to a 10 cm culture dish (Falcon® #3003 works well) containing 8 mL of culture medium.
7. Swirl the culture plate well to mix the cells, then incubate the cells for three days before expansion.

#### ***Important Notes:***

- ***Don't change the culture medium during the first three days. It is normal to see some cells floating after the first 24 hours.***
- ***Don't culture cells to complete confluency. Split cells 4X to 6X every two to three days when the culture reaches 70-90% confluency.***

### **II. Splitting the Cells**

*Note: Avoid forming bubbles as much as possible during this procedure.*

1. Wash cells once with PBS.
2. Add 4 mL of 0.05% Trypsin/0.5 mM EDTA solution to a 10 cm dish and incubate at 37°C for 3-5 min.
3. Remove the cells from the dish surface by tapping the rim of the culture dish.
4. Transfer 10 mL of the culture medium to a 50 mL tube.
5. Using the same pipette with some residual culture medium, wash the dish surface gently three times in 4 mL of the Trypsin/EDTA solution.
6. Gently pipette the cell suspension up and down 7 times and transfer the cell suspension into the 50 mL tube containing 10 mL medium from step 4.
7. Centrifuge the cells for 5 min at 1300-1500 rpm.
8. Discard the supernatant and break the cell pellet by finger tapping.
9. Add a few drops of culture medium with gentle shaking and finger tap the tube a few times.
10. Add 5 mL of culture medium and gently pipet the cell suspension up and down twice.
11. Add 15 mL of culture medium, then count and seed the cells. Typically  $10^7$  cells can be harvested from one 10 cm culture dish.

## **Transfection**

1. Seed  $2 \times 10^6$  cells in a 60 mm culture dish without antibiotics including puromycin and blasticidin one day before transfection.
2. After 16 to 24 hours, start transfection when the culture becomes 70-80% confluent.  
*Note: We suggest transfecting cells with FuGENE® Transfection Reagent (Promega) or Lipofectamine™ Plus (Invitrogen). For example, 3 µg retroviral expression plasmid is mixed with 9 µL FuGENE® Transfection Reagent according to the manufacturer's recommendation. The mixed DNA- FuGENE® complex is added by dropwise into the culture media.*
3. Harvest retroviral supernatant 48 hours after transfection.

## **References**

1. Morita, S., Kojim, T., and Kitamura, T. (2000) *Gene Therapy* **7**: 1063-1066.
2. Takahashi, K. and Yamanaka, S. (2006) *Cell* **126**: 663-676.

## **Recent Product Citations**

1. Torcal Garcia, G. et al. (2023). Carm1-arginine methylation of the transcription factor C/EBPα regulates transdifferentiation velocity. *Elife*. **12**:e83951. doi: 10.7554/eLife.83951.
2. Zhang, Y. et al. (2023). Net39 protects muscle nuclei from mechanical stress during the pathogenesis of Emery-Dreifuss muscular dystrophy. *J Clin Invest*. **133**(13):e163333. doi: 10.1172/JCI163333.
3. Pachmayr, L.O. et al. (2023). Unbiased chemokine receptor screening reveals similar efficacy of lymph node- and tumor-targeted T cell immunotherapy. *Cancer Immunol Immunother*. doi: 10.1007/s00262-023-03472-w.
4. Sato, S. et al. (2023). The circadian clock CRY1 regulates pluripotent stem cell identity and somatic cell reprogramming. *Cell Rep*. **42**(6):112590. doi: 10.1016/j.celrep.2023.112590.
5. Salemme, V. et al. (2023). p140Cap inhibits β-Catenin in the breast cancer stem cell compartment instructing a protective anti-tumor immune response. *Nat Commun*. **14**(1):2350. doi: 10.1038/s41467-023-37824-y.
6. Tabata, K. et al. (2023). Monitoring and assessment of lysosomal membrane damage in cultured cells using the high-content imager. *STAR Protoc*. **4**(2):102236. doi: 10.1016/j.xpro.2023.102236.
7. Jin, J. et al. (2023). CISH impairs lysosomal function in activated T cells resulting in mitochondrial DNA release and inflammaging. *Nat Aging*. **3**(5):600-616. doi: 10.1038/s43587-023-00399-w.
8. Bhatia, V. et al. (2023). Targeting advanced prostate cancer with STEAP1 chimeric antigen receptor T cell and tumor-localized IL-12 immunotherapy. *Nat Commun*. **14**(1):2041. doi: 10.1038/s41467-023-37874-2.
9. Read, K.A. et al. (2023). Aiolos represses CD4+ T cell cytotoxic programming via reciprocal regulation of TFH transcription factors and IL-2 sensitivity. *Nat Commun*. **14**(1):1652. doi: 10.1038/s41467-023-37420-0.
10. Pham, D. et al. (2023). Batf stabilizes Th17 cell development via impaired Stat5 recruitment of Ets1-Runx1 complexes. *EMBO J*. **42**(8):e109803. doi: 10.15252/embj.2021109803.
11. Hahn, A.M. et al. (2023). A monoclonal Trd chain supports the development of the complete set of functional γδ T cell lineages. *Cell Rep*. **42**(3):112253. doi: 10.1016/j.celrep.2023.112253.
12. Liu, Q. et al. (2023). Tcf21 marks visceral adipose mesenchymal progenitors and functions as a rate-limiting factor during visceral adipose tissue development. *Cell Rep*. **42**(3):112166. doi: 10.1016/j.celrep.2023.112166.

13. Diril, M. et al. (2023). Genetic dissection of the Mastl-Arpp19/Ensa-PP2A-B55 $\delta$  pathway in mammalian cells. *TJB*. **48**(2):190-202. doi: 10.1515/tjb-2022-0191.
14. Tang, J. et al. (2023). Runx3-overexpression cooperates with ex vivo AKT inhibition to generate receptor-engineered T cells with better persistence, tumor-residency, and antitumor ability. *J Immunother Cancer*. **11**(2):e006119. doi: 10.1136/jitc-2022-006119.
15. Oh S. et al. (2023). Precision targeting of autoantigen-specific B cells in muscle-specific tyrosine kinase myasthenia gravis with chimeric autoantibody receptor T cells. *Nat Biotechnol*. doi: 10.1038/s41587-022-01637-z.
16. Briukhovetska, D. et al. (2023). T cell-derived interleukin-22 drives the expression of CD155 by cancer cells to suppress NK cell function and promote metastasis. *Immunity*. **56**(1):143-161.e11. doi: 10.1016/j.immuni.2022.12.010.
17. Saltukoglu, D. et al. (2023). Plasma membrane topography governs the 3D dynamic localization of IgM B cell antigen receptor clusters. *EMBO J*. **42**(4):e112030. doi: 10.15252/embj.2022112030.
18. Caravia, X.M. et al. (2022). Loss of function of the nuclear envelope protein LEMD2 causes DNA damage-dependent cardiomyopathy. *J Clin Invest*. **132**(22):e158897. doi: 10.1172/JCI158897.
19. Zenke, S. et al. (2022). Differential trafficking of ligands trogocytosed via CD28 versus CTLA4 promotes collective cellular control of co-stimulation. *Nat Commun*. **13**(1):6459. doi: 10.1038/s41467-022-34156-1.
20. Jain, P. et al. (2022). Discovery and functional characterization of the oncogenicity and targetability of a novel NOTCH1-ROS1 gene fusion in pediatric angiosarcoma. *Cold Spring Harb Mol Case Stud*. **8**(6):a006222. doi: 10.1101/mcs.a006222.
21. Zhong, X. et al. (2022). Decoupling the role of ROR $\gamma$ t in the differentiation and effector function of TH17 cells. *Sci Adv*. **8**(42):eadc9221. doi: 10.1126/sciadv.adc9221.
22. Bresser, K. et al. (2022). Replicative history marks transcriptional and functional disparity in the CD8 $^{+}$  T cell memory pool. *Nat Immunol*. doi: 10.1038/s41590-022-01171-9.
23. Sloat, S.R. & Hoppins, S. (2022). A dominant negative mitofusin causes mitochondrial perinuclear clusters because of aberrant tethering. *Life Sci Alliance*. **6**(1):e202101305. doi: 10.26508/lsa.202101305.
24. Hall, J.A. et al. (2022). Transcription factor ROR $\alpha$  enforces stability of the Th17 cell effector program by binding to a Rorc cis-regulatory element. *Immunity*. doi: 10.1016/j.immuni.2022.09.013.
25. Que, F. et al. (2022). RHOA G17V induces T follicular helper cell specification and involves angioimmunoblastic T-cell lymphoma via upregulating the expression of PON2 through an NF- $\kappa$ B-dependent mechanism. *Oncoimmunology*. **11**(1):2134536. doi: 10.1080/2162402X.2022.2134536.
26. Ren, Y. et al. (2022). Tumorous expression of NAC1 restrains antitumor immunity through the LDHA-mediated immune evasion. *J Immunother Cancer*. **10**(9):e004856. doi: 10.1136/jitc-2022-004856.
27. Hojo, H. et al. (2022). Runx2 regulates chromatin accessibility to direct the osteoblast program at neonatal stages. *Cell Rep*. **40**(10):111315. doi: 10.1016/j.celrep.2022.111315.
28. Koschade, S.E. et al. (2022). Translatome proteomics identifies autophagy as a resistance mechanism to on-target FLT3 inhibitors in acute myeloid leukemia. *Leukemia*. doi: 10.1038/s41375-022-01678-y.
29. Kotov, J.A. et al. (2022). LT $\beta$ R overexpression promotes plasma cell accumulation. *PLoS One*. **17**(8):e0270907. doi: 10.1371/journal.pone.0270907.

30. Ma, S. et al. (2022). Protocol to assess cell-intrinsic regulatory mechanisms using an ex vivo murine T cell polarization and co-culture system. *STAR Protoc.* **3**(3):101543. doi: 10.1016/j.xpro.2022.101543.

### **Warranty**

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

### **License Information**

This licensed product is intended for ACADEMIC, GOVERNMENT AND NON-PROFIT RESEARCH USE ONLY; not for use in diagnostic or therapeutic procedures. The product may not be transferred, sold, or otherwise provided to another laboratory except by an authorized distributor of Cell Biolabs, Inc.

Use of this product by Biotechnology and Pharmaceutical companies requires a license for all fields of use including research. Please contact:

Director of Business Development  
Cell Biolabs, Inc.  
[busdev@cellbiolabs.com](mailto:busdev@cellbiolabs.com)

### **Contact Information**

Cell Biolabs, Inc.  
7758 Arjons Drive  
San Diego, CA 92126  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2008-2023: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.