rainbow® Series Patient Cables, EMS ^{EN} Pour Pulse Oximeters with rainbow Technology

DIRECTIONS FOR USE	MODE D'EMPLOI
Reusable  Not made with natural rubber latex  Non-Sterile	

INDICATIONS

The rainbow® Series Patient Cables, EMS have the same indications for use as the applicable sensors. Please see the compatible sensor's directions for use for indications for use and prescribed information.

DESCRIPTION

The rainbow Series Patient Cables, EMS are used with Masimo® rainbow and M-LNCS® sensors. rainbow Series Patient Cables, EMS are for use only with devices containing Masimo rainbow SET® technology Version 7.1 or higher or licensed to use rainbow compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular sensor and device models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor model.

rainbow Series Patient Cables, EMS were verified using Masimo rainbow SET technology.

WARNING: Masimo sensors and cables are designed for use with devices containing Masimo SET® or Masimo rainbow SET oximetry or licensed to use Masimo sensors.

WARNINGS, CAUTIONS, AND NOTES

- Always refer to the oximeter module operator's manual for complete instructions or additional instructions.
- Ensure the cable is physically intact, with no broken or frayed wires or damaged parts. Visually inspect the cable and discard if cracks or discolored areas are present.
- All sensors and cables are designed for use with specific monitors. Verify the compatibility of the monitor and cable sensor before use, otherwise degraded performance and/or patient injury can result.
- Carefully follow patient cable to reduce the possibility of patient entanglement or strangulation.
- Failure to properly connect the sensor or the oximeter module to the cable will result in interrupted readings, inaccurate results, or no reading.
- To avoid damage to the cable, always hold it by the connector rather than the cable when connecting or disconnecting either it.
- To prevent damage, do not soak or immerse the sensor or cable in any liquid solution. Do not attempt to sterilize the cable or sensor.
- Do not attempt to reprocess, recondition or recycle Masimo sensors or patient cables as these processes may damage the electrical components, potentially leading to patient harm.
- Caution:** Replace the cable when a replace cable message is displayed, or when a low SIQ message is consistently displayed while monitoring consecutive patients, after completing the low SIQ troubleshooting steps identified in the monitoring device operator's manual.
- Note:** The cable is provided with X-Cal® technology to minimize the risk of inaccurate blood flow monitoring. The cable must be replaced every 17,520 hours of patient monitoring time. Replace the cable when the patient monitoring time is exhausted. Some legacy monitoring devices do not distinguish between the cable and the sensor. In these devices, a replace sensor or similar message could pertain either to the sensor or the cable.

CABLE AND SENSOR SELECTION

NOTE: The maximum length for the sensor, cable and extension cable combination is 15 feet. Refer to the table below for correct examples:

Cable	Sensor	
Length	18 inches	3 feet
1 ft	Yes	Yes
4 ft	Yes	No
12 ft	No	No

INSTRUCTIONS

A. Connect the rainbow patient cable, EMS to the Device

- Connect the device connector of the patient cable to the mating connector of the device.
- Ensure the patient cable connector is fully engaged with the device. A tactile or audible click of connection should be heard.

B. Connect the rainbow patient cable, EMS to the Sensors

Connect the rainbow patient cable, EMS to an M-LNCS disposable sensor:

- Refer to Fig. 1. Orient the sensor connector to the patient cable connector as shown.
- Refer to Fig. 2. Insert the sensor connector completely into the patient cable connector until it locks in place as shown.

Connect the rainbow patient cable, EMS to a reusable M-LNCS or a reusable rainbow sensor:

- Refer to Fig. 3. Orient the sensor connector (1) to patient cable connector (2) as shown.
- Refer to Fig. 4. Insert the sensor connector (1) completely into the patient cable connector (2) until it locks in place as shown.

Connect the rainbow patient cable, EMS to a reusable M-LNCS or a reusable rainbow sensor:

- Refer to Fig. 5. Orient sensor connector (1) to patient cable connector (2) as shown.
- Refer to Fig. 6. Insert the sensor connector (1) completely into the patient cable connector (2) until it locks in place.

C. Disconnecting the rainbow patient cable, EMS from the Sensors

NOTE: Pull on the connectors only, do not pull the cables.

Disconnect the M-LNCS disposable sensor:

- Refer to Fig. 2. Pull fully on the sensor connector to remove it from the patient cable.

Disconnect the rainbow disposable sensor:

- Refer to Fig. 4. While holding the sides of cable connector, place thumb on latch (1) and press to disconnect sensor connector from patient cable connector.

Disconnect the reusable M-LNCS or reusable rainbow sensor:

- Refer to Fig. 5. While holding the sides of patient cable connector (2), lift sensor connector cover (3) from patient cable connector (2). Disconnect sensor connector (1) from patient cable connector (2).

CLEANING

- Remove the sensor from the patient and disconnect it from the patient cable.
- Clean the surface of the cable by wiping it with a 70% isopropyl alcohol pad.
- Wipe all surfaces of the cable.
- Saturate another cloth or gauze pad with sterile or distilled water and wipe all surfaces of the cable.

- Dry the cable by wiping all surfaces with a clean cloth or dry gauze pad.

CAUTION:

- Do not immerse the connectors on the cable in any liquid solution.
- Do not sterilize by irradiation, steam, autoclave, or ethylene oxide.
- Do not clean with chemicals not approved above.

ENVIRONMENTALS

Operating Temperature 41°F to 104°F (5°C to 40°C)

Storage Temperature -40°F to 158°F (-40°C to 70°C)

Relative Humidity 10% to 95% non-condensing

COMPATIBILITY

rainbow series patient cables, EMS are intended for use only with devices containing Masimo SET® or Masimo rainbow SET® technology and with compatible sensors. Consult individual oximetry system manufacturers for compatibility of particular devices and sensor models. Each device manufacturer is responsible for determining whether their devices are compatible with each sensor and cable model. Use of this cable with other devices may result in no or improper performance.

For Compatibility Information Reference: www.Masimo.com

WARRANTY

Masimo warrants to the initial buyer only that this product, when used in accordance with the directions provided with the Products by Masimo, will be free of defects in materials and workmanship for a period of 24 months (2 years). Single use products are warranted for single patient use only.

THE FOREGOING IS THE SOLE AND EXCLUSIVE WARRANTY APPLICABLE TO THE PRODUCTS SOLD BY MASIMO TO BUYER. MASIMO EXPRESSLY DISCLAIMS ALL OTHER ALLO, EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. MASIMO'S SOLE OBLIGATION AND BUYER'S EXCLUSIVE REMEDY FOR BREACH OF ANY WARRANTY SHALL BE, AT MASIMO'S OPTION, TO REPAIR OR REPLACE THE PRODUCT.

WARRANTY EXCLUSIONS

This warranty does not extend to any product that has been used in violation of the operating instructions supplied with the product, or has been subject to misuse, neglect, accident or externally created damage. This warranty does not extend to any product that has been connected to any unintended device or system, has been modified, or has been disassembled or reassembled. This warranty does not extend to sensors or patient cables that have been reprocessed, reconditioned or recycled.

IN NO EVENT SHALL MASIMO BE LIABLE TO BUYER OR ANY OTHER PERSON FOR ANY INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES INCLUDING WITHOUT LIMITATION LOST PROFITS, EVEN IF ADVISED OF THE POSSIBILITY THEREOF. IN NO EVENT SHALL MASIMO'S LIABILITY ARISING FROM ANY PRODUCTS SOLD TO BUYER UNDER A CONTRACT, WARRANTY, TORT OR OTHER CLAIM EXCEED THE AMOUNT PAID BY BUYER FOR THE LOT OF PRODUCTS INVOLVED IN SUCH CLAIM. IN NO EVENT SHALL MASIMO BE LIABLE FOR ANY DAMAGES ASSOCIATED WITH A PRODUCT THAT HAS BEEN REPROCESSED, RECONDITIONED OR RECYCLED. THE LIMITATIONS IN THIS SECTION SHALL NOT BE DEEMED TO PRECLUDE ANY LIABILITY THAT, UNDER APPLICABLE PRODUCTS LIABILITY LAW, CANNOT LEGALLY BE PRECLUDED BY CONTRACT.

NO IMPLIED LICENSE

Purchase or possession of this rainbow patient cable does not carry any express or implied license to use this cable with any device that is not an authorized device or separately authorized to use rainbow patient cables.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

If you encounter any serious incident with product, please notify the competent authority in your country and the manufacturer.

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION	
	Follow instructions for use		Consult instructions for use		CE	European Union Conformity Mark
	Manufacturer		Lot		REF	Catalogue number (model number)
	Date of manufacture YYMM-YY		Lot		REF	Maxima reference number
	Caution		Storage handling instruction		LOT	Storage temperature range
	Do not discard		Greater than		REF	Less than
	Fragile, handle with care		Not made with natural rubber latex		REF	Atmospheric pressure limitation
	Non-sterile		Unique device identifier		REF	Instructions/Directions for Use (Manuals are available in electronic format at http://www.Masimo.com/Products . Note: eUDI is not available in all countries.
	Medical device		UDI		REF	Importer
	Importer		Distributor		REF	Distributor

Patents: <http://www.masimo.com/patents.htm>

Masimo rainbow SET, M-LNCS, X-Cal, SET,  and rainbow are the federally registered trademarks of Masimo Corporation.

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