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Isı Sensörlü Foley Catheter 400TM

COVIDIEN™

Covidien, **COVIDIEN™**, the Covidien logo are trademarks of Covidien AG. Other brands are trademarks of a Covidien Company.

Covidien, **COVIDIEN™**, the Covidien logo are trademarks of Covidien AG. Other brands are trademarks of a Covidien Company.
© 2013 Covidien

DEGANIA SILICONE LTD. DEGANIA BET 1513000, ISRAEL

EC REP: BIOMETRIX S.R.O. VICENZY 16, SAMORIN 93101 SLOVAKIA
GERMAN DISTRIBUTOR: COVIDIEN DEUTSCHLAND GmbH, Gewerbepark 1, D-93333 NEUSTADT/DONAU, GERMANY
SWISS DISTRIBUTOR: COVIDIEN SWITZERLAND LTD., ROOSSTRASSE 53, CH-8832 WOLLERAU, SWITZERLAND

ENGLISH

Sterile: Unless packaging has been opened or damaged. Discard after single use, do not re-use. Ethylene Oxide gas sterilized Exposure to elevated temperatures and ultraviolet light should be avoided during storage

This instruction insert applies to the following Covidien product:

SIZE	8Ch	10Ch	12Ch	14Ch	16Ch	18Ch
REF	90053T	90054T	90055T	90056T	90051T	90052T

The temperature is taken by the 400 series thermistor probe which has been inserted into the third channel up to the radiopaque tip of the catheter.

The temperature range is from 25°C. up to 45°C and with a tolerance of ± 0,1°C. Minimum measuring time to obtain required tolerance is 150 seconds.

DIRECTION FOR USE OF BALLOON CATHETER WITH TEMPERATURE SENSOR The device must be connected to a CE Marked medical device for which this catheter is an authorized accessory! The monitoring system must comply with EN 60601-1.

DESCRIPTION: The Mon-a-Therm™ Foley Catheter with Temperature Sensor is packaged individually as a sterile, single-use device and is available in a variety of sizes. The catheter size and balloon volume are also designed on the catheter itself. The catheter is of silicone construction. The balloon inflation system is terminated with a Luer syringe actuated, self-sealing valve. A temperature sensor is located near the distal tip. Each Mon-a-Therm™ Foley Catheter with Temperature Sensor is designed for electrical connection to a specific reusable cable available from Covidien.

INDICATIONS: The Mon-a-Therm™ Foley Catheter with Temperature Sensor is indicated for use in the drainage of the urinary bladder and the simultaneous monitoring of temperature. Foley (balloon) catheterization is usually indicated when continuous bladder drainage is required in the presence of obstruction or paralysis, in the absence of voluntary micturition or after surgery or trauma involving contiguous pelvic organs, and is an aid in evaluating urine output in patients who may have fluid management problems.

CONTRAINDICATIONS: There are no known contraindications.

WARNING: Do not stretch the silicone catheter shaft as this can cause damage/fracture of the temperature sensor.

PRECAUTIONS:
♦ The recommended balloon inflation volume (3ml or 5ml) is printed on the product. In the case of a 5ml balloon, the balloon will withstand some over inflation. However the absolute maximum acceptable inflation will be reached with 10ml of water.
A 3ml balloon must not be inflated beyond the recommended 3ml volume.
♦ Expert clinical judgment should be exercised in the selection of the appropriate size catheter.
♦ Prior to insertion, lubricate the catheter with a suitable water soluble lubricant.
♦ Observe aseptic technique during the insertion process.
♦ Deflate the balloon prior to attempting to remove the catheter.
♦ Puncture of the funnel connector wall could result in leakage. Therefore, aspirating urine through the funnel connector wall with needles should be avoided.
♦ Never use a syringe to puncture the catheter shaft for urine sampling. There is a danger of needle stick injury and the catheter may be compromised.
♦ Do not handle with sharp instruments.
♦ Never clamp the catheter. A catheter plug or valve may be used if required.
♦ Do not attempt to force connection of the sensor. The connectors are keyed by shape. Forcing the connection could result in damage to the sensor or cable connector and a resulting in ability of function.
♦ Take care to insure that the cable connector does not lie in an area where solutions may pool during the surgical procedure.
♦ As with any temperature probe, use acceptable technique during electrosurgical procedures to help reduce the possibility that the probe will serve as an alternate path for radio frequency current and burn the surrounding tissues. Proper technique at minimum should include:
a) Properly locate and connect an adequate electrosurgical dispersive electrode close to the active surgical site.
b) Avoid draping the patient temperature monitor instrument cable over grounded surfaces or intertwining it with the cables of the electrosurgical generator.
♦ The operation of some patient temperature monitor instruments may be temporarily affected during electrosurgical activations. Unusual temperature reading should be checked during a pause between activations.
♦ Precautions should be taken when discarding the device and disposal of the device shall be made in accordance with applicable national regulations or biologically hazardous waste.
♦ Do not re-use or resterilize. Re-use may cause a risk of infection and/or compromise functional reliability.

ADVERSE REACTIONS: The following adverse reactions have been reported to be associated with use of Foley catheters or temperature probes especially during prolonged use. The order listing is alphabetical and does not indicate frequency or severity. Reported adverse reactions include; bladder spasms, electrical burns due to aberrant electrosurgical pathways, encrustation with urinary salts (primarily calcium), kidney or bladder stones, septicemia, urethral discharge, urethritis, urethrocutaneous fistulae and urinary tract infection or irritation.

SUGGESTED DIRECTIONS FOR USE: Do not use if packaging as been opened or damaged.
1. Follow currently accepted aseptic technique for patient preparation and catheter care.
2. Carefully remove the sterile Mon-a-Therm™ Foley Catheter with Temperature Sensor from its protective packaging.
3. Lubricate the catheter with a suitable water lubricant.
4. Insert using accepted medical technique.
5. Be sure the catheter is draining properly before inflating balloon.
6. To inflate balloon, insert a Luer tip syringe filled with sterile water into the valve. Depress the plunger. The recommended balloon inflation volume is printed on the device. In the case of a 5ml balloon, the balloon will withstand some over inflation, however the absolute maximum acceptable inflation will be reached with 10ml of water. A 3ml balloon must not be inflated beyond the recommended 3ml volume. The valve will seal when the syringe tip is withdrawn. Attach drainage tubing and closed system bag.
7. Secure the catheter and drainage tubing in place using accepted medical technique. They should be secured without tension on the meatus/catheter junction.
8. Connect the appropriate reusable cable. To do this, insert the sensor connector into the mating receptacle at the end of the cable until it is firmly seated. The connectors are keyed by shape.
9. Follow the DIRECTIONS FOR USE of the patient temperature monitor instrument. If problems are encountered in monitoring temperature, verify that he necessary cable connections have been made and maintained.
13. Discard the device. Consideration should be given to specific PRECAUTIONS stated in the instruction insert.
14. A patient has to be routinely monitored and according to his needs temperature sensor catheter should be changed or removed. Maximum recommended indwelling time per each particular catheter is up to six weeks.

CONDITIONS OF USE, STORAGE, AND TRANSPORTATION:
♦ Minimum environmental operating range is from +10°C to 40°C and 30% to 75% humidity
♦ Storage and transportation range is from -30°C to +60°C and 85% humidity.

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