


Color-coded
trays



eco-
friendly



NEW!

A new design for a new day

The New epT.I.P.S.® Motion Tip Box System

epT.I.P.S. Motion have been upgraded! To assure the highest level of sample safety when pipetting and to make the epMotion even easier to use epT.I.P.S. Motion will now be available with color-coded trays and a new flexible packaging system. As an Eco-friendly option, epT.I.P.S. Motion are also now available as a Reload-System.

The new epT.I.P.S. Motion racks consist of two options for simple integration into the epMotion liquid handling workstation: A single use box with color coded trays for easy volume identification and a TipHolder Adapter to use the trays as a Reload-System.

With color coded trays you can quickly and easily see the volume range of the tip and assure an even higher level of sample safety. The new Reload trays also offer a more environmentally friendly option and will be delivered in a PET blister packaging with a sealed cap assuring the same well established quality as all Eppendorf tips.

The optional TipHolder Adapter replaces the normal tip box component and reduces waste by up to 40%.

eppendorf

Purity grade

Volume µl	Filter	Racks / Reloads	Tips / unit	eppendorf quality	PCR clean	Sterile	PCR clean & sterile
50	-	10 Racks	960	●			
50	with	10 Racks	960		●		
50	-	10 Racks	960			●	
50	with	10 Racks	960				●
50	-	24 Reload Trays	2304	●			
50	with	24 Reload Trays	2304		●		
300	-	10 Racks	960	●			
300	with	10 Racks	960		●		
300	-	10 Racks	960			●	
300	with	10 Racks	960				●
300	-	24 Reload Trays	2304	●			
300	with	24 Reload Trays	2304		●		
1000	-	10 Racks	960	●			
1000	with	10 Racks	960		●		
1000	-	10 Racks	960			●	
1000	with	10 Racks	960				●
1000	-	24 Reload Trays	2304	●			
1000	with	24 Reload Trays	2304		●		

epT.I.P.S.® Motion pipette tips, without filter Order no.

Eppendorf Quality, 10 × 96 tips in racks	
50 µL, volume range 1–50 µL	0030 014.405
300 µL, volume range 20–300 µL	0030 014.448
1000 µL, volume range 40–1000 µL	0030 014.480
Sterile, 10 × 96 tips in racks	
50 µL, volume range 1–50 µL	0030 015.207
300 µL, volume range 20–300 µL	0030 015.223
1000 µL, volume range 40–1000 µL	0030 015.240
Eppendorf Quality (Reloads), 24 × 96 tips	
50 µL, volume range 1–50 µL	0030 014.421
300 µL, volume range 20–300 µL	0030 014.464
1000 µL, volume range 40–1000 µL	0030 014.502

epT.I.P.S.® Motion pipette tips, with filter Order no.

PCR clean, 10 × 96 tips in racks	
50 µL, volume range 1–50 µL	0030 014.413
300 µL, volume range 20–300 µL	0030 014.456
1000 µL, volume range 40–1000 µL	0030 014.499
PCR clean, sterile, 10 × 96 tips in racks	
50 µL, volume range 1–50 µL	0030 015.215
300 µL, volume range 20–300 µL	0030 015.231
1000 µL, volume range 40–1000 µL	0030 015.258
PCR clean (Reloads), 24 × 96 tips	
50 µL, volume range 1–50 µL	0030 014.430
300 µL, volume range 20–300 µL	0030 014.472
1000 µL, volume range 40–1000 µL	0030 014.510

More safety


- Color-coded trays for direct volume identification
- Double O-ring on dispensing tools assures an optimal fit each time
- Optical sensor automatically identified tip type

Flexible handling

- Available in multiple purity grades
- Easily exchangeable, no additional labware files required
- Also available as refill system

Eco-friendly Reload System

- Easy conversion with TipHolder Adapter
- TipHolder Adapter is autoclavable
- Reduces waste by up to 40%



For more information visit: www.eppendorf.com or contact your local Eppendorf sales representative.

eppendorf
In touch with life

Your local distributor: www.eppendorf.com/worldwide

Eppendorf AG Hamburg · Germany · E-mail: eppendorf@eppendorf.com

Application Support Europe: Tel: +49 1803 666 789 (Preis je nach Tarif im Ausland; 9 ct/min aus dem dt. Festnetz; Mobilfunkhöchstpreis 42 ct/min)
support@eppendorf.com

North America: Tel: +1 516 334 7500 · Toll free phone: +1 800 645 3050 · E-mail: techserv@eppendorf.com

Asia Pacific: Tel: +60 3 8023 6869 · E-mail: support_asiapacific@eppendorf.com