

## Bottle Assemblies with Three-Piece Caps

2

Nominal Filling Capacity (mL)	Nominal Size (mm)	Bottle & Cap Assembly	Bottle Only	Cap Assembly	Cap Only	Plug Only	O-ring Only	Rotors
<b>Bottle and Cap Assembly</b>								
PC — 8.5	16 x 64	355615	355656	355604	335257	335256	870409	Type 50
PC — 10.4	16 x 76	355603	355651	355604	335257	335256	870409	Types 90 Ti, 80 Ti, 75 Ti, 70.1 Ti, 65, 50 Ti, 40, MLA-55
PC — 26.3	25 x 89	355616 <sup>1</sup>	355654	355617	335259	335258	870385	Type 30
PC — 26.3	25 x 89	355618 <sup>1</sup>	355654	355619	338824	335258	870385	Types 70 Ti, 60 Ti, 55.2 Ti, 50.2 Ti, 42.1
PC — 70.0	38 x 102	355620 <sup>2</sup>	355655	355621	334547	334545	870384	Types 35, 21
PC — 70.0	38 x 102	355622 <sup>2</sup>	355655	355623	339165	334545	870384	Type 45 Ti
PA —250.0	60 x 120	334205	355627	362247	362246	332836	812715	Type 19
PP —250.0	62 x 122	356011	358326	358977	357490	N/A	927860	Type 16
PC —250.0	62 x 122	356013	358275	358977	357490	N/A	927860	Type 16

Note: Regarding packaging, bottle and cap assemblies are packaged in units of 6, with the following exceptions: Bottle and Cap assembly 334205 is packaged as one; Bottles 355656 and 355651 are packaged in units of 25; Cap assembly 334543 is packaged as one. Caps, Plugs, and O-rings are packaged separately.

<sup>1</sup> The only difference between these two assemblies is the cap which is Noryl in the 355616 assembly and aluminum in the 355618 assembly. The Noryl cap can be used in the Type 30 Rotor.

<sup>2</sup> The only difference between these two assemblies is the cap which is Noryl in the 355620 assembly and aluminum in the 355622 assembly. The Noryl cap can be used in the Type 35 and 21 Rotors; the aluminum cap is required for the higher speed forces of the Type 45 Ti.