




Vannes de commutation Rheodyne

	Ø passage (mm)	P max (bar)	Θ lim. (°C)	A	B	C
Vanne inox 1/16"	0,6	483	150	7000	7030	7040
Vanne inox "large bore" 1/16"	1,0	340	150	7000L	7030L	7040L
Vanne PEEK 1/16"	0,4	340	50	---	9030	---
Vanne inox 1/8"	1,0	340	50	3000-038	3030-038	3030
Vanne PEEK 1/8"	1,0	276	50	3000	---	---

Vannes de sélection Rheodyne (6 entrées - 1 sortie)

	Ø passage (mm)	P max (bar)	Θ lim. (°C)	Réf.
Vanne de sélection 1/16" Inox	0,4	483	80	7060
Vanne de sélection 1/16" Inox "large bore"	1,0	340	80	7060L
Vanne de sélection 1/16" PEEK	0,4	340	40	9060

Boucles d'injection pour vannes Rheodyne

Volume	Inox 1/16" (7725, 7125 et 7010)	PEEK 1/16" (9725, 9125 et 9010)	Titane 1/16" (7125-081 et 7010-087)	Inox (8125)	Inox 1/8" (3725-038)	PEEK 1/8" (3725)
5 µL	7020	7020-P	---	8020	---	---
10 µL	7021	7021-P	7055-021	8021	---	---
20 µL	7022	7022-P	7055-022	8022	---	---
50 µL	7023	7023-P	7055-023	8023	---	---
100 µL	7024	7024-P	7055-024	7024	---	---
200 µL	7025	7025-P	7055-025	7025	---	---
500 µL	7026	7026-P	7055-026	7026	---	---
1 mL	7027	7027-P	7055-027	7027	---	---
2 mL	7028	7028-P	7055-028	---	3065-018	3055-018
5 mL	7029	7029-P	7055-029	---	3065-019	3055-019
10 mL	---	---	---	---	3065-023	3055-023
20 mL	---	---	---	---	3065-025	3055-025

Description	Réf.
Chargeur de boucles inox	7012
Chargeur de boucles PEEK	9012
Adaptateur pour seringue à aiguille en PEEK	9013
Needle port cleaner	7125-054
Suction needle adapter	9125-076
Equerre de fixation	7160
	7160-010
	7160-029
	5060-007

