



Seringues pour HPLC autosamplers, Gastight® CTC PAL®

Marque : HAMILTON



18969583.jpg



Seringues C-Line

Jauge	Volume	Longueur	Ø ext.	Type de pointe	Type	Cdt	Ref UGAP	Délai de livraison	Ref fabricant	Ref VWR
22s	10 µl	51 mm	0,71 mm	3	1701 N CTC Slimline™	1	2613047	30 jours	203194	549-1318
22s	25 µl	51 mm	0,71 mm	3	1702 N CTC Slimline™	1	2613051	30 jours	203274	548-1046
22	100 µl	51 mm	0,71 mm	3	1710 FN CTC Slimline™	1	2613061	30 jours	203235	549-1321
22	250 µl	51 mm	0,72 mm	3	1725 N CTC	1	2613066	30 jours	203079	549-1322
22	500 µl	51 mm	0,72 mm	3	1750 N CTC	1	2613072	30 jours	203349	549-1323
22s	100 µl	51 mm	0,71 mm	3	1710 N CTC Slimline™	1	2613059	30 jours	203077	549-1320
22	1000 µl	51 mm	0,72 mm	3	1001 LTN CTC	1	2613908	30 jours	203081	549-1330
22	2500 µl	51 mm	0,72 mm	3	1002 LTN CTC	1	2613352	30 jours	203083	549-1331
22	5000 µl	51 mm	0,72 mm	3	1005 LTN CTC	1	2613909	30 jours	203085	549-1332

Seringues S-line

Jauge	Volume	Longueur	Ø ext.	Type de pointe	Type	Cdt	Ref UGAP	Délai de livraison	Ref fabricant	Ref VWR
22	50 µl	51 mm	0,71 mm	3	1705 N CTC	1	2613803	30 jours	67450-01	HAMI67450-01
22s	10 µl	51 mm	0,71 mm	3	1701 N CTC	1	2613469	30 jours	67444-01	HAMI67444-01
22s	25 µl	51 mm	0,71 mm	3	1702 N CTC	1	2613483	30 jours	67446-01	HAMI67446-01
22s	100 µl	51 mm	0,71 mm	3	1710 N CTC	1	2613508	30 jours	67452-01	HAMI67452-01