

Applied Biosystems™

MicroAmp™ 96-Well Tray for VeriFlex™ Blocks



Catalog number: 4379983

Related applications: [PCR](#) | [Real Time PCR \(qPCR\)](#)

Description

The Applied Biosystems® MicroAmp® 96-Well Tray for VeriFlex™ Blocks is specifically designed to snugly hold 0.1 mL MicroAmp® Reaction Tubes with Caps, MicroAmp® Fast Reaction Tubes with Caps, and MicroAmp® Fast 8-Tube Strips.

- Engineered for efficient sample processing
- Designed specifically for 0.1 mL MicroAmp® reaction tubes with caps
- Manufactured to high specification, ensuring a quality product

Convenience and Performance

The MicroAmp® 96-Well Tray for VeriFlex™ Blocks allows you to organize your tubes while adding sample and mix. This handy tray reduces spillage and tube mix-ups by securing and ordering your tubes on the lab bench. When the samples are ready, simply put the tray of tubes directly on the block of the StepOne™ Real-Time PCR System or Veriti™ Thermal Cycler. Not only does the tray eliminate the hassle of transferring the tubes, but it also provides the thermal contact necessary for thermal cycling.

For Research Use Only. Not for use in diagnostic procedures.

Specifications

Color:	Black
For Use With (Equipment):	ProFlex™ PCR System, SimpliAmp™ Thermal Cyclers, StepOnePlus™, Veriti Dx Fast Thermal Cyclers, Veriti Dx Thermal Cyclers, Veriti Fast Thermal Cyclers, Veriti Thermal Cyclers
Plate Compatibility:	96-Well Plates
Product Line:	MicroAmp®
Product Size:	10 trays
Shipping Condition:	Room Temperature
Type:	Tube Tray

Contents & storage

Contains 10 MicroAmp® 96-Well Trays for VeriFlex™ Blocks. Store at room temperature.