

Applied Biosystems™

MicroAmp™ 96-Well Tray/Retainer Set for Veriti™ Systems



Catalog number: 4381850

Related applications: [PCR](#) | [Real Time PCR \(qPCR\)](#)

Description

The Applied Biosystems® MicroAmp® 96-Well Tray/Retainer Set for Veriti™ Systems is specifically designed to snugly hold MicroAmp® 8-Tube Strips (0.2 mL) and MicroAmp® Reaction Tubes without Caps (0.2 mL) in a 96-well format for processing on the Veriti™ 96-Well Thermal Cycler.

- Engineered for efficient sample processing
- Designed specifically for MicroAmp® 8-Tubes Strip (0.2 mL) and MicroAmp® Reaction Tubes without Caps (0.2 mL)
- Manufactured to high specification, ensuring a quality product

Convenience and Performance The MicroAmp® 96-Well Tray/Retainer Set for Veriti™ Systems allows you to organize your tubes while adding sample and mix. This handy tray reduces spillage and tube mix-ups by securing and ordering your tubes on the lab bench. When the samples are ready, simply put the tray of tubes directly on the block of the Veriti™ 96-Well Thermal Cycler. Not only does the tray eliminate the hassle of transferring the tubes, but it also provides the thermal contact necessary for thermal cycling.

Note: See user's manual or package insert for limited label license and trademark information. For Research Use Only. Not for use in diagnostics procedures.

For Research Use Only. Not for use in diagnostic procedures.

Specifications

Color:	Blue
For Use With (Equipment):	ProFlex™ PCR System, SimpliAmp™ Thermal Cycler, Veriti Dx Thermal Cycler, Veriti Thermal Cycler
Plate Compatibility:	96-Well Plates
Product Line:	MicroAmp®
Product Size:	10 trays
Shipping Condition:	Room Temperature
Type:	Tube Tray

Contents & storage

Contains 10 MicroAmp® 96-Well Tray Retainer Sets. Store at room temperature.