

Applied Biosystems™

# MicroAmp™ 12-Cap Strip



**Catalog number:** N8010534

**Related applications:** [PCR](#)

## Description

MicroAmp® Cap Strips are designed to fit on MicroAmp® reaction tubes, tube strips, and 96-well plates. They are available in 8-cap and 12-cap strips in clear or an assorted set of red, orange, blue, and green. They have a domed profile and are suitable for end-point PCR. Attach caps using the MicroAmp® Cap Installing Tool.

- Will not interfere with sample reading
- Provides tight seal to minimize evaporation
- Easy to apply and remove

### Minimize Evaporation

Use MicroAmp® Caps Strips to seal the tubes containing your precious samples. This will reduce the possibility of contamination among sample wells and help ensure reliable results.

### Easy to Use

Use the MicroAmp® Cap Installing Tool to push the caps down and form a tight seal. When ready, just peel away the caps using the tabs at the end.

### Colored for Organization

These caps are available in different colors to help you keep track of your tubes during assay preparation.

*For Research Use Only. Not for use in diagnostic procedures.*

## Specifications

Cap Type:	12-Cap Strip
Color:	Assorted Colors, Clear
For Use With (Equipment):	2720 Thermal Cycler, GeneAmp 9700, ProFlex™ PCR System, SimpliAmp™ Thermal Cycler, Veriti Dx Fast Thermal Cycler, Veriti Dx Thermal Cycler, Veriti Fast Thermal Cycler, Veriti Thermal Cycler
Product Line:	MicroAmp®
Product Size:	200 strips
Reaction Speed:	Standard
Shipping Condition:	Room Temperature
Type:	PCR cap strip

## Contents & storage

Includes 200 MicroAmp® 12-Cap Strips, Clear. Store at room temperature.