

Applied Biosystems™

MicroAmp™ 12-Cap Strip, assorted colors



Catalog number: N8010834

Related applications: [PCR](#) | [Real Time PCR \(qPCR\)](#)

Description

MicroAmp® Cap Strips are designed to fit on MicroAmp® reaction tubes, tube strips, and 96-well plates. They are available in 8-cap and 12-cap strips in clear or an assorted set of red, orange, blue, and green. They have a domed profile and are suitable for end-point PCR. Attach caps using the MicroAmp® Cap Installing Tool.

Will not interfere with sample reading

- Provides tight seal to minimize evaporation
- Easy to apply and remove

Minimizes Evaporation

Use MicroAmp® Caps Strips to seal the tubes containing your precious samples. This will reduce the possibility of contamination among sample wells and help ensure reliable results.

Easy to Use

Use the MicroAmp® Cap Installing Tool to push the caps down and form a tight seal. When ready, just peel away the caps using the tabs at the end.

Colored for Easy Organization

MicroAmp® Caps Strips are available in different colors to help you keep track of your tubes during assay preparation.

Note: See user's manual or package insert for limited label license and trademark information. For Research Use Only. Not for use in diagnostics procedures.

For Research Use Only. Not for use in diagnostic procedures.

Specifications

| | |
|------------------------------|---|
| Cap Type: | 12-Cap Strip |
| Color: | Assorted Colors, Clear |
| For Use With (Equipment): | 2720 Thermal Cycler, GeneAmp 9700, ProFlex™ PCR System, SimpliAmp™ Thermal Cycler, StepOne™, Veriti Dx Fast Thermal Cycler, Veriti Dx Thermal Cycler, Veriti Fast Thermal Cycler, Veriti Thermal Cycler |
| Product Line: | MicroAmp® |
| Product Size: | 200 strips |
| Shipping Condition: | Room Temperature |
| Type: | PCR cap strip |

Contents & storage

Contains 200 MicroAmp® 12-Cap Strips, Assorted Colors. Store at room temperature.