

Applied Biosystems™

MicroAmp™ 8-Cap Strip, clear



Catalog number: N8011535

Related applications: [PCR](#)

Description

MicroAmp® Cap Strips are designed to fit on MicroAmp® reaction tubes, tube strips, and 96-well plates. They are available in 8-cap and 12-cap strips in clear or an assorted set of red, orange, blue, and green. They have a domed profile and are suitable for end-point PCR. For real-time PCR applications, please see our MicroAmp® Optical 8-Cap Strips (cat. no. 4323032), which have a flat-top profile. Attach the caps using the MicroAmp® Cap Installing Tool.

Features:

- Will not interfere with sample reading
- Provide tight seal to minimize evaporation
- Are easy to apply and remove

Minimize Evaporation

Use MicroAmp® Caps Strips to seal the tubes containing your precious samples. This will reduce the possibility of contamination among sample wells and help ensure reliable results.

Easy to Use

Use the MicroAmp® Cap Installing Tool to push the caps down and form a tight seal. When ready, just peel away the caps using the tabs at the end.

Colored for Organization

These caps are available in different colors to help you keep track of your tubes during assay preparation.

For Research Use Only. Not for use in diagnostic procedures.

Specifications

Cap Type:	8-Cap Strip
Color:	Clear
For Use With (Equipment):	2720 Thermal Cycler, GeneAmp 9700, ProFlex™ PCR System, SimpliAmp™ Thermal Cycler, Veriti Dx Fast Thermal Cycler, Veriti Dx Thermal Cycler, Veriti Fast Thermal Cycler, Veriti Thermal Cycler
Product Line:	MicroAmp®
Product Size:	1,500 strips
Type:	PCR cap strip

Contents & storage

Includes: 5 packs of N8010535 (Qty 300 per pack)

Strip of 8 attached caps for use with 0.2 mL individual tubes, strip tubes (0.1 or 0.2 mL), or 48- or 96-well plates. Natural (clear) color.