

Applied Biosystems™

## MicroAmp™ Reaction Tube with Cap, 0.2 mL



**Catalog number:** N8011540

**Related applications:** [PCR](#) | [Real Time PCR \(qPCR\)](#)

### Description

The Applied Biosystems® MicroAmp® Reaction Tubes with Caps are optimally designed for precise PCR and are designed for use with Applied Biosystems 0.2 mL thermal cyclers including the Veriti™ 96-Well Thermal Cycler, GeneAmp® PCR System 9700 and the Applied Biosystems 2720 Thermal Cycler.

#### Features:

- Consistent plastic thickness enables uniform heating and cooling for accurate results.
- Captive lid with positive-click closure assures proper seating.
- Polished surface and conical bottom allows maximum sample recovery.
- Available in economical bulk packages.
- Option of autoclaved tubes provide for a clean, controlled start

#### Optimized Design

MicroAmp® Reaction Tubes provide unmatched temperature accuracy and uniformity for fast, oil-free PCR amplification with our 0.2mL 96-Well thermal cyclers. Like all of our products, MicroAmp® Reaction Tubes undergo stringent quality testing to ensure reproducible results. They ensure the greatest temperature and amplification uniformity for your samples.

#### Convenient

MicroAmp® Reaction Tubes with Caps are available as clear or colored and in bulk packs. Use MicroAmp® Trays for easy sample handling and optimized results on our 0.2 mL 96-well thermal cyclers.

*For Research Use Only. Not for use in diagnostic procedures.*

## Specifications

Cap Type:	Domed
Color:	Clear
For Use With (Equipment):	2720 Thermal Cycler, GeneAmp 9700, ProFlex™ PCR System, SimpliAmp™ Thermal Cycler, Veriti Dx Thermal Cycler, Veriti Thermal Cycler
Product Line:	MicroAmp®
Product Size:	10,000 tubes
Reaction Speed:	Fast
Tube Type:	Individual
Volume (Metric):	0.2 mL

## Contents & storage

Contains 10,000 MicroAmp® Reaction Tubes with Caps, 0.2 ml (10 packs of SKU# N8010540, Qty 1,000 per pack). Store at room temperature.