



21103AST Supelco

Astec® CYCLOBOND® I 2000 RSP Chiral HPLC Guard

Cartridge, 5 µm particle size, L × I.D. 2 cm × 4 mm

eCl@ss 32110501

FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
21103AST	Seulement 4 en stock (d'autres en cours d'arrivage) - A PARTIR DE	204.00	0

Commandes Bulk?

AJOUTER AU PANIER

Propriétés

Related Categories	Analytical/Chromatography , Astec CYCLOBOND Chiral HPLC Columns (Native and Derivatized Cyclodextrins) , CYCLOBOND I 2000 RSP , Chiral HPLC Columns , Chiral HPLC Guard Columns and Accessories , Plus...
description	HPLC guard column
packaging	pkg of 1 ea
application(s)	HPLC: suitable
L × I.D.	2 cm × 4 mm
matrix	silica particle platform
matrix active group	cyclodextrin, beta- phase
particle size	5 µm
suitability	suitable for L45 per USP
mode of chromatography	chiral

montrer moins de résultats

Description

Recommended products
Discover LiChropur reagents ideal for **HPLC** or **LC-MS** analysis

Legal Information
Astec is a registered trademark of Sigma-Aldrich Co. LLC
CYCLOBOND is a registered trademark of Sigma-Aldrich Co. LLC

Informations Sécurité

Aucune information relative à la sécurité de ce produit est actuellement disponible.

Documents

Certificat d'Analyse

Entrez le Lot N°

[Devis/ Commande de produits Bulk](#)

Customers Also Viewed

21013AST Astec® CYCLOBOND® I 2000 RSP Chiral HPLC Guard Column, 5 µm particle size, L × I.D. 2 cm × 1 mm	20324AST Astec® CYCLOBOND® I 2000 RSP Chiral HPLC Column 5 µm particle size, L × I.D. 25 cm × 4.6 mm	 21150AST Supelguard™ Guard Cartridge Holder Stand-Alone (Swivel-type), for use with Supelguard cartridges (2 cm L. x 2 to 4.6 mm I.D.)	20024AST Astec® CYCLOBOND® I 2000 Chiral HPLC Column 5 µm particle size, L × I.D. 25 cm × 4.6 mm	50179-U Nucleosil® C18 HPLC Column 3 µm particle size, L × I.D. 7.5 cm × 4.6 mm
---	---	--	---	--

Recently Viewed

21100AST Astec® CYCLOBOND® I 2000 Chiral HPLC Guard Cartridge, 5 µm particle size, L × I.D. 2 cm × 4 mm	21066-U Purge/Trap K Vocab® 3000 configured for Tekmar Velocity XPT, LSC-1, LSC-2, LSC-2000, 4000 (Straight, 12 in./30.5 cm, Swagelok fitting), HT-3	21062-U Purge/Trap D configured for Tekmar Velocity XPT, LSC-1, LSC-2, LSC-2000, 4000 (Straight, 12 in./30.5 cm, Swagelok fitting), HT-3	21061-U Purge/Trap C configured for Tekmar Velocity XPT, LSC-1, LSC-2, LSC-2000, 4000 (Straight, 12 in./30.5 cm, Swagelok fitting), HT-3	21060-U Purge/Trap B configured for Tekmar Velocity XPT, LSC-1, LSC-2, LSC-2000, 4000 (Straight, 12 in./30.5 cm, Swagelok fitting), HT-3
--	---	---	---	---

Protocoles et articles

Related Content

[Supelco HPLC Calculator for Isocratic or Gradient Method Transfer](#)
Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...
Keywords: High performance liquid chromatography

Documentation référencée

Did you use this product in your Paper? If so click here.

[read abstract](#)

[Separation of oxazepam, lorazepam, and temazepam enantiomers by HPLC on a derivatized cyclodextrin-bonded phase: application to the determination of oxazepam in plasma.](#)
Chuong Pham-Huy et. al
Journal of biochemical and biophysical methods, 54(1-3), undefined (2003-1-25)
The enantioselective high-performance liquid chromatography (HPLC) of three racemic 3-hydroxybenzodiazepines, oxazepam (Oxa), lorazepam (Lor), and temazepam (Tem), is a difficult operation because of the spontaneous chiral inversion in polar solvent....[Read More](#)

[read abstract](#)

[Chiral stationary phases in HPLC for the stereoselective determination of methadone.](#)
S Rudaz and J L Veuthey
Chirality, 11(4), undefined (1999-5-29)
An extensive study of the behavior of three chiral stationary phases (CSP) is presented for the stereoselective determination of methadone. The following chromatographic columns were selected: a cellulose, Chiralcel OJ; a modified cyclodextrin. Cyclo...[Read More](#)

[read abstract](#)

[Enantiomeric separations of illicit drugs and controlled substances using cyclofructan-based \(LARIHC\) and cyclobond I 2000 RSP HPLC chiral stationary phases.](#)
Nilusha L T Padivitage et. al
Drug testing and analysis, 6(6), undefined (2013-10-12)
Recently a novel class of chiral stationary phases (CSPs) based on cyclofructan (CF) has been developed. Cyclofructans are cyclic oligosaccharides that possess a crown ether core and pendent fructofuranose moieties. Herein, we evaluate the applicabil...[Read More](#)

[read abstract](#)

[\(+\)-Catechin is more bioavailable than \(-\)-catechin: relevance to the bioavailability of catechin from cocoa.](#)
Jennifer L Donovan et. al
Free radical research, 40(10), undefined (2006-10-4)
Catechin is a flavonoid present in fruits, wine and cocoa products. Most foods contain the (+)-enantiomer of catechin but chocolate mainly contains (-)-catechin, in addition to its major flavanol, (-)-epicatechin. Previous studies have shown poor...[Read More](#)

[read abstract](#)

[Solid-phase extraction of methadone enantiomers and benzodiazepines in biological fluids by two polymeric cartridges for liquid chromatographic analysis.](#)
Hua He et. al
Journal of chromatography. B, Analytical technologies in the biomedical and life sciences, 814(2), undefined (2005-1-11)
The aim of this work was to present the advantages of two polymeric cartridges (Oasis HLB from Waters and Absolut Nexus from Varian) for the solid-phase extraction of methadone enantiomers and its major metabolite, 2-ethylidene-1,5-dimethyl-3,3-diphe...[Read More](#)


Produits associés



**Brighter
LC-MS**

LC-MS Grade Solvents and Reagents
for accurate, brilliant results

[LEARN MORE](#)



**Milli-Q® water
purification
solutions**

Ultrapure & Pure Systems
for Every Lab

[Learn More](#)

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)