




24100AST Supelco

Astec® CYCLOBOND® I 2000 HP-RSP Chiral HPLC Guard

Cartridge, 5 µm particle size, L × I.D. 2 cm × 4 mm

eCl@ss 32110501



FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
24100AST	Expédition estimée le 18.09.18	204.00	0

Commandes Bulk?

AJOUTER AU PANIER

Propriétés

Related Categories	Analytical/Chromatography, Astec CYCLOBOND Chiral HPLC Columns (Native and Derivatized Cyclodextrins), CYCLOBOND I 2000 HP-RSP, Chiral HPLC Columns, Chiral HPLC Guard Columns and Accessories, HPLC Less...
description	HPLC guard column
packaging	pkg of 1 ea
application(s)	HPLC: suitable
L × I.D.	2 cm × 4 mm
matrix	silica particle platform
matrix active group	cyclodextrin, beta- phase
particle size	5 µm
pore size	100 Å
suitability	suitable for L45 per USP
mode of chromatography	chiral

[montrer moins de résultats](#)

Description

Recommended products

Discover LiChropur reagents ideal for **HPLC** or **LC-MS** analysis

Legal Information

Astec is a registered trademark of Sigma-Aldrich Co. LLC

CYCLOBOND is a registered trademark of Sigma-Aldrich Co. LLC

Informations Sécurité

Aucune information relative à la sécurité de ce produit est actuellement disponible.

Documents

Certificat d'Analyse

Entrez le Lot N°

Devis/ Commande de produits Bulk

FDS

Customers Also Viewed

 21150AST Supelguard™ Guard Cartridge Holder Stand-Alone (Swivel-type), for use with Supelguard cartridges (2 cm L. x 2 to 4.6 mm I.D.)	24024AST Astec® CYCLOBOND® I 2000 HP-RSP Chiral HPLC Column 5 µm particle size, L × I.D. 25 cm × 4.6 mm	24098-U SPB®-1 Capillary GC Column L × I.D. 15 m × 0.32 mm, d _f 1.00 µm	24085-U SUPELCOWAX® 10 Capillary GC Column L × I.D. 60 m × 0.32 mm, d _f 0.50 µm	24080-U SUPELCOWAX® 10 Capillary GC Column L × I.D. 30 m × 0.32 mm, d _f 0.25 µm
--	---	--	--	--

Recently Viewed

24098-U SPB®-1 Capillary GC Column L × I.D. 15 m × 0.32 mm, d _f 1.00 µm	24085-U SUPELCOWAX® 10 Capillary GC Column L × I.D. 60 m × 0.32 mm, d _f 0.50 µm	24080-U SUPELCOWAX® 10 Capillary GC Column L × I.D. 30 m × 0.32 mm, d _f 0.25 µm	24045-U SPB®-1 Capillary GC Column L × I.D. 30 m × 0.32 mm, d _f 1.00 µm	24024AST Astec® CYCLOBOND® I 2000 HP-RSP Chiral HPLC Column 5 µm particle size, L × I.D. 25 cm × 4.6 mm
--	--	--	--	---

Protocoles et articles

Related Content

[Supelco HPLC Calculator for Isocratic or Gradient Method Transfer](#)
Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...
Keywords: High performance liquid chromatography

Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Produits associés

Brighter
LC-MS
LC-MS Grade Solvents and Reagents
for accurate, brilliant results



[LEARN MORE](#)

Milli-o® water purification solutions
Ultrapure & Pure Systems for Every Lab



[Learn More](#)