



54258-U Supelco

# Large Volume SPE Reservoir

polypropylene body, for use with 6 mL polypropylene SPE tubes, volume 25 mL, pk of 30



FDS

Similar Products

| Conditionnement - SKUDisponibilité |                                                                   | Prix (EUR) | Quantité                       |
|------------------------------------|-------------------------------------------------------------------|------------|--------------------------------|
| 54258-U                            | Seulement 5 en stock (d'autres en cours d'arrivage) - A PARTIR DE | 74.00      | <input type="text" value="0"/> |

Commandes Bulk?

AJOUTER AU PANIER

## Propriétés

Related Categories

Analytical/Chromatography, Large Volume SPE Reservoir, Large Volume SPE Reservoirs, SPE Frit Insertion Tool, Tube Adapters, Large Vol Reservoirs and Single Tube Processors, SPE Tube Components, Plus...

|               |                                           |
|---------------|-------------------------------------------|
| material      | polypropylene body                        |
| packaging     | pk of 30                                  |
| volume        | 25 mL                                     |
| compatibility | for use with 6 mL polypropylene SPE tubes |

## Description

**General description**  
Large volume SPE reservoirs are designed to increase the head space volume of standard polypropylene SPE tubes. Because these reservoirs are designed to connect directly to the mouth of the SPE tube, they are ideal for gravity applications where increased headspace volume is required.

## Informations Sécurité

Aucune information relative à la sécurité de ce produit est actuellement disponible.

## Documents

Certificat d'Analyse

Entrez le Lot N°

Devis/ Commande de produits Bulk FDS

### Customers Also Viewed



**54259-U**  
**Large Volume SPE Reservoir**  
PTFE body, for use with 6 mL polypropylene SPE tubes, volume 25 mL, pk of 3



**57267**  
**Cartridge Adapter for H300 Cartridges**  
pkg of 6 ea



**57020-U**  
**SPE Tube Adapter**  
for use with 1, 3, & 6 mL tubes, pk of 12



**57022**  
**Empty SPE Tube (no frits)**  
volume 60 mL, pk of 12



**55253-U**  
**Supelclean™ Sulfoxide SPE Tube**  
PE frit, bed wt. 3 g, volume 6 mL, pk of 30

## Recently Viewed

|                                                                                                                               |                                                                                    |                                                                                                                             |                                                                                                                                    |                                                                                                                            |
|-------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|
|                                              |  |                                            |                                                  |                                           |
| 54256-U<br><b>Ascentis® Express C18, 2.7 Micron Capillary HPLC Column</b><br>2.7 µm particle size, L × I.D.<br>15 cm × 100 µm | 54254-U<br><b>Supelpak™ -4</b><br>bed wt. 400 mg, size 1 mL                        | 54229-U<br><b>Ascentis® Express C8 Capillary, 2.7 micron HPLC Column</b><br>2.7 µm particle size, L × I.D.<br>15 cm × 75 µm | 54226-U<br><b>Discovery® Ag-Ion SPE Tube</b><br>polypropylene hardware (Rezorian Cartridge), bed wt. 750 mg, volume 1 mL, pk of 10 | 53345-U<br><b>Ascentis® Express Phenyl-Hexyl, 2.7 Micron HPLC Column</b><br>2.7 µm particle size, L × I.D.<br>10 cm × 3 mm |

## Protocoles et articles

### Articles

#### [AOAC First Action Official Method<sup>SM</sup> 2015.04](#)

Determination of Monofluoroacetate in Powdered Nutritional Products by Derivatization with 2-Nitrophenylhydrazine and LC-MS/MS  
Keywords: Derivatizations, High performance liquid chromatography, Liquid chromatography mass spectrometry, Mass spectrometry

## Related Content

#### [SPE Cartridges \(Tubes\) - Phase Selection Guide](#)

The SPE Tube Configuration Guide will help you choose the appropriate bed weight and cartridge dimensions.

Keywords: Agriculture, Calorimetry, Environmental, Fractionation, High performance liquid chromatography, Liquid chromatography mass spectrometry, Mass spectrometry, Pesticides, Pharmaceutical, Solid phase extractions

## Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

## Produits associés

**Milli-o® water purification solutions**

Ultrapure & Pure Systems for Every Lab

[Learn More](#)



© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)