




55050-U Supelco

RheBuild® Kit

for use with 8125/8126 valves



FDS

[Similar Products](#)

Conditionnement - SKUDisponibilité		Prix (EUR)	Quantité
55050-U	Seulement 1 en stock (d'autres en cours d'arrivage) - A PARTIR DE	246.00	0

[Commandes Bulk?](#)

[AJOUTER AU PANIER](#)

Propriétés

Related Categories	Analytical/Chromatography, For 8000 Series, HPLC, HPLC Accessories, RheBuild Kits, Plus...
compatibility	for use with 8125/8126 valves

Description

General description
These kits include all of the tools and parts you need to repair your Rheodyne injector. The parts included are those most likely to be damaged, worn, or lost. Kits for front-loading injectors include: rotor seal, stator face assembly, isolation seal, needle guide, needle port cleaner, 2 hex keys, operating instructions, and a mini-manual. Parts may vary, depending on valve model.

Recommended products
Discover LiChropur reagents ideal for [HPLC](#) or [LC-MS](#) analysis

Legal Information
RheBuild is a registered trademark of IDEX Health & Science, LLC

Informations Sécurité

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Entrez le Lot N°

[Devis/ Commande de produits Bulk](#) [FDS](#)

Customers Also Viewed

57955 Rheodyne® 8125 Rotor Seal for use with 8125 Injector	54612 Rheodyne® Stator Face Assy for 8125/8126 (#8125-074)	 57950-U Rheodyne® Model 8125 Low Dispersion Injector	504602 RheBuild® Kit for use with 7010 with Stator	55052 RheBuild® Kit for use with 9010/9040 Injector
---	--	--	---	--

Recently Viewed

		
---	--	---

55009-U Super Flangeless Ferrules for 1/8 in. tubing	55008-U Super Flangeless Ferrules for 1/16 in. tubing	55007-U Tube Fitting, Sealtight™ ferrule pack	55006-U Tube Fitting, Sealtight™ Extra-long nut, 10-32, with ferrule	55003-U Tube Fitting, Sealtight™ Long nut, 10-32, with ferrule
---	--	--	--	---

Protocoles et articles

Related Content

[Supelco HPLC Calculator for Isocratic or Gradient Method Transfer](#)

Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...

Keywords: High performance liquid chromatography

Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Produits associés



© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)