




56710-U Supelco

Stainless Steel 1/16 in. Capillary Tubing

L × O.D. × I.D. 30 cm × 1/16 in. × 0.005 in.

eCl@ss 32119290



FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
56710-U	Disponible pour expédition le 31.08.18 - A PARTIR DE	14.70	0

[Commandes Bulk?](#)

AJOUTER AU PANIER

Propriétés

Related Categories	30 cm length, Analytical/Chromatography, HPLC, HPLC Accessories, Stainless Steel Couplings, Unions, Tees, Tubing, Plus...
L × O.D. × I.D.	30 cm × 1/16 in. × 0.005 in.
color	red band

Description

General description
316 stainless steel, 1/16 in. O.D. 1 m length is coiled, all other dimensions are straight.

Recommended products
Discover LiChropur reagents ideal for [HPLC](#) or [LC-MS](#) analysis

Informations Sécurité

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Entrez le Lot N°

[Devis/ Commande de produits Bulk](#) [FDS](#)

Customers Also Viewed

56711 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 0.5 m × 1/16 in. × 0.005 in.	56709 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 20 cm × 1/16 in. × 0.005 in.	56722 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.010 in.	56716 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.007 in.	56740-U Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.040 in.
---	---	---	---	---

Recently Viewed

56666-U Micro-Mixer Static Mixers (2 - 25 µL) housing, Binary	56665-U Micro-Mixer Static Mixers (2 - 25 µL) housing, In-Line	56664-U Micro-Mixer Static Mixers (2 - 25 µL) 25 µL cartridge, stainless steel	56661-U Micro-Mixer Static Mixers (2 - 25 µL) 2 µL cartridge, stainless steel	56634-U LC-MS Post Column Resistor Sets Binary, Split Ratio = 3:1
---	--	--	---	---

Protocoles et articles

Related Content


[Supelco HPLC Calculator for Isocratic or Gradient Method Transfer](#)

Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...
Keywords: High performance liquid chromatography

Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Produits associés



© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)