




56720-U Supelco

Stainless Steel 1/16 in. Capillary Tubing

L × O.D. × I.D. 10 cm × 1/16 in. × 0.010 in.

eCl@ss 32119290



FDS

[Similar Products](#)

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
56720-U	Expédition estimée le 19.09.18	7.00	0

[Commandes Bulk?](#)

[AJOUTER AU PANIER](#)

Propriétés

Related Categories	10 cm length, Analytical/Chromatography, HPLC, HPLC Accessories, Stainless Steel Couplings, Unions, Tees, Tubing, Plus...
L × O.D. × I.D.	10 cm × 1/16 in. × 0.010 in.
color	blue band

Description

General description
316 stainless steel, 1/16 in. O.D. 1 m length is coiled, all other dimensions are straight.

Recommended products
Discover LiChropur reagents ideal for [HPLC](#) or [LC-MS](#) analysis

Informations Sécurité

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Entrez le Lot N°

[Devis/ Commande de produits Bulk](#) [FDS](#)

Customers Also Viewed

56721 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 20 cm × 1/16 in. × 0.010 in.	56714 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 10 cm × 1/16 in. × 0.007 in.	56722 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.010 in.	56732 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 10 cm × 1/16 in. × 0.030 in.	56719 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 5 cm × 1/16 in. × 0.010 in.
---	---	---	---	--

Recently Viewed

--	--	--	--	--

567209-U Discovery® BIO Wide Pore C18 HPLC Column 10 µm particle size, L × I.D. 25 cm × 10 mm	567208-U Discovery® BIO Wide Pore C18 HPLC Column 10 µm particle size, L × I.D. 15 cm × 10 mm	567207-U Discovery® BIO Wide Pore C18 HPLC Column 10 µm particle size, L × I.D. 5 cm × 10 mm	567205-U Discovery® BIO Wide Pore C18 HPLC Column 3 µm particle size, L × I.D. 15 cm × 4.6 mm	567204-U Discovery® BIO Wide Pore C18 HPLC Column 3 µm particle size, L × I.D. 10 cm × 4.6 mm
--	--	---	--	--

Protocoles et articles

Related Content

[Supelco HPLC Calculator for Isocratic or Gradient Method Transfer](#)

Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...

Keywords: High performance liquid chromatography

Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Produits associés

