

SIGMA-ALDRICH®


56734-U Supelco

Stainless Steel 1/16 in. Capillary Tubing

L × O.D. × I.D. 30 cm × 1/16 in. × 0.030 in.

Synonym: 0.030IN ID SS TUBING, 1/16IN OD, 30CM

eCl@ss 32119290



FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
56734-U	Expédition estimée le 13.09.18	22.50	0

Commandes Bulk?

AJOUTER AU PANIER

Propriétés

Related Categories	30 cm length, Analytical/Chromatography, HPLC, HPLC Accessories, Stainless Steel Couplings, Unions, Tees, Tubing, Plus...
material	white band
L × O.D. × I.D.	30 cm × 1/16 in. × 0.030 in.

Description

General description
316 stainless steel, 1/16 in. O.D. 1 m length is coiled, all other dimensions are straight.

Recommended products
Discover LiChropur reagents ideal for HPLC or LC-MS analysis

Informations Sécurité

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Entrez le Lot N°

Devis/ Commande de produits Bulk FDS

Customers Also Viewed

56740-U Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.040 in.	56732 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 10 cm × 1/16 in. × 0.030 in.	56710-U Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.005 in.	56716 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.007 in.	56728 Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 30 cm × 1/16 in. × 0.020 in.
---	---	---	---	---

Recently Viewed

56733-U Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 20 cm × 1/16 in. × 0.030 in.	56730-U Stainless Steel 1/16 in. Capillary Tubing L × O.D. × I.D. 1.0 m × 1/16 in. × 0.020 in.	567284-U Discovery® BIO Wide Pore C8 Supelguard™ Cartridge 10 µm particle size, L × I.D. 1 cm × 10 mm	567282-U Discovery® BIO Wide Pore C18 Supelguard™ Cartridge 10 µm particle size, L × I.D. 1 cm × 10 mm	567281-U Discovery® BIO Wide Pore C5 Supelguard™ Cartridge 3 µm particle size, L × I.D. 2 cm × 4 mm, pkg of 1 kit
--	--	---	--	---

Protocoles et articles

Related Content

[Supelco HPLC Calculator for Isocratic or Gradient Method Transfer](#)

Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...
Keywords: High performance liquid chromatography

Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Produits associés



© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)