


France Accueil FFA57289-U - SPME fiber assembly Polydimethylsiloxane (PDMS)



FFA57289-U Supelco

SPME fiber assembly Polydimethylsiloxane (PDMS)

df 30 µm (nonbonded), needle size 23 ga, for use with multi fiber exchanger



FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
FFA57289-U	Expédition estimée le 20.09.18	480.00	<input type="text" value="0"/>

Commandes Bulk?

AJOUTER AU PANIER

Propriétés

Related Categories	Analytical/Chromatography, SPME Fast Fit Fiber Assemblies (FFA), SPME Fast Fit Fiber Assemblies (FFA) & Accessories, SPME Fibers and Holders, Sample Preparation & Purification, Solid Phase Microextraction (SPME) Less...
material	fused silica fiber
	yellow hub plain
needle size	23 ga
packaging	pkg of 3 × ea
df	30 µm (nonbonded)
matrix active group	Polydimethylsiloxane (PDMS) coating
compatibility	for analyte group volatiles (MW 60-275)
	for use with multi fiber exchanger

[montrer moins de résultats](#)

Informations Sécurité

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Entrez le Lot N°

Devis/ Commande de produits Bulk FDS

Customers Also Viewed

<p>FFA57301</p> <p>SPME fiber assembly Polydimethylsiloxane (PDMS)</p> <p>df 100 µm (nonbonded phase, for use with multi fiber exchanger, needle size 24 ga</p>	<p>FFA57293-U</p> <p>SPME Fiber Assembly Poly dimethylsiloxane/Divinylben zene (PDMS/DVB)</p> <p>df 65 µm, needle size 23 ga, Stableflex, for use with multi fiber exchanger</p>	<p>FFA57305</p> <p>SPME fiber assembly polyacrylate (PA)</p> <p>df 85 µm (partially crosslinked phase, for use with multi fiber exchanger, needle size 24 ga</p>	<p>FFA57295-U</p> <p>SPME fiber assembly Carbo xen/Polydimethylsiloxane (CAR/PDMS)</p> <p>df 85 µm, needle size 23 ga, Stableflex, for use with multi fiber exchanger</p>	<p>FFA57362-U</p> <p>SPME fiber assortment kit 5</p> <p>needle size 23 ga, for use with multi fiber exchanger</p>
---	--	--	---	---

Recently Viewed

--	--	--	--	--

FFA57291-U SPME fiber assembly Polydimethylsiloxane (PDMS) df 7 µm (bonded, needle size 23 ga, for use with multi fiber exchanger	FFA57287-U SPME fiber assortment kit 4 for use with multi fiber exchanger, needle size 23 ga	FFA57286-U SPME fiber assortment kit 2 needle size 23 ga, for use with multi fiber exchanger	FFA57285-U SPME fiber assortment kit 1 for use with multi fiber exchanger, needle size 23 ga	FFA57284-U SPME StableFlex™ fiber assortment kit needle size 23 ga, for use with multi fiber exchanger
---	--	--	--	--

Protocoles et articles

Related Content

Fundamentos basicos de microextraccion en fase solida (SPME) y sus aplicaciones [Webinar]
En este webinar presentamos los fundamentos básicos de microextracción en fase sólida (SPME) y sus aplicaciones. SPME es una técnica para análisis de compuestos volátiles y semi-volátiles, fácil de u...
Keywords: Solid phase microextractions

SUPELCO® SPME Fast-Fit Fiber Assemblies [VIDEO]
Ideal for SPME method development, SPME Fast-Fit Fiber Assemblies (FFA) from Supelco are easy to install, feature a non-threaded attachment, and are equipped with a traceable color-coded barcode lab...
Keywords: Solid phase microextractions

Documentation référencée

Did you use this product in your Paper? If so click here.

[read abstract](#)
Analysis of free fatty acids in whey products by solid-phase microextraction.
R M Tomaino et. al
Journal of agricultural and food chemistry, 49(8), undefined (2001-8-22)
To evaluate the impact of Cheddar cheese starter cultures on the level of free fatty acids in liquid whey, a solid-phase microextraction (SPME) technique was utilized. The determination of response factors relative to an internal standard and the ver...[Read More](#)

[read abstract](#)
Solid-phase microextraction to predict bioavailability and accumulation of organic micropollutants in terrestrial organisms after exposure to a field-contaminated soil.
Leon van der Wal et. al
Environmental science & technology, 38(18), undefined (2004-10-19)
The risk posed by soil contaminants strongly depends on their bioavailability. In this study, a partition-based sampling method was applied as a tool to estimate bioavailability in soil. The accumulation of organic micropollutants was measured in two...[Read More](#)

[read abstract](#)
Measurement of trihalomethanes in potable and recreational waters using solid phase micro extraction with gas chromatography-mass spectrometry.
M A Stack et. al
Chemosphere, 41(11), undefined (2000-11-1)
Solid phase micro extraction (SPME) was applied to the determination of selected trihalomethanes (THMs), chloroform, bromodichloromethane, dibromochloromethane, bromoform, in potable and recreational waters. The selected samples were environmentally ...[Read More](#)

[read abstract](#)
Evidence for selectivity of absorption of volatile organic compounds by a polydimethylsiloxane solid-phase microextraction fibre.
S Niedziella et. al
Journal of chromatography. A, 885(1-2), undefined (2000-8-15)
Solid-phase microextraction using a 30 microns polydimethylsiloxane fibre has been used to sample the volatile organic compounds from standard mixtures and from mixtures produced by the decomposition of organic compounds. This method of sampling has ...[Read More](#)

[read abstract](#)
2-Methyl-3-furanthiol and methional are possible off-flavors in stored orange juice: aroma-similarity, NIF/SNIF GC-O, and GC analyses.
Y Bezman et. al
Journal of agricultural and food chemistry, 49(11), undefined (2001-11-21)
The occurrence of methional in fresh orange juice, and possible occurrence of beta-damascenone in heated orange juice, has been previously suggested. Here we report on the occurrence of 2-methyl-3-furanthiol in the headspace, collected by solid-phase...[Read More](#)

[read abstract](#)
Influence of vine training and sunlight exposure on the 3-alkyl-2-methoxypyrazines content in musts and wines from the Vitis vinifera variety cabernet sauvignon.
Cristina Sala et. al
Journal of agricultural and food chemistry, 52(11), undefined (2004-5-27)
The influence of vine training and sunlight exposure on the 3-alkyl-2-methoxypyrazines contents in musts and wines was studied by means of two previously reported methods based on headspace solid-phase micro-extraction. Experimental samples were moni...[Read More](#)

[read abstract](#)
High-throughput analysis of bergamot essential oil by fast solid-phase microextraction-capillary gas chromatography-flame ionization detection.
Peter Quinto Tranchida et. al
Journal of chromatography. A, 1103(1), undefined (2005-12-21)
The advantages of using a narrow-bore column in headspace solid-phase microextraction-gas chromatographic (HS-SPME-GC) analysis are investigated. An automated rapid HS-SPME-GC method for the determination of volatile compounds in a complex sample (be...[Read More](#)

[read abstract](#)
[Preparation of a novel activated carbon coating fiber for solid phase micro-extraction and its application for halocarbon compound analysis in water].
Shutao Wang et. al
Se pu = Chinese journal of chromatography, 22(5), undefined (2005-2-15)
A novel activated carbon coating fiber used for solid phase micro-extraction (SPME) was prepared using activated carbon powder and silica resin adhesive. The extraction properties of the novel activated carbon coating fiber were investigated. The res...[Read More](#)

[read abstract](#)
The scent of Mycobacterium tuberculosis.
Mona Syhre and Stephen T Chambers
Tuberculosis (Edinburgh, Scotland), 88(4), undefined (2008-2-26)
Worldwide, tuberculosis (TB) kills nearly 2 million people annually, yet rapid diagnosis still relies on a 100-year-old method of sputum staining for acid-fast bacilli. The advent of solid phase micro-extraction and gas chromatography/mass spectromet...[Read More](#)

[read abstract](#)
Experimental design for the optimization of the extraction conditions of polycyclic aromatic hydrocarbons in milk with a novel diethoxydiphenylsilane solid-phase microextraction fiber.
F Bianchi et. al
Journal of chromatography. A, 1196-1197, undefined (2008-5-6)

An innovative solid-phase microextraction coating based on the use of diethoxydiphenylsilane synthesized by sol-gel technology was used for the determination of polycyclic aromatic hydrocarbons at trace levels in milk. The effects of time and tempera...[Read More](#)

[read abstract](#)

Coating of peanuts with edible whey protein film containing alpha-tocopherol and ascorbyl palmitate.

J H Han et. al

Journal of food science, 73(8), undefined (2008-11-21)

Physical properties of whey protein isolate (WPI) coating solution incorporating ascorbic palmitate (AP) and alpha-tocopherol (tocopherol) were characterized, and the antioxidant activity of dried WPI coatings against lipid oxidation in roasted peanu...[Read More](#)

[read abstract](#)

Simultaneous determination of UV filters and polycyclic musks in aqueous samples by solid-phase microextraction and gas chromatography-mass spectrometry.

Hongtao Liu et. al

Journal of chromatography. A, 1217(43), undefined (2010-6-29)

A simple, precise and accurate method for the simultaneous determination of four UV filters and five polycyclic musks (PCMs) in aqueous samples was developed by solid-phase microextraction coupled with gas chromatography-mass spectrometry (SPME-GC-MS...[Read More](#)

[read abstract](#)

Solid-phase micro-extraction procedure for the determination of 1,3-dichloro-2-propanol in water by on-fibre derivatisation with bis(trimethylsilyl)trifluoroacetamide.

Antonia Maria Carro et. al

Analytical and bioanalytical chemistry, 394(3), undefined (2009-4-11)

The headspace solid-phase micro-extraction technique with on-fibre derivatisation followed by gas chromatography-tandem mass spectrometry has been evaluated for the analysis of 1,3-dichloro-2-propanol in water. An asymmetric factorial design has been...[Read More](#)

[read abstract](#)

Solid-phase micro-extraction (SPME) in the early detection of potentially active volatile compounds from organic wastes used for the management of soil-borne pathogens.

Esther De La Fuente et. al

Journal of environmental science and health. Part A, Toxic/hazardous substances & environmental engineering, 44(10), undefined (2009-10-16)

The complex molecular assemblages were analysed in the soil gas phase after applying pine forest wastes (PFW) or sugarbeet vinasses (SBV) for soil-borne crop pests' management. For this purpose, solid-phase micro-extraction (SPME) and gas chromatogra...[Read More](#)

[read abstract](#)

Optimisation of the headspace-solid phase microextraction for organomercury and organotin compound determination in sediment and biota.

Alejandra Delgado et. al

Journal of separation science, 31(4), undefined (2008-2-2)

Headspace solid-phase microextraction was optimised for the simultaneous preconcentration of methylmercury (MeHg⁺), monobutyltin, dibutyltin, tributyltin, monophenyltin (MPhT), diphenyltin (DPhT), and triphenyltin (TPhT) from sediments and biota. Ext...[Read More](#)

[Show more references \(2013 remain\)](#)

Produits associés



© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)