



Accueil > Life Science Research > Protein Detection and Quantification > Western Blotting > Western Blotting Kits and Reagents > Bløk - CH Noise Cancelling Reagents for Chemiluminescence Detection, 500 mL

## WBAVDCH01 | Bløk - CH Noise Cancelling Reagents for Chemiluminescence Detection, 500 mL



**Chemiluminescent detection kit used for Western blotting, dot/slot blotting, mass spectrometry**

### WBAVDCH01

1 Sufficient reagent for 19 mini blots (8x7 cm)

Le prix n'a pas pu être récupéré

La quantité minimale doit être un multiple de

Les prix peuvent être modifiés sans préavis

[Aperçu](#) [Documentation complémentaire](#) [Produits & Applications associés](#)

### Aperçu

Description  
Product Information  
Applications  
Dimensions  
Product Usage Statements  
Storage and Shipping Information  
Packaging Information  
Tableau de caractéristiques principal

#### Tableau de caractéristiques principal

| Key Applications    |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| WB                  |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Catalogue Number    | WBAVDCH01                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Trade Name          | <ul style="list-style-type: none"> <li>bløk</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                        |
| Description         | Bløk - CH Noise Cancelling Reagents for Chemiluminescence Detection, 500 mL                                                                                                                                                                                                                                                                                                                                                                   |
| Overview            | Bløk Noise Cancelling Reagents are designed especially for Western blotting applications. This family of protein-free, ready-to-use buffers has been optimized to reduce background levels when using chemiluminescence, fluorescent, or phosphotyrosine detection. Bløk buffers are ideal for use with Immobilon®-P Western blotting membranes and the SNAP i.d.® system, and are also compatible with standard Western blotting techniques. |
| Product Information |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Detection method    | Chemiluminescence                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Format              | HRP                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Applications        |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Application         | Chemiluminescent detection kit used for Western blotting, dot/slot blotting, mass spectrometry                                                                                                                                                                                                                                                                                                                                                |
| Key Applications    | <ul style="list-style-type: none"> <li>Western Blotting</li> </ul>                                                                                                                                                                                                                                                                                                                                                                            |
| Application Notes   | For use on immunoblots.                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Dimensions          |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Membrane Coverage   | 1000 cm <sup>2</sup>                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Volume              | 500 mL                                                                                                                                                                                                                                                                                                                                                                                                                                        |

| Product Usage Statements |                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Quality Assurance        | For optimum results, use Immobilon-P membrane for chemiluminescence (luminol) or visual detection (e.g. TMB), and Immobilon-FL membrane for fluorescent detection applications.                                                                                                                                                                                                                                                                      |
| Usage Statement          | <ul style="list-style-type: none"> <li>Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.</li> </ul> |

| Storage and Shipping Information |                               |
|----------------------------------|-------------------------------|
| Storage Conditions               | Stable for 2 years at 20–25°C |
| Storage Temperature              | 20 to 25 °C                   |

| Packaging Information |                                               |
|-----------------------|-----------------------------------------------|
| Material Size         | 1                                             |
| Material Package      | Sufficient reagent for 19 mini blots (8x7 cm) |

[Nous contacter](#)

[Groupe Merck](#) | 
 [Mentions légales](#) | 
 [Conditions d'utilisation](#) | 
 [Respect de la vie privée](#) | 
 [Conditions de vente](#)

© 2018 Merck KGaA, Darmstadt, Allemagne et/ou ses filiales. Tous droits réservés.

© Merck KGaA, Darmstadt, Allemagne, 2014. Toutes les références à Merck désignent Merck KGaA, Darmstadt, Allemagne.