

Z227307

Supelco

PEEK Tubing

L × O.D. × I.D. 10 ft × 1/16 in. × 0.005 in., red, parameter 8000 max. pressure

eCl@ss 32119290

FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (GBP)	Quantité
Z227307	Disponible pour expédition le 07.06.19 - A PARTIR DE	45.50	0

Bulk orders?

AJOUTER AU PANIER

Product Recommendations



Z226661

PEEK Tubing

L × O.D. × I.D. 10 ft × 1/16 in. × 0.010 in., blue, parameter 6000 psi max. pressure



Z226688

PEEK Tubing

L × O.D. × I.D. 10 ft × 1/16 in. × 0.007 in., yellow, parameter 8000 max. pressure



Z227293

PEEK Tubing

L × O.D. × I.D. 10 ft × 1/16 in. × 0.020 in., orange, parameter 6000 psi max. pressure

Properties

Related Categories	Analytical/Chromatography , HPLC , HPLC Accessories , PEEK Tubing , PEEK Tubing , PEEKsil Tubing , PTFE Tubing Plus...
parameter	8000 max. pressure
L × O.D. × I.D.	10 ft × 1/16 in. × 0.005 in.
color	red

Description

General description

Polyetheretherketone (PEEK) tubing has become a popular replacement for stainless steel tubing in various places in the HPLC system. It is especially useful when contact between the sample and metal components must be avoided, such as when working with metal complexing agents or certain biochemical compounds. PEEK tubing also offers other advantages: it is flexible, easy to cut, has excellent mechanical stability and chemical compatibility, and can be used to 100 °C. Tetrahydrofuran (THF), dimethyl sulfoxide (DMSO), methylene chloride, and concentrated nitric and sulfuric acids should not be used with PEEK tubing. Our 1/16 in. PEEK tubing is color coded for easy identification of the internal diameter. 10 ft./3 m length.

Recommended products

Discover LiChropur reagents ideal for [HPLC](#) or [LC-MS](#) analysis

Safety Information

Safety Information for this product is unavailable at this time.

Documents

Certificat d'Analyse

Enter Lot No.

[Devis/ Commande de produits Bulk](#)[FDS](#)

Customers Also Viewed

 Z226955 PEEK Tubing L × O.D. × I.D. 10 ft × 1/16 in. × 0.030 in., green, parameter 5000 psi max. pressure	 54994 PEEK Tubing L × O.D. × I.D. 10 ft × 1/16 in. × 0.055 in., natural, parameter 1000 psi max. pressure	 54995 PEEK Tubing L × O.D. × I.D. 10 ft × 1/8 in. × 0.062 in., natural	 Z227250 Supelco® PEEK Fitting One-piece fingertight, pkg of 5 ea	 55067U Upchurch PEEK Fingertight Fitting fitting 10-32 (one-piece), pkg of 10 ea
--	--	---	--	---

Recently Viewed

 SU860197 2in1 kit universal vial, short thread autosampler vials pkg of 100 ea, volume 0.3 mL, blue polypropylene cap (with 6 mm center hole), clear	 29104U Vials, open-top screw cap, 2 mL standard opening with closures, convenience pack, unassembled clear glass vial, PTFE/silicone septum, pkg of × 100 ea	 27001U Vials, screw top with solid green Thermoset cap with PTFE liner, preassembled, pkg of 100 volume 4 mL, amber glass vial, O.D. × H 15 mm ×	 854996 Crimp seals with PTFE/silicone septa silver aluminum seal, open center (10 mm opening), diam. × thickness 20 mm × 3.25 mm, tan PTFE/silicone, pkg of	 27273 Assembled screw cap with hole with PTFE/silicone septum black polypropylene hole cap, thread 10-425, PTFE/silicone, septum thickness 1.5 mm, for use
---	---	---	---	---

Milli-o® water purification solutions

Ultrapure & Pure Systems for Every Lab

Learn More



SAME PRODUCT
NEW PACKAGE

Our new packaging features 2D barcodes to provide you with more information about our products.

LEARN MORE



Protocols & Articles

Related Content

Supelco HPLC Calculator for Isocratic or Gradient Method Transfer

Calculates conditions for transfer of an isocratic or gradient method from one HPLC column to another. Allows method scaling from microbore through preparative column range. Calculates time and solve...
Keywords: High performance liquid chromatography

Peer-Reviewed Papers

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Related Products

