




52604-U Supelco

Discovery® DSC-18 SPE Tube

bed wt. 500 mg, volume 6 mL, pk of 30



FDS

Similar Products

Conditionnement - SKU	Disponibilité	Prix (GBP)	Quantité
52604-U	Seulement 5 en stock (d'autres en cours d'arrivage) - A PARTIR DE	93.50	0

Bulk orders?

AJOUTER AU PANIER

Product Recommendations

52599U

Discovery® DSC-18 SPE Tube

bed wt. 500 mg, volume 6 mL, pk of 1000

52615U

Discovery® DSC-18Lt SPE Tube

bed wt. 500 mg, volume 6 mL, pk of 30

52603U

Discovery® DSC-18 SPE Tube

bed wt. 500 mg, volume 3 mL, pk of 54

Properties

Related Categories	β-Glucuronidase Drug Analysis Bundle, Analytical and Industrial Enzymes, Analytical/Chromatography, Biochemicals and Reagents, Discovery DSC-18 SPE Products, Plus...
material	PE frit (20 μm)
	polypropylene
feature	endcapped
composition	bed wt., 500 mg
packaging	pk of 30
extent of labeling	18% C loading
volume	6 mL
	montrer plus (14)

Description

General description
Retention Mechanism: Reversed-phase
Sample Matrix Compatibility: Aqueous solutions (biological fluids, water)

- Polymerically bonded, octadecyl (18 %C), endcapped
- Higher 18 %C loading for increased binding capacities and higher recoveries
- The least selective phase: retains most organic analytes from aqueous matrices
- Can also be used for desalting aqueous matrices
- Beneficial for extracting structurally diverse analytes from the same sample

Legal Information
Discovery is a registered trademark of Sigma-Aldrich Co. LLC

Safety Information

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Enter Lot No.

Devis/ Commande de produits Bulk FDS

Customers Also Viewed



52606U
Discovery® DSC-18 SPE Tube
bed wt. 1 g, volume 6 mL, pk of 30



57054
Supelclean™ LC-18 SPE Tube
bed wt. 500 mg, volume 6 mL, pkg of 30 ea



57064
Supelclean™ ENVI™ -18 SPE Tube
bed wt. 500 mg, volume 6 mL, pkg of 30 ea



52688U
Discovery® DSC-SCX SPE Tube
bed wt. 500 mg, volume 6 mL, pk of 30



52609U
Discovery® DSC-18 SPE Tube
bed wt. 10 g, volume 60 mL, pk of 16

Recently Viewed



52654U
Discovery® DSC-Si SPE Tube
bed wt. 500 mg, volume 3 mL, pk of 54



27003
Vials, screw top with solid green Thermoset cap with PTFE liner, preassembled, pkg of 100
volume 15 mL, amber glass vial, O.D. × H



57299U
SPME fiber assembly Di vinylbenzene/Carboxen/ Polydimethylsiloxane (D VB/CAR/PDMS)
needle size 23 ga, StableFlex, for use with manual



SU860083
Vials, screw top, amber glass, 1.5 mL, 8-425 thread
O.D. × H × I.D. 11.6 mm × 32 mm × 4.6 mm, pkg of 100 ea



33326
Vials for Waters® 96-position carousel, 8 x 40 mm
volume 1 mL, clear glass vial, vial O.D. × H 8 mm × 40 mm, pkg of 1000 ea

Milli-Q® water purification solutions

Ultrapure & Pure Systems for Every Lab

Learn More



Protocols & Articles

Articles

[Solid Phase Extraction: Reversed-Phase Methodology](#)

Reversed-phase SPE is considered the least selective retention mechanism when compared to normalphase or ion-exchange SPE. In other words, it may be difficult for a reversed-phase method or bonded ch...

Keywords: Environmental, Evaporation, Gas chromatography, High performance liquid chromatography, Homogenization, Metabolites, Phase transitions, Solid phase extractions, Solvents

Related Content

[SPE Cartridges \(Tubes\) - Phase Selection Guide](#)

The SPE Tube Configuraton Guide will help you choose the appropriate bed weight and cartridge dimensions.

Keywords: Agriculture, Calorimetry, Environmental, Fractionation, High performance liquid chromatography, Liquid chromatography mass spectrometry, Mass spectrometry, Pesticides, Pharmaceutical, Solid phase extractions

[Solid Phase Extraction: Fundamentals and Recent Developments \[Webinar\]](#)

Solid Phase Extraction (SPE) has become a popular technique for sample preparation due to its convenience, effectiveness, relatively low solvent usage and versatility. The focus of this webinar is on...

Keywords: Sample preparations, Solid phase extractions

Peer-Reviewed Papers

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Related Products

© 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Reproduction of any materials from the site is strictly forbidden without permission. Sigma-Aldrich Products are sold exclusively through Sigma-Aldrich, Inc. [Site Use Terms](#) | [Privacy](#)