



Z203688 Aldrich

PTFE syringe tubing

gauge 18, L 36 in.

Synonym: ® fun



FDS

[Similar Products](#)

Conditionnement - SKU	Disponibilité	Prix (EUR)	Quantité
Z203688-1EA	Seulement 1 en stock (d'autres en cours d'arrivage) - A PARTIR DE	35.50	0

[Bulk orders?](#)

[AJOUTER AU PANIER](#)

Product Recommendations



Z203653
PTFE syringe tubing
gauge 20, L 36 in.



Z117366
PTFE syringe tubing
gauge 12, L 12 in.



Z117374
PTFE syringe tubing
gauge 20, L 24 in.

Properties

Related Categories	Aldrich , Aldrich Syringes & Needles , Analytical/Chromatography , Labware , Manufacturer Browser , Plus...
L	36 in.
gauge	18

Description

General description
With KEL-F Luer hub.

A septum to be pierced with these PTFE needles must be punctured with a hollow SS needle. The PTFE needle is then threaded through. After threading, the SS needle is withdrawn, allowing the septum to seal around the PTFE needle.

For 20 ga PTFE use Z126756 (14 ga SS); 18 ga PTFE use Z126764 (13 ga SS); 16 ga PTFE use Z126772 (12 ga SS).

Legal Information
QuadraSii is a registered trademark of Johnson Matthey Finland Oy

Safety Information

Safety Information for this product is unavailable at this time.

Documents

Certificat d'Analyse

[Devis/ Commande de produits Bulk](#) [FDS](#)
[Needle Gauge Guide \(29 KB \)](#)

Customers Also Viewed

 Z117315 PTFE syringe tubing gauge 20, L 12 in.	 Z203726 PTFE syringe tubing gauge 12, L 36 in.	 Z203718 PTFE syringe tubing gauge 14, L 36 in.	 Z117382 PTFE syringe tubing gauge 18, L 24 in.	 Z203696 PTFE syringe tubing gauge 16, L 36 in.
--	--	--	---	--

Recently Viewed

 Z117382 PTFE syringe tubing gauge 18, L 24 in.	 Z117323 PTFE syringe tubing gauge 18, L 12 in.	 Z203696 PTFE syringe tubing gauge 16, L 36 in.	 Z117390 PTFE syringe tubing gauge 16, L 24 in.	 Z117331 PTFE syringe tubing gauge 16, L 12 in.
---	---	---	--	---

Protocols & Articles

Related Content

[Syringe Needle Gauge Chart](#)
From our library of Related Content, Sigma-Aldrich presents Syringe Needle Gauge Chart

Peer-Reviewed Papers

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Related Products