




AP1003700 | Tampon absorbant, 37 mm



37 mm, cellulose support pads are used to reinforce filters in Millipore's monitors for contamination analysis. When saturated with growth medium, they can also be used for culturing microorganisms [More>>](#)

AP1003700

100

Le prix n'a pas pu être récupéré

La quantité minimale doit être un multiple de

Les prix peuvent être modifiés sans préavis

[Aperçu](#) [Documentation complémentaire](#) [Produits & Applications associés](#)

Aperçu

Description
Informations produit
Applications
Informations biologiques
Informations physico-chimiques
Dimensions
Informations sur les matériaux
Informations sur l'emballage

Description	
Référence	AP1003700
Description	Tampon absorbant, 37 mm
Informations générales	Cellulose support pads are used to reinforce filters in Millipore's monitors for contamination analysis. When saturated with growth medium, they can also be used for culturing microorganisms on membranes. To increase flow rate and throughput, our woven mesh spacers can be placed between filters in a serial filtration to prevent the downstream screen filter from "blinding" the upstream filter pores. Features & Benefits • Avoid membrane damage during high pressure or fast flow conditions • Spacers enable you to combine multiple stages of filtration into a single process step, saving time Applications: Air Monitoring, Environmental Monitoring, Aerosol Contamination Monitoring

Informations produit

Code du filtre	AP10
Couleur du filtre	Blanc
Type de filtre	Filtre en profondeur

Applications

Application	37 mm, cellulose support pads are used to reinforce filters in Millipore's monitors for contamination analysis. When saturated with growth medium, they can also be used for culturing microorganisms on membranes. Comes in a pack of 100.
-------------	---

Informations biologiques

Mouillabilité	Hydrophile
---------------	------------

Informations physico-chimiques

Informations physico-chimiques	
Dimension de pores	Not Applicable
Débit d'eau	• 5.8 mL/min x cm²

Dimensions	
Surface du filtre	Unie
Diamètre du filtre (ø)	37 mm

Informations sur les matériaux	
Chimie	• Fibre de cellulose

Informations sur l'emballage	
Quantité	100

[Nous contacter](#)

[Groupe Merck](#) | [Mentions légales](#) | [Conditions d'utilisation](#) | [Respect de la vie privée](#) | [Conditions de vente](#)

© 2018 Merck KGaA, Darmstadt, Allemagne et/ou ses filiales. Tous droits réservés.

© Merck KGaA, Darmstadt, Allemagne, 2014. Toutes les références à Merck désignent Merck KGaA, Darmstadt, Allemagne.