



Biosigma s.r.l.

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |
Fax ++39 0426 302228 | SMS ++39 348 4077376 |
E-mail info@biosigmaeu.com | <http://www.biosigma.com> |

M521D

Last Updating
18th January 2013

TECHNICAL DATA SHEET

REF: BSA022



CONE UTILIZABLE ON EPPENDORF, WITH CAP, GRADUATED, IN PP NEUTRAL

Box 16 X 500 Pcs.



CONETTI		Vol. ml	Dim. mm
BSA022	Cone Utilizable on Eppendorf 1,5 ml PP NEUTRAL	1,5	Ø 10 x 40

IN VITRO DIAGNOSTIC MEDICAL DEVICE CE FOR PROFESSIONAL USE ONLY

- Use only for laboratory analysis
- Use gloves: contamination risk if contents leaks
- Keep out of flame or heat sources which might damage the product
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place, temperature range: min. - 10°C max 50°C
- Can be centrifuged up to 16.400 rpm or 25.000 x g
- Temperature of use: -80°C +121 °C,
- Autoclavable open at +121°C for 20 minutes,
- Disposal: use appropriate personal protective equipment and act according to applicable regulations