

## **HARD TISSUE IMPACT GRINDING MIX W/ 1/4" CERAMIC (30 ML TUBES) - 50 PACK**

**Sku: 19-6350Z**



### **GENERAL INFO**

#### **Hard Tissue Impact Grinding Mix w/ 1/4" Ceramic (30 mL Tubes) - 50 Pack**

Bead Ruptor bead beating tubes are pre-filled with an optimized amount of lysing matrix to provide efficient and high yield sample homogenates.

#### **Product Features:**

- 30 mL polypropylene screw cap tubes pre-filled with 2.4 mm metal beads and one 1/4" ceramic cylinder
- Recommended for: larger pieces of hard plant or animal tissues, dry grinding, roots and seeds

**Ordering Information:** Unit includes 50 bead beating tubes pre-filled with 2.4 mm metal beads and one 1/4" ceramic cylinder

### **SPECIFICATIONS**

<b>Tube material:</b>	Polypropylene
<b>Tube volume:</b>	30.0 mL
<b>Bead material:</b>	2.4 mm metal and one 1/4" ceramic cylinder
<b>Processing range:</b>	2 mL - 25 mL
<b>Temperature Range:</b>	-80 to 121°C
<b>Force:</b>	5,000 x g