



HARD TISSUE IMPACTION HOMOGENIZING MIX W/ 1/4" CERAMIC (30 ML TUBES) - 50 PACK

Sku: **19-6358Z**

GENERAL INFO

Hard Tissue Impaction Homogenizing Mix w/ 1/4" Ceramic (30 mL Tubes) - 50 Pack

Bead Ruptor bead beating tubes are pre-filled with an optimized amount of lysing matrix to provide efficient and high yield sample homogenates.

Product Features:

- 30 mL polypropylene screw cap tubes pre-filled with 2.8 mm ceramic beads and one 1/4" ceramic cylinder
- Recommended for: larger pieces of tough and fibrous plant material, hard tissues, seeds and feed

Ordering Information: Unit includes 50 bead beating tubes pre-filled with 2.8 mm ceramic beads and one 1/4" ceramic cylinder

SPECIFICATIONS

Tube material:	Polypropylene
Tube volume:	30.0 mL
Bead material:	2.8 mm ceramic tough and one 1/4" ceramic cylinder
Processing range :	2 mL - 25 mL
Temperature Range:	-80 to 121°C
Force:	5,000 x g