

LANDANGER	FICHE TECHNIQUE PRODUIT <i>PRODUCT TECHNICAL SHEET</i>	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	---	---

Éditée le / *Publishing date*

2022-04-28

Page

1

REFERENCE	DESIGNATIONS - <i>DESCRIPTIONS</i>	CODE EAN
B25320	CISEAUX POINTUS COURBES L:120MM CATARACTE	3661219004097
	SHARP SCISSORS CURVED L:120MM CATARACT	

Catalogues / *catalogs*

Landanger 2015

Chapitre-page

03-012

DESCRIPTIF DU DISPOSITIF/ *PRODUCT DESCRIPTION*

LANDANGER

Ciseaux délicats - Delicate scissors - Tijeras delicadas

STRABISME STRABISMUS ESTRABISMO <table border="0"> <tr> <td>Droit Straight Recta</td> <td>Courbe Curved Curva</td> </tr> </table> 12 cm - 4 3/4"  B25265 B25270	Droit Straight Recta	Courbe Curved Curva	<table border="0"> <tr> <td>Droit Straight Recta</td> <td>Courbe Curved Curva</td> </tr> </table>  B253680 B253730	Droit Straight Recta	Courbe Curved Curva	CATARACTE CATARACT CATARATA <table border="0"> <tr> <td>Droit Straight Recta</td> <td>Courbe Curved Curva</td> </tr> </table>  B25315 B25320	Droit Straight Recta	Courbe Curved Curva
Droit Straight Recta	Courbe Curved Curva							
Droit Straight Recta	Courbe Curved Curva							
Droit Straight Recta	Courbe Curved Curva							

DAUPHIN

Mousse / Mousse Blunt / Blunt Roma / Roma  B25270	Pointu / Mousse Sharp / Blunt Aguda / Roma  B253730	Pointu / Pointu Sharp / Sharp Aguda / Aguda  B25320
---	---	---

 B25265 12 cm - 4 3/4"	 B253680 12 cm - 4 3/4"	 B25315 12 cm - 4 3/4"
---	--	---

III - 12

III - Ciseaux
III - Scissors
III - Tijeras

Copyright © LANDANGER - 2009 - 2015

LANDANGER	FICHE TECHNIQUE PRODUIT <i>PRODUCT TECHNICAL SHEET</i>	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	---	---

Éditée le / *Publishing date*

2022-04-28

Page

2

REFERENCE	DESIGNATIONS - <i>DESCRIPTIONS</i>	CODE EAN
B25320	CISEAUX POINTUS COURBES L:120MM CATARACTE	3661219004097
	SHARP SCISSORS CURVED L:120MM CATARACT	 3 661219 004097

INFORMATIONS REGLEMENTAIRES / *REGULATORY INFORMATION*

Classe CE :
CE Class :

I

Code GMDN :
GMDN code :

38727

PERFORMANCES REVENDIQUEES / *INTENDED PERFORMANCE*

Dispositif médical de bloc opératoire :
Operating room medical device :

REUTILISABLE/REUSABLE

CARACTERISTIQUES SPECIFIQUE / *SPECIFIC FEATURES*

Type matière :
Material type :

ACIER INOXYDABLE / STAINLESS STEEL

Nuance matière :
Material grade :

X39Cr13(1.4031)

Dimensions : Voir désignations produits
Outer dimensions : *See products description*
Chapitre-pages

Landanger 2015

03-012

Nettoyage/Pré-désinfection : Standard pour instrumentation de bloc opératoire réutilisable.
Cleaning / Pre-disinfection : *Standard conditions for reusable surgical instruments*

Sterilisation :

La stérilisation doit être effectuée selon un procédé de stérilisation à la vapeur conformément à la norme ISO 17665-1 (procédé sous vide fractionné) avec un stérilisateur conformément à la norme EN 285.

Méthode de stérilisation à la vapeur saturée	Température	Durée	Temps de séchage	Pression absolue
Pré-vide	132 °C (270 °F)	4 minutes	30 minutes	2,9 bars, 2,8 atm, 41,6 psi
Pré-vide	134 °C (273 °F)	3 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	5 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	18 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi

Sterilization process :

The Sterilization shall be performed according to ISO 17665-1 validated steam sterilisation procedure (fractional vacuum procedure) in an autoclave in compliance with EN 285.

Saturated steam sterilization method	Temperature	Time	Drying time	Absolute pressure
Pre-vacuum	132 °C (270 °F)	4 minutes	30 minutes	2.9 bar, 2.8 atm, 41.6 psi
Pre-vacuum	134 °C (273 °F)	3 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	5 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	18 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi

EMBALLAGE / *PACKAGING*

Conditionnement : Pochette plastique en polyéthylène (PE).
Packaging : *Polyethylene plastic pouch (PE).*