

|                                                                                  |                     |                                                                                    |
|----------------------------------------------------------------------------------|---------------------|------------------------------------------------------------------------------------|
|  | SPECIFICATION SHEET | FAQ.026<br>Indice : 1<br>Page : 1<br>Date : 30/04/2002<br><br>Révisée : 19/02/2007 |
|                                                                                  | CE LABELLING        |                                                                                    |

## NAVARRO THORACIC RETRACTOR

| <u>ARTICLE CODE</u> | <u>DESIGNATION</u>                                                                                          |
|---------------------|-------------------------------------------------------------------------------------------------------------|
| <b>DC30650-00</b>   | NAVARRO THORACIC RETRACTOR L: 190MM W: 55MM MAX OPENING: 130MM INCLUDED DC30651-30 (PAIR) DC30652-30 (PAIR) |
| <b>DC30651-20</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:55MM D:20MM                                                          |
| <b>DC30651-30</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:55MM D:30MM                                                          |
| <b>DC30651-40</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:55MM D:40MM                                                          |
| <b>DC30652-20</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:27MM D:20MM                                                          |
| <b>DC30652-30</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:27MM D:30MM                                                          |
| <b>DC30653-30</b>   | FIXATION + ACCESSORY BLADE FOR NAVARRO RETRACTOR W:55MM D:30MM                                              |



### PERFORMANCES

| <u>ARTICLE CODE</u> | <u>PERFORMANCES</u>                                |
|---------------------|----------------------------------------------------|
| <b>DC30650-00</b>   | THORACIC RETRACTOR FOR A CONICAL FIELD             |
| <b>DC30651-20</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:55MM D:20MM |
| <b>DC30651-30</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:55MM D:30MM |
| <b>DC30651-40</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:55MM D:40MM |
| <b>DC30652-20</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:27MM D:20MM |
| <b>DC30652-30</b>   | PAIR OF BLADES FOR NAVARRO RETRACTOR W:27MM D:30MM |

|                                                                                  |                     |  |                                                                     |
|----------------------------------------------------------------------------------|---------------------|--|---------------------------------------------------------------------|
|  | SPECIFICATION SHEET |  | FAQ.026                                                             |
|                                                                                  | CE LABELLING        |  | Indice : 1<br>Page : 2<br>Date : 30/04/2002<br>Révisée : 19/02/2007 |

**DC30653-30**

*FIXATION + ACCESSORY BLADE FOR NAVARRO RETRACTOR W:55MM  
D:30MM*

### **SPECIFICATIONS**

Components :      *Metallic parts :*

- *Austenitic stainless steel Z06CN18.09 (1.4301)*
- *Martensitic stainless steel Z30C13 (1.4028), tempered.*

Dimensions :      *According to product descriptions.*

### **INSTRUCTIONS FOR USE**

*See Instructions for use if provided.*

### **PREDISINFECTION**

*Use neutral pre-disinfection bath.  
Incompatible with chloride.  
Compatible with soda process.*

### **CLEANING**

*Recommended : autowasher cleaning  
Compatible with ultrasound cleaning process.*

### **STERILIZATION**

Sterilization :      *Saturated steam at a maximum 135°C / 20 minutes/ 2.2 bars.*