

LANDANGER	FICHE TECHNIQUE PRODUIT PRODUCT TECHNICAL SHEET	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	--	---

Éditée le / Publishing date

2022-04-28

Page

1

REFERENCE	DESIGNATIONS - DESCRIPTIONS	CODE EAN
DC41103-66TI	CLAMP BULLDOG DE BAKEY MORS 66MM L:125MM DROIT TITANE	3661219265603
	DE BAKEY BULLDOG CLAMP JAWS : 66MM L:125MM STRAIGHT TI	

Catalogues / catalogs

Delacroix 2020

Chapitre-page

04.1-012

DESCRIPTIF DU DISPOSITIF/ PRODUCT DESCRIPTION



DeBakey Bulldog Clamps / Clamps Bulldog DeBakey / Clamps Bulldog DeBakey

Pressure / Pression / Presión : 250 g



DeBakey

Straight
Droit
Recta





Curved
Courbe
Curva

DC41100-23 

L : 80 mm - 3 1/8"

Jaws / Mors / Bocas : 23 mm

DC41101-29 

L : 88 mm - 3 1/2"

Jaws / Mors / Bocas : 29 mm

DC41102-50 

L : 110 mm - 4 1/4"

Jaws / Mors / Bocas : 50 mm

 **DC41103-66** 

L : 125 mm - 4 7/8"

Jaws / Mors / Bocas : 66 mm

DC41104-70 

L : 130 mm - 5 1/8"

Jaws / Mors / Bocas : 70 mm

DC41105-21

L : 77 mm - 3"

Jaws / Mors / Bocas : 21 mm

DC41106-23 

L : 80 mm - 3 1/8"

Jaws / Mors / Bocas : 23 mm

DC41107-32 

L : 90 mm - 3 1/2"

Jaws / Mors / Bocas : 32 mm

DC41108-49 

L : 105 mm - 4 1/8"

Jaws / Mors / Bocas : 49 mm

 **DC41109-64** 

L : 123 mm - 4 7/8"

Jaws / Mors / Bocas : 64 mm

Ⓟ : Add TI for Titanium Model / Ajoutez TI pour le Modèle Titane / Añadir TI para Modelo Titánico

12



Copyright LANDANGER - 2009 - 2020

LANDANGER	FICHE TECHNIQUE PRODUIT <i>PRODUCT TECHNICAL SHEET</i>	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	---	---

Éditée le / *Publishing date*

2022-04-28

Page

2

REFERENCE	DESIGNATIONS - <i>DESCRIPTIONS</i>	CODE EAN
DC41103-66TI	CLAMP BULLDOG DE BAKEY MORS 66MM L:125MM DROIT TITANE	3661219265603
	DE BAKEY BULLDOG CLAMP JAWS : 66MM L:125MM STRAIGHT TI	 3 661219 265603

INFORMATIONS REGLEMENTAIRES / *REGULATORY INFORMATION*

Classe CE :

I

Code GMDN :

62470

CE Class :

GMDN code :

PERFORMANCES REVENDIQUEES / *INTENDED PERFORMANCE*

Dispositif médical de bloc opératoire :

REUTILISABLE/REUSABLE

Operating room medical device :

CARACTERISTIQUES SPECIFIQUE / *SPECIFIC FEATURES*

Type matière :

TITANE / TITANIUM

Material type :

Nuance matière :

TA6V(TITANE)

Material grade :

Dimensions :

Voir désignations produits

Delacroix 2020

Outer dimensions :

See products description

Chapitre-pages

04.1-012

Nettoyage/Pré-désinfection : Standard pour instrumentation de bloc opératoire réutilisable.

Cleaning / Pre-disinfection : *Standard conditions for reusable surgical instruments*

Stérilisation :

La stérilisation doit être effectuée selon un procédé de stérilisation à la vapeur conformément à la norme ISO 17665-1 (procédé sous vide fractionné) avec un stérilisateur conformément à la norme EN 285.

Méthode de stérilisation à la vapeur saturée	Température	Durée	Temps de séchage	Pression absolue
Pré-vide	132 °C (270 °F)	4 minutes	30 minutes	2,9 bars, 2,8 atm, 41,6 psi
Pré-vide	134 °C (273 °F)	3 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	5 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	18 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi

Sterilization process :

The Sterilization shall be performed according to ISO 17665-1 validated steam sterilisation procedure (fractional vacuum procedure) in an autoclave in compliance with EN 285.

Saturated steam sterilization method	Temperature	Time	Drying time	Absolute pressure
Pre-vacuum	132 °C (270 °F)	4 minutes	30 minutes	2.9 bar, 2.8 atm, 41.6 psi
Pre-vacuum	134 °C (273 °F)	3 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	5 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	18 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi

EMBALLAGE / *PACKAGING*

Conditionnement :

Pochette plastique en polyéthylène (PE).

Packaging :

Polyethylene plastic pouch (PE).