

L A N D A N G E R	FICHE TECHNIQUE PRODUIT PRODUCT TECHNICAL SHEET	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-------------------	--	--

Éditée le / Publishing date

2022-04-28

Page

1

REFERENCE	DESIGNATIONS - DESCRIPTIONS	CODE EAN
F36410	MANDRIN DILATATEUR HEGAR DOUBLE 3X4MM	3661219069713
	HEGAR UTERINE DILATOR DOUBLE EDGED STAINLESS STEEL 3X4MM	 3 661219 069713

Catalogues / catalogs

Landanger 2015

Chapitre-page

08-084

DESCRIPTIF DU DISPOSITIF/ PRODUCT DESCRIPTION

L A N D A N G E R

Dilatateurs utérins - Uterine dilators - Dilataores uterinos

HEGAR	HANK	PRATT																																																																			
Modèle double Double model Modelo doble																																																																					
 <table border="0"> <tr><td></td><td>Ø</td></tr> <tr><td>F36400</td><td>1 x 2 mm</td></tr> <tr><td>F36410</td><td>3 x 4 mm</td></tr> <tr><td>F36420</td><td>5 x 6 mm</td></tr> <tr><td>F36430</td><td>7 x 8 mm</td></tr> <tr><td>F36440</td><td>9 x 10 mm</td></tr> <tr><td>F36450</td><td>11 x 12 mm</td></tr> <tr><td>F36460</td><td>13 x 14 mm</td></tr> <tr><td>F36470</td><td>15 x 16 mm</td></tr> <tr><td>F36480</td><td>17 x 18 mm</td></tr> <tr><td>F36485</td><td>19 x 20 mm</td></tr> <tr><td>F36490</td><td>21 x 22 mm</td></tr> <tr><td>F36494</td><td>23 x 24 mm</td></tr> <tr><td>F36498</td><td>25 x 26 mm</td></tr> </table>		Ø	F36400	1 x 2 mm	F36410	3 x 4 mm	F36420	5 x 6 mm	F36430	7 x 8 mm	F36440	9 x 10 mm	F36450	11 x 12 mm	F36460	13 x 14 mm	F36470	15 x 16 mm	F36480	17 x 18 mm	F36485	19 x 20 mm	F36490	21 x 22 mm	F36494	23 x 24 mm	F36498	25 x 26 mm	 <table border="0"> <tr><td></td><td>Ø</td><td>Fig. US</td></tr> <tr><td>F365050</td><td>4,5 x 5 mm</td><td>9/10</td></tr> <tr><td>F365060</td><td>5,5 x 6 mm</td><td>11/12</td></tr> <tr><td>F365070</td><td>6,5 x 7 mm</td><td>13/14</td></tr> <tr><td>F365080</td><td>7,5 x 8 mm</td><td>15/16</td></tr> <tr><td>F365090</td><td>8,5 x 9 mm</td><td>17/18</td></tr> <tr><td>F365100</td><td>9,5 x 10 mm</td><td>19/20</td></tr> </table>		Ø	Fig. US	F365050	4,5 x 5 mm	9/10	F365060	5,5 x 6 mm	11/12	F365070	6,5 x 7 mm	13/14	F365080	7,5 x 8 mm	15/16	F365090	8,5 x 9 mm	17/18	F365100	9,5 x 10 mm	19/20	 <table border="0"> <tr><td></td><td>Ø</td></tr> <tr><td>F365200</td><td>4,5 x 5,0 mm</td></tr> <tr><td>F365201</td><td>5,6 x 6,3 mm</td></tr> <tr><td>F365202</td><td>7,0 x 7,6 mm</td></tr> <tr><td>F365203</td><td>8,3 x 9,0 mm</td></tr> <tr><td>F365204</td><td>9,6 x 10,3 mm</td></tr> <tr><td>F365205</td><td>11,0 x 11,6 mm</td></tr> <tr><td>F365206</td><td>12,3 x 13,0 mm</td></tr> <tr><td>F365207</td><td>13,6 x 14,3 mm</td></tr> </table>		Ø	F365200	4,5 x 5,0 mm	F365201	5,6 x 6,3 mm	F365202	7,0 x 7,6 mm	F365203	8,3 x 9,0 mm	F365204	9,6 x 10,3 mm	F365205	11,0 x 11,6 mm	F365206	12,3 x 13,0 mm	F365207	13,6 x 14,3 mm
	Ø																																																																				
F36400	1 x 2 mm																																																																				
F36410	3 x 4 mm																																																																				
F36420	5 x 6 mm																																																																				
F36430	7 x 8 mm																																																																				
F36440	9 x 10 mm																																																																				
F36450	11 x 12 mm																																																																				
F36460	13 x 14 mm																																																																				
F36470	15 x 16 mm																																																																				
F36480	17 x 18 mm																																																																				
F36485	19 x 20 mm																																																																				
F36490	21 x 22 mm																																																																				
F36494	23 x 24 mm																																																																				
F36498	25 x 26 mm																																																																				
	Ø	Fig. US																																																																			
F365050	4,5 x 5 mm	9/10																																																																			
F365060	5,5 x 6 mm	11/12																																																																			
F365070	6,5 x 7 mm	13/14																																																																			
F365080	7,5 x 8 mm	15/16																																																																			
F365090	8,5 x 9 mm	17/18																																																																			
F365100	9,5 x 10 mm	19/20																																																																			
	Ø																																																																				
F365200	4,5 x 5,0 mm																																																																				
F365201	5,6 x 6,3 mm																																																																				
F365202	7,0 x 7,6 mm																																																																				
F365203	8,3 x 9,0 mm																																																																				
F365204	9,6 x 10,3 mm																																																																				
F365205	11,0 x 11,6 mm																																																																				
F365206	12,3 x 13,0 mm																																																																				
F365207	13,6 x 14,3 mm																																																																				
 F36440 20 cm - 8"	 F365050 27,5 cm - 10 3/4"	 F365203 28 cm - 11"																																																																			

VIII - Voies basses
VIII - Gynae Utero Procto
VIII - Gynae Utero Procto

VIII - 84

Copyright © LANDANGER - 2009 - 2015

LANDANGER	FICHE TECHNIQUE PRODUIT <i>PRODUCT TECHNICAL SHEET</i>	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	---	---

Éditée le / *Publishing date*

2022-04-28

Page

2

REFERENCE	DESIGNATIONS - <i>DESCRIPTIONS</i>	CODE EAN
F36410	MANDRIN DILATATEUR HEGAR DOUBLE 3X4MM	3661219069713
	HEGAR UTERINE DILATOR DOUBLE EDGED STAINLESS STEEL 3X4MM	 3 661219 069713

INFORMATIONS REGLEMENTAIRES / *REGULATORY INFORMATION*

Classe CE :

I-MDR

Code GMDN :

11267

CE Class :

GMDN code :

PERFORMANCES REVENDIQUEES / *INTENDED PERFORMANCE*

Dispositif médical de bloc opératoire :

REUTILISABLE/REUSABLE

Operating room medical device :

CARACTERISTIQUES SPECIFIQUE / *SPECIFIC FEATURES*

Type matière :

ACIER INOXYDABLE / STAINLESS STEEL

Material type :

Nuance matière :

X2CrNi19-11(1.4306)

Material grade :

Dimensions :

Voir désignations produits

Landanger 2015

Outer dimensions :

See products description

Chapitre-pages

08-084

Nettoyage/Pré-désinfection : Standard pour instrumentation de bloc opératoire réutilisable.

Cleaning / Pre-disinfection : *Standard conditions for reusable surgical instruments*

Stérilisation :

La stérilisation doit être effectuée selon un procédé de stérilisation à la vapeur conformément à la norme ISO 17665-1 (procédé sous vide fractionné) avec un stérilisateur conformément à la norme EN 285.

Méthode de stérilisation à la vapeur saturée	Température	Durée	Temps de séchage	Pression absolue
Pré-vide	132 °C (270 °F)	4 minutes	30 minutes	2,9 bars, 2,8 atm, 41,6 psi
Pré-vide	134 °C (273 °F)	3 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	5 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	18 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi

Sterilization process :

The Sterilization shall be performed according to ISO 17665-1 validated steam sterilisation procedure (fractional vacuum procedure) in an autoclave in compliance with EN 285.

Saturated steam sterilization method	Temperature	Time	Drying time	Absolute pressure
Pre-vacuum	132 °C (270 °F)	4 minutes	30 minutes	2.9 bar, 2.8 atm, 41.6 psi
Pre-vacuum	134 °C (273 °F)	3 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	5 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	18 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi

EMBALLAGE / *PACKAGING*

Conditionnement :

Pochette plastique en polyéthylène (PE).

Packaging :

Polyethylene plastic pouch (PE).