

L A N D A N G E R	FICHE TECHNIQUE PRODUIT PRODUCT TECHNICAL SHEET	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-------------------	--	--

Editée le / Publishing date

2022-04-29

Page

1

REFERENCE	DESIGNATIONS - DESCRIPTIONS	CODE EAN
J163159	SET COMPLET DE 6 CURETTES ANNEAU L:320MM BLACK MAGNUM MANCHE DIAM 28MM L:198MM	3661219443834
	COMPLETE SET WITH 6 RING CURETTES L:320MM BLACK MAGNUM HANDLE DIAM 28MM L:198MM	 3 661219 443834

Catalogues / catalogs

OrthoVertebre 202x

Chapitre-page

J16314x

DESCRIPTIF DU DISPOSITIF/ PRODUCT DESCRIPTION

L A N D A N G E R

Curettes anneaux - Ring curettes - Curetas de anillo



Magnum

Dimensions du manche
Handle dimensions
Dimensiones del mango:

	Ø 28mm L: 198mm	Ø 32mm L: 221mm
Ø : 6 mm	J163143	J163140
Ø : 8 mm	J163144	J163141
Ø : 10 mm	J163145	J163142



Ø 28mm
L: 198mm

Ø 32mm
L: 221mm

	Ø 28mm L: 198mm	Ø 32mm L: 221mm
Ø : 6 mm	J163153	J163150
Ø : 8 mm	J163154	J163151
Ø : 10 mm	J163155	J163152



Set complet
Complete set
Juego completo

Dimensions du manche
Handle dimensions
Dimensiones del mango:

	Ø 28mm - L: 198mm	Ø 32mm - L: 221mm
	J163159	J163158

Copyright LANDANGER - 2019

347

IX - Divers
IX - Miscellaneous
IX - Diversos

LANDANGER	FICHE TECHNIQUE PRODUIT <i>PRODUCT TECHNICAL SHEET</i>	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	---	---

Éditée le / *Publishing date*

2022-04-29

Page

2

REFERENCE	DESIGNATIONS - <i>DESCRIPTIONS</i>	CODE EAN
J163159	SET COMPLET DE 6 CURETTES ANNEAU L:320MM BLACK MAGNUM MANCHE DIAM 28MM L:198MM	3661219443834
	COMPLETE SET WITH 6 RING CURETTES L:320MM BLACK MAGNUM HANDLE DIAM 28MM L:198MM	 3 661219 443834

INFORMATIONS REGLEMENTAIRES / *REGULATORY INFORMATION*

Classe CE :
CE Class :

I

Code GMDN :
GMDN code :

31335

PERFORMANCES REVENDIQUEES / *INTENDED PERFORMANCE*

Dispositif médical de bloc opératoire :
Operating room medical device :

REUTILISABLE/REUSABLE

CARACTERISTIQUES SPECIFIQUE / *SPECIFIC FEATURES*

Type matière :
Material type :

ACIER INOXYDABLE + ALUMINIUM / STAINLESS STEEL + ALUMINIUM

Nuance matière :
Material grade :

X20Cr13(1.4021) + X8CrNiS18-9(1.4305) + AlMgSi1

Dimensions : Voir désignations produits
Outer dimensions : See products description
Chapitre-pages

OrthoVertebre 202x

J16314x

Nettoyage/Pré-désinfection : Standard pour instrumentation de bloc opératoire réutilisable.
Cleaning / Pre-disinfection : Standard conditions for reusable surgical instruments

Sterilisation :

La stérilisation doit être effectuée selon un procédé de stérilisation à la vapeur conformément à la norme ISO 17665-1 (procédé sous vide fractionné) avec un stérilisateur conformément à la norme EN 285.

Méthode de stérilisation à la vapeur saturée	Température	Durée	Temps de séchage	Pression absolue
Pré-vide	132 °C (270 °F)	4 minutes	30 minutes	2,9 bars, 2,8 atm, 41,6 psi
Pré-vide	134 °C (273 °F)	3 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	5 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	18 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi

Sterilization process :

The Sterilization shall be performed according to ISO 17665-1 validated steam sterilisation procedure (fractional vacuum procedure) in an autoclave in compliance with EN 285.

Saturated steam sterilization method	Temperature	Time	Drying time	Absolute pressure
Pre-vacuum	132 °C (270 °F)	4 minutes	30 minutes	2.9 bar, 2.8 atm, 41.6 psi
Pre-vacuum	134 °C (273 °F)	3 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	5 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	18 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi

EMBALLAGE / *PACKAGING*

Conditionnement : Pochette plastique en polyéthylène (PE).
Packaging : Polyethylene plastic pouch (PE).