

Picric acid aqueous 1.2% solution

IVD In-vitro diagnostic medical device **CE**

CND Code: W01030799

Catalog number	Unit size
05-M05027	500 ml

Packaging Primary container: white bottle in High Density Polyethylene (HDPE). Useful capacity 500 ml. HDPE cap. Tamper evident cap.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Expected aim Product for the preparation of cyto-histological samples for optical microscopy.

Application Component used in the preparation of several fixative mixtures, the Bouin fixative included.

Components

Components	CAS	CE	Index
Picric acid	88-89-1	2018659	609-009-00-X
Deionized water			

Warning and precaution

The product must be used exclusively by specialized technical operators.
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release.
Do not use if the primary container is damaged.

Storage Store the preparation at room temperature. Keep the containers tightly closed.

Stability After the first opening, the product is usable until the expiry date, if correctly stored. Product validity: 2 years.

Disposal Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

Date of issue: May 2013