

EN

Tubing maintenance sol.

REF IK705 6 x 20 ml

INTENDED USE

Solution dedicated to The Binding Site Optilite® analyzer.
To be used according to Operation Manual instructions, during monthly maintenance procedure.

PRODUCT INFORMATION

6 x 20 ml
Ready to use liquid.

Composition

Detergent, Proclin

Precautions

Exercise the normal precautions required for handling all laboratory reagents.
EUH208: Contains Isothiazolinone compounds. May produce an allergic reaction.

Storage and Stability

The Tubing maintenance solution is stable at 2...25 °C until the expiration date printed on the label.

MANUFACTURER

Thermo Fisher Scientific Oy
Ratastie 2, FI-01620 Vantaa, Finland
Tel. +358 10 329 200
www.thermofisher.com

Date of revision (yyyy-mm-dd)

2019-07-26

Changes from previous version

Detergent changed.
Manufacturer contact information updated.

Template: D01630_4 © 2019 Thermo Fisher Scientific Inc. All rights reserved. In the European Union Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, UK.. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.



The Binding Site Group Ltd

B15 1QT, United Kingdom.

8 Calthorpe Road, Edgbaston, Birmingham,

Tel: +44 (0) 121 456 9500
info@bindingsite.co.uk

Fax: +44 (0) 121 456 9749
www.bindingsite.com

