

Terumo® Surflo® Winged Infusion Set


Product Information Card



Description:

The Surflo® winged infusion set (= Scalp Vein or SV-set) is intended to introduce or withdraw fluids from the vascular system after vein puncture.

Components:



Description	Raw Materials
1. Cannula	Stainless steel
2. Wing hub	Polyvinyl chloride + master batch
3. Tube	Polyvinyl chloride
4. Adapter: BL item code EL/BLK item code	Polypropylene PMMA
5. Adapter cap: BL item codes EL/BLK item code	Polyethylene + master batch HDPE + master batch
Protector	Polyethylene



Product Range:

Product Code	Colour	External Cannula Ø	Cannula Length	Dead Space Volume	Tubing Length
SV*18BLK03	Pink	18G 1.20 mm	19 mm	0.46 ml	300 mm
SV-19BL	Ivory	19G 1.10 mm	19 mm	0.43 ml	300 mm
SV-21BL	Green	21G 0.80 mm	19 mm	0.42 ml	300 mm
SV-22BL	Black	22G 0.70 mm	19 mm	0.42 ml	300 mm
SV-23BL	Blue	23G 0.60 mm	19 mm	0.41 ml	300 mm
SV-25BL	Orange	25G 0.50 mm	19 mm	0.41 ml	300 mm
SV*27EL	Grey	27G 0.40 mm	19 mm	0.27 ml	200 mm

Terumo® Surflo® Winged Infusion Set

Packaging:

Single unit: blister pack (thermoformed blister film)

Product Code	Unit Box	Shipping Carton	Gross Weight	Dimensions (LxWxH)
SV*18BLK03	50	500	3.2 kg	430 x 220 x 265 mm
SV-19BL	50	500	2.6 kg	438 x 228 x 276 mm
SV-21BL	50	500	2.6 kg	438 x 228 x 276 mm
SV-22BL	50	500	2.6 kg	438 x 228 x 276 mm
SV-23BL	50	500	2.6 kg	438 x 228 x 276 mm
SV-25BL	50	500	2.6 kg	438 x 228 x 276 mm
SV*27EL	100	1000	n.a.	475 x 255 x 465 mm

Sterilization:

With ethylene oxide to a Sterility Assurance Level of at least 10^{-6} for all BL/EL-items
With gamma radiation to a Sterility Assurance Level of at least 10^{-6} for BLK-item

Storage:

Shelf life of 5 years for all BL-items / 3 years for all EL/BLK-items

* The lot number and expiry date are embossed in the film of the thermoformed blister pack

* Do not store at extreme temperature and humidity

Precautions for Use:

Check integrity of the unit package before use
For single use only
Dispose of safely after single use to avoid risk of infection

Status of Regulation:

Class IIa, medical device following MDD 93/42/EEC

Notified body: TÜV Rheinland 0197
Manufacturer: TERUMO EUROPE N.V.
Interleuvenlaan 40
B - 3001 Leuven
Belgium

TERUMO CORPORATION JAPAN (for EL/BLK item code)



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

European Offices

**Terumo Europe N.V.
Benelux Sales Division**
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 39 25 80
Fax: +32 16 39 25 99
The Netherlands:
Tel.: 0800 0220396
Fax: 0800 0220414

Laboratoires Terumo France S.A.
Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex, France
Tel.: +33 130 96 13 00
Fax: +33 130 43 60 85

Terumo Deutschland GmbH
Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel.: +49 6196 80 230
Fax: +49 6196 80 23 200

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel.: +41 56 419 10 10
Fax: +41 56 419 10 11

Terumo Europe España SL
Avda. Juan Carlos I, Nº 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares (Madrid),
Spain
Tel.: +34 9021 01 298
Fax: +34 9021 01 358

Terumo Italia SRL
Via Simone Martini, 143 / 145
00142 Roma RM,
Italy
Tel.: +39 0651 96 14 20
Fax: +39 0650 30 407

Terumo Sweden AB
Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
426 71 Västra Frölunda,
Sweden
Tel.: +46 3174 85 880
Fax: +46 3174 85 890

**Terumo Denmark
Filial of Terumo Sweden AB**
Tel.: +45 7020 93 80
Fax: +45 7020 94 80

Terumo UK Ltd
Tamesis
The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel.: +44 1784 476 220
Fax: +44 1784 476 201

**Terumo Europe N.V.
Moscow Representative Office**
BC Northern Tower
13Fr., 10 Testovskaya Street
Moscow 123317, Russia
Tel.: +7 495 988 4740
Fax: +7 495 988 4739

**Terumo Corporation Dubai Branch
(Middle East and Africa)**
Al Masraf Tower, 22nd Floor
P. O. Box 20291 Dubai
U.A.E.
Tel.: +971 4 221 2220
Fax: +971 4 221 3330



**Addresses subject to change,
please consult our Website.**

HO 29GB-0307TE-III(1.12)E

© Registered Trademark
Published by Terumo Europe N.V.